

## VWR® STAINLESS-STEEL TEST-TUBE BASKETS

As per the manufacturer, the below product meets the following criteria:

<b>VWR Catalog Number</b>	89259-822
<b>Lot Number</b>	1455
<b>Description</b>	VWR Basket S Steel Test Tube
<b>Date of Manufacture</b>	December 07, 2017
<b>Country of Origin</b>	PAKISTAN

### Quality System Compliance

Products are manufactured under the ISO 9001:2008 and ISO13485:2003 & F-899 standard. Products are Inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP

### QC Testing

Following tests performs as per procedure:

- Material Chemical Analysis 7-element Test.
- Boil Test
- Passivation

### Product Specifications

<b>Material:</b>	AISI 410 TYPE
<b>Sterilization Process:</b>	N/A
<b>BPA statement:</b>	N/A.
<b>BSE/TSE Statement:</b>	All Products are not of animal origin and do not contain, make use of, or involve, at any point of the manufacturing processes
<b>Non-Pyrogenic:</b>	N/A
<b>DNase &amp; RNase Free:</b>	N/A
<b>Cytotoxicity:</b>	N/A
<b>REACH statement:</b>	N/A
<b>Latex-Free Statement:</b>	This product is latex free.
<b>Storage Condition:</b>	Keep Dry Place.

Signed:



Jamie Ethier  
VP Global Quality  
VWR, Part of Avantor

**Date:** December 31, 2019