

Part of Thermo Fisher Scientific

Material Safety Data Sheet

Product and company identification

Product name : Dulbecco's Phosphate Buffered Saline

Synonym : Wash Buffer Powder

Supplier : Thermo Fisher Scientific Manufacturer

Pierce Biotechnology P.O. Box 117 Rockford, IL 61105 United States 815.968.0747 or 800.874.3723

Rockford, IL 61105 United States 815.968.0747 800.874.3723

P.O. Box 117

: Thermo Fisher Scientific

Pierce Biotechnology

Product No. : 0028374 1858094 1887840

MSDS # : 1093 Validation date : 11/4/2009. Print date : 11/4/2009.

Responsible name : MSDS (Regulatory Affairs)

In case of emergency : CHEMTREC: Use of

800.424.9300 Substance/Preparation

OUTSIDE US: 202.483.7616

 Refer to the instruction booklet for proper and intended use.
 Otherwise, contact supplier for specific applications.

2. Hazards identification

Physical state : Solid. [Crystalline powder.]

OSHA/HCS status : This material is considered hazardous by the OSHA Hazard Communication Standard

(29 CFR 1910.1200).

Emergency overview : WARNING!

CAUSES RESPIRATORY TRACT. EYE AND SKIN IRRITATION. MAY BE HARMFUL

IF SWALLOWED. CAN CAUSE TARGET ORGAN DAMAGE.

Harmful if swallowed. Irritating to eyes, respiratory system and skin. Handling and/or processing of this material may generate a dust which can cause mechanical irritation of the eyes, skin, nose and throat. Do not ingest. Avoid breathing dust. Avoid contact with eyes, skin and clothing. Can cause target organ damage. Use only with adequate ventilation. Keep container tightly closed and sealed until ready for use. Wash

thoroughly after handling.

Routes of entry : Eye contact. Inhalation. Ingestion.

Potential acute health effects

Inhalation : Irritating to respiratory system.

Ingestion : Harmful if swallowed.

Skin : Irritating to skin.

Eyes : Irritating to eyes.

Potential chronic health effects

Chronic effects : Can cause target organ damage.

Carcinogenicity: No known significant effects or critical hazards.

11/4/2009.

ences Division 3747 N. Meridian Rd. Rockford, IL (815) 968-0747 www.thermo.com/pierce



Dulbecco's Phosphate Buffered Saline

2. Hazards identification

 Mutagenicity
 : No known significant effects or critical hazards.

 Teratogenicity
 : No known significant effects or critical hazards.

 Developmental effects
 : No known significant effects or critical hazards.

 Fertility effects
 : No known significant effects or critical hazards.

 Target organs
 : Causes damage to the following organs: heart.

Contains material which causes damage to the following organs: mucous membranes,

gastrointestinal tract, skin, eye, lens or cornea.

Contains material which may cause damage to the following organs: stomach.

Over-exposure signs/symptoms

Inhalation : Adverse symptoms may include the following:

respiratory tract irritation

coughing

Ingestion : No specific data.

Skin : Adverse symptoms may include the following:

irritation redness

Eyes : Adverse symptoms may include the following:

pain or irritation watering redness

Medical conditions aggravated by over-

: Pre-existing digestive disorders and disorders involving any other target organs mentioned in this MSDS as being at risk may be aggravated by over-exposure to this

exposure prod

The preparation is not classified as dangerous according to Directive 1999/45/EC and its amendments.

Classification : Not classified.

See toxicological information (section 11)

3. Composition/information on ingredients

United States

 Name
 CAS number
 %

 Sodium Chloride
 7647-14-5
 80 - 95

 Sodium Phosphate, Dibasic
 7558-79-4
 10 - 20

 Potassium Chloride
 7447-40-7
 1 - 3

 Potassium Phosphate, Monobasic
 7778-77-0
 1 - 3

Europe

Substance/preparation : Preparation

 Ingredient name
 CAS number
 %
 EC number
 Classification

 Sodium Phosphate, Dibasic
 7558-79-4
 10 - 20
 231-448-7
 Xi;
 [1]

 R36/37/38

There are no ingredients or additional ingredients present which, within the current knowledge of the supplier and in the concentrations applicable, are classified as hazardous to health or the environment and hence require reporting in this section.

[1] Substance classified with a health or environmental hazard

[2] Substance with a workplace exposure limit

11/4/2009.

2/9

 Biosciences Division
 3747 N. Meridian Rd.
 Rockford, IL
 (815) 968-0747
 www.thermo.com/pi

 Pierce
 PO Box 117
 61105
 (815) 968-7316 Fax



Composition/information on ingredients

Occupational exposure limits, if available, are listed in section 8.

First aid measures

Inhalation

: Move exposed person to fresh air. If it is suspected that fumes are still present, the rescuer should wear an appropriate mask or self-contained breathing apparatus. Keep person warm and at rest. If not breathing, if breathing is irregular or if respiratory arrest occurs, provide artificial respiration or oxygen by trained personnel. It may be dangerous to the person providing aid to give mouth-to-mouth resuscitation. Get medical attention. If unconscious, place in recovery position and get medical attention immediately. Maintain an open airway. Loosen tight clothing such as a collar, tie, belt or waistband.

Ingestion

: Wash out mouth with water. Remove dentures if any. Move exposed person to fresh air. Keep person warm and at rest. If material has been swallowed and the exposed person is conscious, give small quantities of water to drink. Stop if the exposed person feels sick as vomiting may be dangerous. Do not induce vomiting unless directed to do so by medical personnel. If vomiting occurs, the head should be kept low so that vomit does not enter the lungs. Get medical attention. Never give anything by mouth to an unconscious person. If unconscious, place in recovery position and get medical attention immediately. Maintain an open airway. Loosen tight clothing such as a collar, tie, belt or waistband.

Skin contact

: Flush contaminated skin with plenty of water. Remove contaminated clothing and shoes. Continue to rinse for at least 10 minutes. Get medical attention. Wash clothing before reuse. Clean shoes thoroughly before reuse.

Eye contact

: Immediately flush eyes with plenty of water for at least 15 minutes, occasionally lifting the upper and lower eyelids. Check for and remove any contact lenses. Get medical

Protection of first-aiders

: No action shall be taken involving any personal risk or without suitable training. If it is suspected that fumes are still present, the rescuer should wear an appropriate mask or self-contained breathing apparatus. It may be dangerous to the person providing aid to give mouth-to-mouth resuscitation.

See section 11 for more detailed information on health effects and symptoms.

5. Fire-fighting measures

Flammability of the product : No specific fire or explosion hazard.

Extinguishing media Suitable

: Use an extinguishing agent suitable for the surrounding fire.

Not suitable

Special exposure hazards

: Promptly isolate the scene by removing all persons from the vicinity of the incident if there is a fire. No action shall be taken involving any personal risk or without suitable

Hazardous combustion

products

: Decomposition products may include the following materials:

phosphorus oxides halogenated compounds

metal oxide/oxides

Special protective equipment for fire-fighters : Fire-fighters should wear appropriate protective equipment and self-contained breathing apparatus (SCBA) with a full face-piece operated in positive pressure mode.

Thermo

Dulbecco's Phosphate Buffered Saline

6. Accidental release measures

Personal precautions

No action shall be taken involving any personal risk or without suitable training. Evacuate surrounding areas. Keep unnecessary and unprotected personnel from entering. Do not touch or walk through spilled material. Avoid breathing dust. Provide adequate ventilation. Wear appropriate respirator when ventilation is inadequate. Put on appropriate personal protective equipment (see section 8).

Environmental precautions

Avoid dispersal of spilled material and runoff and contact with soil, waterways, drains and sewers. Inform the relevant authorities if the product has caused environmental pollution (sewers, waterways, soil or air).

Large spill

: Move containers from spill area. Approach release from upwind. Prevent entry into sewers, water courses, basements or confined areas. Vacuum or sweep up material and place in a designated, labeled waste container. Avoid creating dusty conditions and prevent wind dispersal. Dispose of via a licensed waste disposal contractor. Note: see section 1 for emergency contact information and section 13 for waste disposal.

Small spill

: Move containers from spill area. Vacuum or sweep up material and place in a designated, labeled waste container. Dispose of via a licensed waste disposal

Handling and storage

Handling

Put on appropriate personal protective equipment (see section 8). Eating, drinking and smoking should be prohibited in areas where this material is handled, stored and processed. Workers should wash hands and face before eating, drinking and smoking. Do not breathe dust. Do not ingest. Avoid contact with eyes, skin and clothing. Use only with adequate ventilation. Wear appropriate respirator when ventilation is inadequate. Keep in the original container or an approved alternative made from a compatible material, kept tightly closed when not in use. Empty containers retain product residue and can be hazardous. Do not reuse container

Storage

: Store between the following temperatures: 20 to 25°C (68 to 77°F). Store in accordance with local regulations. Store in original container protected from direct sunlight in a dry, cool and well-ventilated area, away from incompatible materials (see section 10) and food and drink. Keep container tightly closed and sealed until ready for use. Containers that have been opened must be carefully resealed and kept upright to prevent leakage. Do not store in unlabeled containers. Use appropriate containment to avoid environmental contamination.

Exposure controls/personal protection

Europe

No exposure limit value known.

Consult local authorities for acceptable exposure limits. Recommended monitoring

procedures

: If this product contains ingredients with exposure limits, personal, workplace atmosphere or biological monitoring may be required to determine the effectiveness of the ventilation or other control measures and/or the necessity to use respiratory protective equipment.

Engineering measures

: Use only with adequate ventilation. If user operations generate dust, fumes, gas, vapor or mist, use process enclosures, local exhaust ventilation or other engineering controls to keep worker exposure to airborne contaminants below any recommended or statutory

3/9 4/9 11/4/2009. 11/4/2009 3747 N. Meridian Rd. PO Box 117 Rockford, IL 61105



Exposure controls/personal protection

Hygiene measures

; Wash hands, forearms and face thoroughly after handling chemical products, before eating, smoking and using the lavatory and at the end of the working period. Appropriate techniques should be used to remove potentially contaminated clothing. Wash contaminated clothing before reusing. Ensure that eyewash stations and safety showers are close to the workstation location.

Personal protection

Respiratory

: Use a properly fitted, air-purifying or air-fed respirator complying with an approved standard if a risk assessment indicates this is necessary. Respirator selection must be based on known or anticipated exposure levels, the hazards of the product and the safe working limits of the selected respirator.

Hands

: Chemical-resistant, impervious gloves complying with an approved standard should be worn at all times when handling chemical products if a risk assessment indicates this is

Eyes

Skin

: Safety eyewear complying with an approved standard should be used when a risk assessment indicates this is necessary to avoid exposure to liquid splashes, mists, gases or dusts. If operating conditions cause high dust concentrations to be produced,

: Personal protective equipment for the body should be selected based on the task being performed and the risks involved and should be approved by a specialist before handling this product.

Environmental exposure

controls

: Emissions from ventilation or work process equipment should be checked to ensure they comply with the requirements of environmental protection legislation. In some cases, fume scrubbers, filters or engineering modifications to the process equipment will be necessary to reduce emissions to acceptable levels.

Physical and chemical properties

: Solid. [Crystalline powder.] **Physical state**

Color : White.

pН : 7.2 to 7.6 [Conc. (% w/w): 1%]

Solubility : Easily soluble in the following materials: cold water and hot water.

10. Stability and reactivity

Chemical stability

: The product is stable. Under normal conditions of storage and use, hazardous polymerization will not occur.

Conditions to avoid Incompatible materials : No specific data. : No specific data.

Hazardous decomposition products

: Under normal conditions of storage and use, hazardous decomposition products should

not be produced.

Possibility of hazardous

: Will not occur.

reactions

11/4/2009. 5/9 www.thermo.com/pierce

Thermo

Dulbecco's Phosphate Buffered Saline

11. Toxicological information

Acute toxicity

	Product/ingredient name	Result	Species	Dose	Exposu
	sodium chloride	LD50 Dermal	Rabbit	10000 mg/kg	-
		LD50 Intraperitoneal	Rat	2600 mg/kg	-
		LD50 Oral	Rat	3000 mg/kg	-
		LDLo Intraperitoneal	Rat	3.72 g/kg	-
		LDLo Subcutaneous	Rat	3500 mg/kg	-
		LC50 Inhalation Dusts and mists	Rat	>42 gm/m3	1 hours
	sodium phosphate, dibasic	LD50 Oral	Rat	17000 mg/kg	-
	, , , , , , , , , , , , , , , , , , , ,	LDLo Intraperitoneal	Rat	1 g/kg	-
		LDLo Intramuscular	Rat	1 g/kg	-
		LDLo Subcutaneous	Rat	1 g/kg	-
	potassium chloride	LD50 Intraperitoneal	Rat	660 mg/kg	-
		LD50 Intravenous	Rat	142 mg/kg	-
		LD50 Intravenous		142 mg/kg	-
		LD50 Oral	Rat - Male	2600 mg/kg	-
		TDLo Intracerebral	Rat	272.7 mg/kg	-
	phosphoric acid, monopotassium salt	LD50 Dermal	Rabbit	>4640 mg/kg	-
	F			1010	

LD50 Dermal LDLo Oral

Conclusion/Summary : Not available.

Chronic toxicity

Conclusion/Summary : Not available.

Carcinogenicity

Conclusion/Summary : Not available.

Classification

Product/ingredient name **OSHA ACGIH** IARC **EPA** NIOSH sodium chloride None. None. sodium phosphate, dibasic None. None. potassium chloride None. None. phosphoric acid, monopotassium salt None. None.

Rat

4640 mg/kg

6/9

Mutagenicity

Conclusion/Summary : Not available.

Teratogenicity

Conclusion/Summary : Not available.

Reproductive toxicity

Conclusion/Summary : Not available.

Europe

Chronic effects : Repeated or prolonged inhalation of dust may lead to chronic respiratory irritation.

11/4/2009. 3747 N. Meridian Rd. PO Box 117 Rockford, IL 61105



11. Toxicological information

Carcinogenicity : No known significant effects or critical hazards.

Mutagenicity : No known significant effects or critical hazards.

Teratogenicity : No known significant effects or critical hazards.

Developmental effects : No known significant effects or critical hazards.

Fertility effects : No known significant effects or critical hazards.

12. Ecological information

Environmental effects: No known significant effects or critical hazards.

United States

Aquatic ecotoxicity

Conclusion/Summary

: Not available.

Other adverse effects

: No known significant effects or critical hazards.

13. Disposal considerations

 Waste disposal
 : The generation of waste should be avoided or minimized wherever possible. Dispose of surplus and non-recyclable products via a licensed waste disposal contractor. Disposal of this product, solutions and any by-products should at all times comply with the

or this product, solutions and any by-products similar at miles comply with the requirements of environmental protection and waste disposal legislation and any regional local authority requirements. Avoid dispersal of spilled material and runoff and contact

with soil, waterways, drains and sewers.

Hazardous waste : Within the present knowledge of the supplier, this product is not regarded as

hazardous waste, as defined by EU Directive 91/689/EEC.

Disposal should be in accordance with applicable regional, national and local laws and regulations.

Refer to Section 7: HANDLING AND STORAGE and Section 8: EXPOSURE CONTROLS/PERSONAL PROTECTION for additional handling information and protection of employees.

14. Transport information

Regulatory information	UN number	Proper shipping name	Classes	PG*
DOT Classification	Not regulated.	-	-	-
IATA-DGR Class	Not regulated.	-	-	-

PG* : Packing group

15 . Regulatory information

United States

HCS Classification

: Irritating material Target organ effects

U.S. Federal regulations

: United States inventory (TSCA 8b): All components are listed or exempted.

Thermo

Dulbecco's Phosphate Buffered Saline

15 . Regulatory information

SARA 302/304/311/312 extremely hazardous substances: No products were found. SARA 302/304 emergency planning and notification: No products were found. SARA 302/304/311/312 hazardous chemicals: potassium chloride; sodium phosphate, dibasic; sodium chloride

SARA 311/312 MSDS distribution - chemical inventory - hazard identification: potassium chloride: Immediate (acute) health hazard, Delayed (chronic) health hazard; sodium phosphate, dibasic: Immediate (acute) health hazard; sodium chloride: Immediate (acute) health hazard, Delayed (chronic) health hazard

Clean Water Act (CWA) 307: No products were found.
Clean Water Act (CWA) 311: sodium phosphate, dibasic

Clean Air Act (CAA) 112 accidental release prevention: No products were found.
Clean Air Act (CAA) 112 regulated flammable substances: No products were found.
Clean Air Act (CAA) 112 regulated toxic substances: No products were found.

Canada

WHMIS (Canada)

: Class D-2B: Material causing other toxic effects (Toxic).

Canadian lists : CEPA Toxic substances: None of the components are listed.

Canadian ARET: None of the components are listed.

Canadian NPRI: None of the components are listed.

Alberta Designated Substances: None of the components are listed. Ontario Designated Substances: None of the components are listed. Quebec Designated Substances: None of the components are listed.

Canada inventory

EU regulations

: Canada inventory: All components are listed or exempted.

Risk phrases : This product is not classified according to EU legislation.

:

International regulations

International lists

: Australia inventory (AICS): All components are listed or exempted. China inventory (IECSC): All components are listed or exempted. Korea inventory (KECI): All components are listed or exempted. Philippines inventory (PICCS): All components are listed or exempted. Japan inventory (ENCS): All components are listed or exempted.

16 . Other information

Label requirements

: CAUSES RESPIRATORY TRACT, EYE AND SKIN IRRITATION. MAY BE HARMFUL IF SWALLOWED. CAN CAUSE TARGET ORGAN DAMAGE.

Hazardous Material Information System (U.S.A.)



The customer is responsible for determining the PPE code for this material.

National Fire Protection Association (U.S.A.)

11/4/2009.

3747 N. Meridian Rd. PO Box 117

Rockford, IL 61105 (815) 968-0747 (815) 968-7316 Fax ww.thermo.com/pierce

8/9

1/4/2009.

7/9

11/4/2009.

Biosciences Division

PO Box 117



16. Other information



 Date of printing
 : 11/4/2009.

 Date of issue
 : 11/4/2009.

 Date of previous issue
 : 2/13/2008.

 Version
 : 1.01

✓ Indicates information that has changed from previously issued version.

Full text of R-phrases : R36/37/38- Irritating to eyes, respiratory system and skin.

referred to in sections 2 and

3 - Europe

Full text of classifications : Xi - Irritant

referred to in sections 2 and

3 - Europe

Notice to reader

To the best of our knowledge, the information contained herein is accurate. However, neither the above-named supplier, nor any of its subsidiaries, assumes any liability whatsoever for the accuracy or completeness of the information contained herein.

Final determination of suitability of any material is the sole responsibility of the user. All materials may present unknown hazards and should be used with caution. Although certain hazards are described herein, we cannot guarantee that these are the only hazards that exist.

This product has been classified according to the hazard criteria of the CPR and the MSDS contains all the information required by the CPR.

11/4/2009. 9/9