



Interior

Front View



Side View


115VAC/60Hz/1PH
NEMA-5-15P

Cat# 76210-112/HCUCBI0404SSLH

Cu.Ft.	Defrost	Door	Int Door	Shelves	Interior Dimensions			Refrigerant	H.P.	Amps	Weight
					W"	D"	H"				
4.5	Cycle	1 Solid LH	0	2 Adj./1 Fixed	19 1/4	17 1/2	24 1/2	R600a	1/8	1.74	100 lbs.

Note: Interior dimensions are estimates and subject to change without notice