


Anti-OR1D4 + OR1D5 Antibody (A97066)

Specifications:

Name:	Anti-OR1D4 + OR1D5 Antibody
Description:	Rabbit polyclonal antibody to OR1D4 + OR1D5.
Applications:	WB, IF, ELISA
Recommended Dilutions:	WB: 1:500-1:1000, IF: 1:100-1:500, ELISA: 1:20000
Reactivity:	Human
Immunogen:	Synthetic peptide derived from human OR1D4/5 (amino acids 200-249).
Host:	Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugate:	Unconjugated
Purification:	Purified from rabbit serum by antigen affinity chromatography using the immunizing peptide.
Molecular Weight:	35kDa
Product Form:	Liquid
Formulation:	Supplied in Phosphate Buffered Saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, with 150mM NaCl, 0.02% Sodium Azide, and 50% Glycerol.
Storage:	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Disclaimer:	This product is for research use only. It is not intended for diagnostic or therapeutic use.

Images:




Western blot analysis of lysates from HeLa, HUVEC, Jurkat, HepG and MCF-7 cells using Anti-OR1D4 + OR1D5 Antibody. The right hand lane represents a negative control, where the antibody is blocked by the immunising peptide.

Anti-OR1D4 + OR1D5 Antibody (A97066)

Images continued:



Immunofluorescence analysis of MCF7 cells using Anti-OR1D4 + OR1D5 Antibody. The right hand panel represents a negative control, where the antibody was pre-incubated with the immunising peptide.



Western blot analysis of various cells using Anti-OR1D4 + OR1D5 Antibody.