



Certificate of Analysis

Apr 1, 2015 (JST)

TOKYO CHEMICAL INDUSTRY CO.,LTD.
4-10-1 Nihonbashi-Honcho, Chuo-ku, Tokyo 103-0023 Japan

Chemical Name: Dibutyl Succinate		
Product Number: S0102 CAS: 141-03-7	Lot: GI01	

Tests	Results	Specifications
Refractive index n ₂₀ /D	1.4292	1.4280 to 1.4310
Purity(GC)	99.8 %	min. 99.0 %
Specific gravity (20/20)	0.9780	0.9760 to 0.9800
Acid Value	0.1	max. 2.0

TCI Lot numbers are 4-5 characters in length.
Characters listed after the first 4-5 characters are control numbers for internal purpose only.

Customer service:

TCI AMERICA
Tel: +1-800-423-8616 / +1-503-283-1681
Fax: +1-888-520-1075 / +1-503-283-1987
E-mail: Sales-US@TCIchemicals.com