

FREEZE DRYING **FLASK w/SNAP-ON CAP**

Complete “fast-freeze” flask with snap-on rubber cap and plain or stopcock PTFE adapter. Flasks are designed for ease of handling and faster loading. Wide mouth provides easy access to inner walls. Snap-on black rubber cap is grease free; one size fits all capacity flasks. Plain straight adapter for connecting flask cap to lyophilizer manifold is removable.

Wide Mouth Series Flasks

- Compatible with all major brands of freeze dryers and manifolds
- Uniform wall weight for good temperature transfer
- Wide mouth provides ease in getting product in and out of the flask
- Black rubber cap matches O.D. and I.D. on the flask
- PTFE adapter tubes available plain or with stopcock
- Caps fit standard lyophilizer manifolds



Contact us for...

Custom glass and repair services.

FREEZE DRYING

FLASK w/TWISTLOR SEAL

These are just a few of the more widely used glassware components to replace your existing glass. All our rotary evaporator replacement glass is manufactured in the USA, to exacting specifications and standards.

Twistlor Series Flasks

- O-Ring for vacuum-tight seal
- Dual purpose: Rotary evaporator and freeze drying flask
- Adapter heads will fit all capacity flasks
- Adapter head comes with nylon retainers and FETFE O-Ring
- Note: These flasks are NOT for use in elevated temperatures



Also Available...

Contact us for all your freeze drying and rotavap glassware needs, for both large and bench scale applications.

avantor™
delivered by **vwr™**

VWR.COM

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada unless otherwise noted, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2022 Avantor, Inc. All rights reserved.