

# ADAPTER KITS

# ACE ASSORTMENTS™

Assortment kits ensure you'll have the right adapter or stopper on hand when you need it most. Choose from our best selling reducing and enlarging standard taper glass joint adapters, our ubiquitous Ace-Thred™ probe adapters and our most popular hollow glass stoppers.

## Available Kits:

- Standard taper enlarging and reducing adapters
- Standard taper reducing adapters
- Ace-Thred™ (#7) to standard taper adapters
- Ace-Thred™ (#11 & #15) to standard taper adapters
- Glass stoppers
- Replacement Ferrules and O-Rings for adapter kits



## O-Ring Kits

- Various popular sizes
- Various materials available



Contact us for custom applications.

# ACE-THRED™ GLASS ADAPTERS

Ace-Thred™ adapters come with a variety of joints at the bottom and a threaded nylon bushing at the top. The bushing tightens into an Ace-Thred™ to form an O-Ring compression seal with electrodes, thermometers, thermocouples, bleed tubes, etc. The nylon bushing comes complete with FETFE O-Ring. These adapters are suitable for vacuum work.

## Ace-Thred™ Adapter Advantages:

- Grease-free connections
- Clamp-free
- Vacuum-tight seal
- Convert virtually any connection to Ace-Thred™
- Reduce or enlarge your connection
- Air sensitive applications

## Ace-Thred™ Adapter Features:

- Available in various sizes: #7, #11, #15, #18, #25, #36 & #50
- Transition to standard taper joint, spherical joint, beaded pipe or NW flange connections
- Works with nylon or PTFE bushing
- Works with various O-Ring materials
- Custom applications available

*For complete listing of glass adapters go to [aceglass.com](http://aceglass.com)*

