

Quality Certificate

VWR® SINGLE SCALE GRADUATED CYLINDERS, POLYPROPYLENE

As per the manufacturer, the below product meets the following criteria:

North American Catalog No:	65000-002
Lot Number:	J335456
Description:	VWR CYLINDER POLYPRO 25ML
Date of Manufacture:	September 16, 2019
Country of Origin:	United States of America

Quality Conformance

The product is in conformance with the Purchase Order Contract, drawings, specifications, and other applicable documentation.

Product Specifications

Material: Polypropylene

Dimensions: 25 ml cylinder has an overall height of 7.48", Base dia. @ max 1.96", Cylinder dia. .728" OD and .61" ID. Scale is in 0.5ml subdivisions. Tolerance: +/- 0.020"

Other Compliances:

- Graduated cylinders have a total contained, total delivered scale on one side with raised numerals.
- Large spout facilitates pouring and hexagon base provides stability.
- Non-wetting surface eliminates concave meniscus and gives more accurate readings.
- Cylinders are chemical resistant and steam autoclavable at 121°C +/- 2 ml / 5% accuracy. Class B

Signed:



Jamie Ethier
VP Global Quality
VWR, Part of Avantor

Date: October 12, 2020