

## Quality Certificate

### VWR® ADVANCED PROTECTION CLASS 100 NITRILE GLOVES

As per the manufacturer, the below product meets the following criteria:

<b>North American Catalog No:</b>	89130-836
<b>Lot No:</b>	1804N2154
<b>Description:</b>	VWR GLOVES ADVP NITRL CS1000 SM
<b>Manufacturing Date:</b>	April-2018
<b>Product Shelf Life :</b>	5 years from the Manufacturing Date
<b>Country of Origin :</b>	Malaysia

**Material:** Acrylonitrile-Butadiene Copolymer Latex

**QC Testing:** Representative production samples are collected and inspected in accordance with current applicable Product specification.

**Accuracy:**

Water tight test:	Inspection level = G1;	AQL 1.5
Visual inspection:	Inspection level = G1;	AQL 1.5 (critical defects)
		AQL 2.5 (major defects)
		AQL 4.0 (minor defects)
Dimensional:	Inspection level = S2;	AQL 4.0
Physical Property:	Inspection level = S2;	AQL 4.0
Glove count:	Inspection level = 5 bags;	Tolerance 100+/-2 pieces

**Sterility** N/A

**Non-Pyrogenic:** N/A

**RNase/DNase Free:** N/A

**Regulatory:** ISO 9001, ISO 14001

**Signed:**



Linda K. Shine  
Manager - Supplier Qualification  
Global Compliance  
VWR International LLC

**Date:** April 26, 2018