

Quality Certificate

VWR® SPINBAR® STIR BARS, CYLINDRICAL

As per the manufacturer, the below product meets the following criteria:

North American Catalog No:	58949-108
Lot Number:	J256515
Description:	VWR Stir Bar Ptfе 5/16x1-1/2in.
Date of Manufacture:	March 26, 2020
Product Shelf Life:	Not Applicable
Country of Origin:	United States of America

Quality Conformance

The product is in conformance with the Purchase Order Contract, drawings, specifications, and other applicable documentation.

Product Specifications

- Material:** Dupont Virgin White Teflon #807N X; Magnet, Alnico 5
- Dimensions:** 1 ½" x 5/16"
- Other Compliances:**
- Autoclavable at 121 °C
 - Materials Meet FDA Regulations 21CFR177.1550, 21CFR177.1520
 - Teflon used to completely encapsulate the magnetic bar meets US Pharmacopia Class VI
 - Based on the information provided the product does not contain or come in contact with lead, latex, Mercury, Formaldehyde or other contaminants.
 - The product does not contain or come in contact with any animal origin.
 - All Items are carefully inspected to conform with ISO 9001:2008
 - There is no expiration for this product.

Signed:



Jamie Ethier
VP Global Quality
VWR, Part of Avantor

Date: October 15, 2020