

Culture medium	Protein yields using FectoPRO®
ExpiCHO™ (Gibco™)	+++
FreeStyle™ CHO (Gibco™)	+++
ProCHO™ 4 (Lonza®)	+++
CD-CHO™ (Gibco™)	+++
FreeStyle™ F17 (Gibco™)	+++
BalanCD® Transfectory™ CHO (Irvine Scientific™)	+++
CD-FortiCHO™ (Gibco™)	+++
CD-OptiCHO™ (Gibco™)	-
HyClone™ CDM4 CHO (GE Healthcare™)	-
PowerCHO™ 2 (Lonza®)	-

FectoPRO® is compatible with different synthetic serum free media formulations for CHO cells.