

## Quality Certificate

### VWR® NITRILE EXTENDED CUFF EXAM GLOVES

As per the manufacturer, the below product meets the following criteria:

<b>North American Catalog Number:</b>	89521-902
<b>Lot Number:</b>	GT19K313
<b>Description:</b>	VWR Glove Nitrile EC PF Exam XL PK50
<b>Manufacturing Date:</b>	August 2019
<b>Product Shelf Life:</b>	5 years from the Manufacturing Date
<b>Country of Origin:</b>	Malaysia

**Material:** Nitrile Butadiene Rubber, 100% Latex Free

**QC Testing:** Representative production samples are collected and inspected in accordance with current applicable Product specification.

<b>Accuracy:</b>	Water tight test:	Inspection level = G1;	AQL 0.65
	Visual inspection:	Inspection level = G1;	AQL 0 (critical defects)
			AQL 1.5 (major defects)
			AQL 2.5 (minor defects)
	Dimensional:	Inspection level = S2;	AQL 4.0
	Physical Property:	Inspection level = S2;	AQL 4.0
	Powder Content:	Inspection level = 5 pieces;	Average ≤ 2 mg/glove
	Moisture Content:	Inspection level = 5 pieces;	Average ≤ 1%
	Glove count:	Inspection level = 5 dispensers;	Tolerance 50+/-1 pieces

**Sterility** N/A

**Non-Pyrogenic:** N/A

**RNase/DNase Free:** N/A

**Regulatory:** FDA 510K, ISO 9001, ISO 13485, GMP Certificate

**Signed:**



Jamie Ethier  
VP Global Quality  
VWR, Part of Avantor

**Date:** January 11, 2020