

VWR® disposable glove guide

01. CLEANROOM GLOVES

02. EXAM OR SAFETY
GLOVES



VWR disposable glove guide

Visit
**VWR.COM/
VWRGLOVES**
to search
online!

The VWR Collection line of disposable gloves offers a variety of choices depending on your application. Whether your needs are production, laboratory, or industrial focused, we have you covered. Use the guide below to understand the features of our VWR® Disposable Gloves and choose the right level of protection from various material, thickness, and length options.



CLEANROOM GLOVES

	Product Line	Packed ISO Class	Static Dissipative	Finger Thickness	Length	Color	Cuff	Texture	Hand Specific/ Ambidextrous	Packaging	Cat No. Series
Latex											
	VWR® CERTICLEAN® Class 10 Latex Gloves	ISO Class 4 (FED-STD-209E Class 10/M2.5)	No	8.7	12 in	Natural	Beaded	Full	Ambidextrous	1000/cs	40101-146
	VWR® CERTICLEAN® Class 100 Latex Gloves	ISO Class 5 (FED-STD-209E Class 100/M3.5)	No	8.7	12 in	Natural	Beaded	Full	Ambidextrous	1000/cs	40101-158
	VWR Cleanroom Latex, Non-Sterile, Class 10	ISO Class 4 (FED-STD-209E Class 10/M2.5)	No	7.3	12 in	White	Beaded	Full	Ambidextrous	1000/cs	76579-674
	VWR Cleanroom Latex, Sterile, Class 10	ISO Class 4 (FED-STD-209E Class 10/M2.5)	No	10.0	12 in	White	Beaded	Fingers	Hand Specific	200pr/cs	76533-186

	Product Line	Packed ISO Class	Static Dissipative	Finger Thickness	Length	Color	Cuff	Texture	Hand Specific/ Ambidextrous	Packaging	Cat No. Series
Nitrile											
	VWR® CERTICLEAN® Class 10 Nitrile Gloves	ISO Class 4 (FED-STD-209E Class 10/M2.5)	Yes	6.3	12 in	White	Beaded	Fingertips	Ambidextrous	1000/cs	40101-352
	VWR® Maximum Protection Class 10 Nitrile Gloves	ISO Class 4 (FED-STD-209E Class 10/M2.5)	Yes	5.5	12 in	White	Beaded	Fingers	Ambidextrous	1000/cs	89130-846
	VWR® Advanced Protection Class 100 Nitrile Gloves	ISO Class 5 (FED-STD-209E Class 100/M3.5)	Yes	5.5	12 in	White	Beaded	Fingers	Ambidextrous	1000/cs	89130-836
	VWR Cleanroom Nitrile Non-Sterile Class 10 Non-Tacky	ISO Class 4 (FED-STD-209E Class 10/M2.5)	Yes	6.0	12 in	White	Beaded	Fingers	Ambidextrous	1000/cs	76533-198
	VWR Cleanroom Nitrile Non-Sterile Class 10 Tacky	ISO Class 4 (FED-STD-209E Class 10/M2.5)	Yes	6.0	12 in	White	Beaded	Fingers	Ambidextrous	1000/cs	76533-210
	VWR Cleanroom Nitrile, Sterile, Class 10	ISO Class 4 (FED-STD-209E Class 10/M2.5)	Yes	6.5	12 in	White	Beaded	Fingers	Hand Specific	200pr/cs	76533-224
	VWR Sterile Cleanroom Nitrile Glove, Class 100	ISO Class 5 (FED-STD-209E Class 100/M3.5)	Yes	5.5	12 in	White	Beaded	Fingers	Ambidextrous	1000/cs	76423-448
	VWR Sterile Class 100 HS Nitrile Gloves	ISO Class 5 (FED-STD-209E Class 100/M3.5)	Yes	5.7	12 in	White	Beaded	Fingers	Hand Specific	200pr/cs	76458-120
	VWR Cleanroom Nitrile, Non-Sterile, ISO 4, 400mm (16")	ISO Class 4 (FED-STD-209E Class 10/M2.5)	Yes	5.3	16 in	White	Beaded	None	Ambidextrous	1000/cs	76557-792
	VWR Cleanroom Nitrile, Sterile, ISO 4, 400mm (16")	ISO Class 4 (FED-STD-209E Class 10/M2.5)	Yes	5.3	16 in	White	Beaded	None	Ambidextrous	200pr/cs	76557-804

EXAM GLOVES

	Product Line	Feature	Material	Chemo Rated	Finger Thickness	Length	Texture	Cuff	Color	Packaging	Cat No. Series
	VWR® Blue Nitrile Examination Gloves	Everyday Glove, Ideal for Many Applications	Nitrile	Yes	4.7	9 in	Finger	Beaded	Blue	100pk/1000cs	76518-338
	VWR® Premium Thickness Nitrile Examination Gloves	Increased Thickness for Enhanced Protection	Nitrile	Yes	6.0	9.5 in	Full	Beaded	Blue	100pk/1000cs	76582-340
	VWR® Black Nitrile Examination Gloves	Exceptional Barrier Protection, Sleek Look, Fentanyl Approved	Nitrile	Yes	5.5	9.5 in	Finger	Beaded	Black	100pk/1000cs	76307-464
	VWR® Ultra Light Black Nitrile Glove	Space-Saving and Waste-Reducing Packaging	Nitrile	Yes	2.5	9.5 in	Finger	Beaded	Black	250pk/2500cs	76457-132
	VWR® Nitrile Extended Cuff Examination Gloves	Extended Cuff High Risk Protection	Nitrile	Yes	6.0	12 in	Finger	Beaded	Blue	100pk/1000cs	76582-350
	VWR Vinyl Nitrile Hybrid Exam Glove	Benefits of nitrile and durability of Vinyl	Vinyl/ Nitrile hybrid	No	4.0	9.5 in	None	Beaded	Blue	100pk/1000cs	76456-894
	VWR® Powder-Free Latex Examination Gloves	Natural Rubber Latex For Everyday Use, High Barrier Protection	Latex	Yes	5.5	9.5 in	Finger	Beaded	Natural	100pk/1000cs	76319-668
	VWR Vinyl Examination Glove	Durable Basic and Food Services Gloves	Vinyl	No	4.5	9.5 in	None	Rolled	Transparent	100pk/1000cs	76459-214