

ONE / TWO-PIECE PRESSURE REACTORS

Designed to perform low to moderate positive pressure reactions, synthesis and catalysis. Jacketed for heating and cooling. Reactor capacities range from 500mL to 5000mL and use Ace-Threds to achieve a leak-tight system. One-piece and two-piece reactor systems available. *Rated up to 45psig at 100°C.*

Two-Piece Base Systems include:

- Jacketed Reactor Body
- Head and Clamp
- Bearing, Stir Shaft and Agitator
- Pressure Relief Manifold
- Addition Funnel, Condenser & Thermowell
- Stir Motor and Controller

Bottom Outlet System replaces:

- Jacketed Reactor Body w/bottom outlet

One-Piece Reactor Body

- Base systems similar to above except in a one-piece flask design
- With or without bottom outlet

JACKETED



Contact us for custom or large scale applications.

ONE / TWO-PIECE PRESSURE REACTORS

Designed to perform low to moderate positive pressure reactions, synthesis and catalysis. Reactor capacities range from 500mL to 5000mL and use Ace-Threads to achieve a leak-tight system. One-piece and two-piece reactor systems available. *Rated up to 45psig at 100°C.*

Two-Piece Base Systems include:

- Unjacketed Reactor Body
- Head and Clamp
- Bearing, Stir Shaft and Agitator
- Pressure Relief Manifold
- Addition Funnel, Condenser & Thermowell
- Stir Motor and Controller
- Heating Mantle and Temperature Controller

Bottom Outlet System replaces:

- Unjacketed Reactor Body w/bottom outlet

One-Piece Reactor Body

- Base systems similar to above except in a one-piece flask design
- With or without bottom outlet



UNJACKETED



Shown with optional Stand and Support Shelf.



1.800.932.5000 | VWR.COM

Prices and product details are current when published and subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2019 Avantor, Inc. All rights reserved.