

| | |
|---------------|--|
| Product Name: | PF051SSS/0M VWR #76271-310 |
| Description: | <p>White Diamond Series Undercounter Freezer Built-In</p> <p>4.2 Cu. Ft. Capacity</p> <p>One stainless steel swing door, right hinged*</p> <p>Microprocessor temperature controller</p> <p>Digital temperature display</p> <p>Audible and visual high and low temperature alarms</p> <p>2 probes; 1 in air and 1 in sample bottle</p> <p>Remote alarm contacts</p> <p>Adjustable operating temperature range: -15°C – -25°C</p> <p>Two fixed shelves</p> <p>Keyed door lock</p> <p>Recessed door handle</p> <p>Magnetic door gasket for positive seal</p> <p>Leveling legs</p> <p>Probe access port (3/8")</p> <p>Stainless steel exterior</p> <p>High density urethane foam cabinet and door insulation</p> <p>HFC-free refrigerant (R600a), foam insulation and packaging</p> <p>Three year parts and labor warranty, plus an additional two year compressor parts only warranty</p> <p>Pyxis®, Omnicell® and AcuDose RX® compatible</p> <p>Exterior dimensions: 23 ¾ "W x 24 ½" D x 33 3/8 "H</p> <p>Manual defrost</p> <p>115V, 60 Hz, 1.5 Amps, 1/8 HP</p> <p>UL/C-UL listed</p> <p>Shipping weight: 131 lbs.</p> <p>*This product is available left hinged, product number PF051SSSLH/0M</p> |