

Certificate of Quality

VWR Catalogue Number	See List
Description	VWR® 18 mm Screw-Thread Headspace Vials and Caps
Country of Origin	Manufactured in Germany
Date of Issue (yyyy-mm-dd)	2021-06-03

Vials

VWR Catalogue Number	Capacity [mL]	Colour
10119-906	10	Clear
10119-908	10	Amber
10119-902	20	Clear
10119-904	20	Amber

Caps

VWR Catalogue Number	Septa thickness [mm]	Septa Type
10119-914	1.3	White Silicone/Red PTFE
10119-916	1.3	Blue Silicone/White PTFE
10119-910	1.5	White Silicone/Blue PTFE
10119-918	1.5	White Silicone/Red PTFE, Star Slit
10119-912	1.6	Red Butyl/Grey PTFE

Quality System Compliance

Products are manufactured under the **ISO 9001** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

Product Specifications

Material	Type 1A 33 expansion Hydrolytic glass.
Sterilization	Not applicable.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	The product is Latex free.
BPA Statement	Not applicable.
DEHP Statement	Not applicable.
RoHS	Not applicable.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 01/19/2021 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.
Storage Conditions	Store at room temperature.
Shelf Life	2 years.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.