



Quality Certificate

VWR* ASPIRATING PIPETS

This document certifies that these pipets meet the following criteria:

Product Number	414004-266
Lot Number	110413
Description	VWR PIPET ASPIRATING 5ML CS200

- Material:** 100% Virgin High Clarity Polystyrene
- QC Testing:** Representative production samples are collected and inspected in accordance with current applicable product specifications.
- Accuracy:** Pipets are accurate $\pm 3\%$ at full volume.
- Sterility:** Product labelled as sterile is gamma irradiated and dorsimetrically released upon US Association for the Advancement of Medical Instrumentation (AAMI) recommended practices in effect at the time of validation. Products labelled sterile meets a minimum requirement of 10^{-6} SAL.
- Non-Pyrogenic:** The acceptance level of the product is 0.5 EU/ml or less than 20EU/device. LAL Gel Clot Method.
- Date of Sterilization:** June 12, 2011
- Shelf-life:** 4 years from sterilization date
- BSE TSE Statement:** This document certifies that the above product is not derived from and does not contain any source containing BSE/TSE. It is of chemical synthesis origin.

Signed:

A handwritten signature in blue ink, appearing to read 'Ken Crossley'.

Ken Crossley
Manager
Quality Assurance

Date of Issue: 22 June 2011