

# VWR® Absorbents

**vwr**™  
part of avantor

01. ABSORB LEAKS AND SPILLS

02. HAZMAT, GENERAL-PURPOSE &  
AND OIL-ONLY SPILL CONTROL

03. CHOOSE FROM PADS, ROLLS,  
PILLOWS AND SOCKS



# VWR® Hazmat Absorbents

## IDEAL FOR SPILL CLEANUP

High-capacity absorbents stay on the job longer and are perfect for stocking spill response kits.

## HANDLE HAZARDOUS CHEMICALS

Our Hazmat absorbents are designed for aggressive fluids such as acids, caustics, bases and unknown liquids, which can cause regular cellulose mats, socks, and pillows to degrade or react.

## EASY TO PICK OUT AT A GLANCE

Yellow color is easily distinguishable to ensure workers use the correct mat during a spill emergency.

## THE RIGHT ABSORBENT FOR EVERY JOB

Choose from pads for lining benches, rolls to cover large areas, pillows for high-volume spills, and socks to contain spills and keep them from spreading. Our Hazmat spill response kit includes absorbents plus PPE so you have everything you need to tackle chemical spills.

**Make this VWR absorbents collection the solution to your everyday needs. Contact your VWR Sales Representative or call 1 800 932 5000 today for more information!**



Specially treated to absorb highly corrosive liquids.



82028-024

### VWR® HAZMAT PADS

Great for lining benches and fume hoods, our fast-wicking dimpled design absorbs liquids quickly to leave area clean and dry for improved safety.

Description	Size, W x L, cm	Cat. No.	Unit
Hazmat Pad	38.1 x 48.3	<b>82028-024</b>	Pack of 100



82028-026

### VWR® HAZMAT ROLLS

Rolls are ideal for covering large areas and soaking up leaks and overspray. Dimple pattern speeds wicking across surface for fast absorption.

Description	Size, W x L, cm	Cat. No.	Unit
Hazmat Rolls	76.2 x 45.7	<b>82028-026</b>	Each



89169-002

### VWR® HAZMAT PILLOWS

Large surface area and sorbent filler handles large spills quickly and efficiently. Easily pulls liquids from the floor, catches leaks, and absorbs spills.

Description	Size, W x L x D, cm	Cat. No.	Unit
Hazmat Pillows	5.7 x 20.3 x 5.1	<b>89169-002</b>	Case of 20



89169-000

### VWR® HAZMAT SOCKS

Mold easily for hugging corners and walls. Made of 100% chemical-resistant skin and fast-wicking filler. Contain and absorb to keep fluid from spreading to larger areas.

Description	Size, D x L, cm	Cat. No.	Unit
Hazmat Socks	7.6 x 122	<b>89169-000</b>	Case of 30




89168-998

### VWR® HAZMAT LAB SPILL RESPONSE KIT

Includes 20 pads, three socks, one pair nitrile gloves, one pair safety goggles, and one water-resistant yellow bag with tie. Supports compliance with 40 CFR 263.30(a) safety standard.

Description	Absorbency	Cat. No.	Unit
Hazmat Lab Spill Response Kit	Up to 34.1 L	<b>89168-998</b>	Each





## VWR® GENERAL-PURPOSE PADS AND ROLLS

General-purpose pads and rolls absorb water-based liquids, petroleum products, and solvents for everyday maintenance and cleanup. Soak up leaks and spills fast, keeping floors at your workplace clean, dry, and safe.

Description	Size, W x L, cm	Cat. No.	Unit
General-Purpose Pad	38.1 x 48.3	<b>82028-028</b>	Pack of 100
General-Purpose Roll	76.2 x 45.7	<b>82028-018</b>	Each



## VWR® OIL-ONLY PADS AND ROLLS

Absorbs oil and fuels, but not water. Ideal for catching drips and soaking up spills. Clearly shows saturation level during spill response.

Description	Size, W x L, cm	Cat. No.	Unit
Oil-Only Pads	38.1 x 48.3	<b>82028-020</b>	Pack of 100
Oil-Only Rolls	76.2 x 45.7	<b>82028-022</b>	Each



800 932 5000 | VWR.COM

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit [vwr.com](http://vwr.com) to view our privacy policy, trademark owners, and additional disclaimers. © 2019 Avantor, Inc. All rights reserved.