

The following list contains the Material Safety Data Sheets you requested. Please scroll down to view the requested MSDS(s).

Product	MSDS	Distributor	Format	Language	Quantity
8508600	2756549	Hach Company	ROWGHS	English	1
8508600	2283449	Hach Company	ROWGHS	English	1
8508600	2283549	Hach Company	ROWGHS	English	1
8508600	2283649	Hach Company	ROWGHS	English	1

Total Enclosures: 4

World Headquarters
Hach Company
P.O.Box 389
Loveland, CO USA 80539
(970) 669-3050

MSDS No: M01702

SAFETY DATA SHEET

1. CHEMICAL PRODUCT AND COMPANY IDENTIFICATION

Product Name: pH Storage Solution

Catalog Number: 2756549

Hach Company
P.O.Box 389
Loveland, CO USA 80539
(970) 669-3050

Emergency Telephone Numbers:
(Medical and Transportation)
(303) 623-5716 24 Hour Service
(515)232-2533 8am - 4pm CST

MSDS Number: M01702

Chemical Name: Not applicable

CAS Number: Not applicable

Additional CAS No. (for hydrated forms): Not applicable

Chemical Formula: Not applicable

Chemical Family: Not applicable

Intended Use: Laboratory Reagent Electrode storage solution

2. HAZARDS IDENTIFICATION

This mixture is not classified as hazardous per GHS (UN publication ST/SG/AC.10/36/Add.3)

GHS Classification:

Hazard categories: . Serious Eye Damage/Eye Irritation: Eye Irrit. 2A

GHS Label Elements:

WARNING



Hazard statements: . Causes serious eye irritation.

Precautionary statements: Wear eye protection. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention. If skin irritation occurs: Get medical advice/attention.

HMIS:

Health: 1

Flammability: 0

Reactivity: 0

Protective Equipment: X - See protective equipment, Section 8.

NFPA:

Health: 1

Flammability: 0

Reactivity: 0

Symbol: Not applicable

WHMIS Hazard Classification: Class D, Division 2, Subdivision B - Toxic material (other toxic effects)

WHMIS Symbols: Other Toxic Effects

3. COMPOSITION / INFORMATION ON INGREDIENTS

Hazardous Components according to GHS:

Potassium Chloride

CAS Number: 7447-40-7

Chemical Formula: KCl

GHS Classification: Acute Tox. 5 -Orl, H303; Skin Irrit. 3, H316; Eye Irrit. 2, H319; Aquatic Acute 3, H402

Percent Range (Trade Secret): 15.0 - 25.0

Percent Range Units: weight / weight

PEL: 15 mg/m³ as inhalable dust; 5 mg/m³ as respirable dust

TLV: 10 mg/m³ as inhalable dust; 3 mg/m³ as respirable dust

WHMIS Symbols: Not applicable

Potassium Phosphate, Monobasic

CAS Number: 7778-77-0

Chemical Formula: KH₂PO₄

GHS Classification: Acute Tox. 4 - Orl, H302; Eye Irrit. 2, H319

Percent Range (Trade Secret): < 0.5

Percent Range Units: weight / weight

PEL: 15 mg/m³ as total dust; 5 mg/m³ as respirable dust

TLV: 10 mg/m³ as inhalable fraction; 3 mg/m³ as respirable fraction

WHMIS Symbols: Not applicable

Glutaric dialdehyde

CAS Number: 111-30-8

Chemical Formula: C₅H₈O₂

GHS Classification: Acute Tox. 3 - Inh, H331; Acute Tox. 3 - Orl, H301; Skin Corr. 1B, H314; Resp. Sens. 1, H334; Skin Sens. 1, H317; Aquatic Acute 1, H400

Percent Range (Trade Secret): < 0.1

Percent Range Units: weight / weight

PEL: Not established

TLV: Not established

WHMIS Symbols: Acute Poison

Sodium Phosphate, Dibasic

CAS Number: 7558-79-4

Chemical Formula: Na₂HPO₄

GHS Classification: Skin Irrit. 2, H315; Eye Irrit. 2A, H319

Percent Range (Trade Secret): < 0.1

Percent Range Units: weight / weight

PEL: 15 mg/m³ as inhalable dust; 5 mg/m³ as respirable dust

TLV: 10 mg/m³ as inhalable dust; 3 mg/m³ as respirable dust

WHMIS Symbols: Not applicable

Hazardous Components according to GHS: No

Demineralized Water

CAS Number: 7732-18-5

Chemical Formula: H₂O

GHS Classification: Not a dangerous substance according to GHS.

Percent Range (Trade Secret): 75.0 - 85.0

Percent Range Units: weight / weight

PEL: Not established

TLV: Not established

WHMIS Symbols: Not applicable

4. FIRST AID MEASURES

General Information: In the event of exposure, show this Material Safety Data Sheet and label (where possible) to a doctor.

Advice to doctor: Treat symptomatically.

Eye Contact: Flush eyes with water. Call physician if irritation develops.

Skin Contact (First Aid): Wash skin with plenty of water. Call physician if irritation develops. Remove contaminated clothing.

Inhalation: Remove to fresh air. Give artificial respiration if necessary. If you feel unwell, contact a physician

Ingestion (First Aid): Never give anything by mouth to an unconscious person. Rinse mouth with plenty of water. Give 1-2 glasses of water. If you feel unwell, contact a physician.

5. FIRE FIGHTING MEASURES

Flammable Properties: Material will not burn. Material is not classified as flammable according to GHS criteria. During a fire, irritating and highly toxic gases may be generated by thermal decomposition.

Fire Fighting Instruction: As in any fire, wear self-contained breathing apparatus pressure-demand and full protective gear. Evacuate area and fight fire from a safe distance.

Extinguishing Media: Use media appropriate to surrounding fire conditions

Extinguishing Media NOT To Be Used: Not applicable

Fire / Explosion Hazards: This product will not burn or explode. May react violently with: strong acids strong bases strong oxidizers

Hazardous Combustion Products: This material will not burn.

6. ACCIDENTAL RELEASE MEASURES

Spill Response Notice:

Only persons properly qualified to respond to an emergency involving hazardous substances may respond to a spill according to federal regulations (OSHA 29 CFR 1910.120(a)(v)) and per your company's emergency response plan and guidelines/procedures. See Section 13, Special Instructions for disposal assistance. Outside of the US, only persons properly qualified according to state or local regulations should respond to a spill involving chemicals.

Containment Technique: Stop spilled material from being released to the environment. Absorb spilled liquid with non-reactive sorbent material. Releases of this material may contaminate the environment.

Clean-up Technique: If permitted by regulation, Dilute with a large excess of water. Flush the spilled material to the drain with a large excess of water. Otherwise, Dispose of in accordance with local, state and federal regulations or laws. Decontaminate the area of the spill with a soap solution.

Evacuation Procedure: Evacuate as needed to perform spill clean-up. If conditions warrant, increase the size of the evacuation.

DOT Emergency Response Guide Number: Not applicable

7. HANDLING AND STORAGE

Handling: Avoid contact with eyes skin clothing Wash thoroughly after handling. Use with adequate ventilation. Do not breathe mist or vapors. Maintain general industrial hygiene practices when using this product.

Storage: Keep container tightly closed when not in use. Keep this product in its original container when not in use. Store in a cool, dry, well-ventilated place. Protect from: heat extreme temperatures Keep away from: acids bases oxidizers

Flammability Class: Not applicable

8. EXPOSURE CONTROLS / PERSONAL PROTECTION

Engineering Controls: Maintain general industrial hygiene practices when using this product.

Personal Protective Equipment:

Eye Protection: safety glasses with top and side shields

Skin Protection: disposable latex gloves In the EU, the selected gloves must satisfy the specifications of EU Directive 89/686/EEC and standard EN 374 derived from it.

Inhalation Protection: adequate ventilation

Precautionary Measures: Avoid contact with: eyes skin clothing Wash thoroughly after handling. Use with adequate ventilation. Do not breathe: mist/vapor Protect from: heat Keep away from: acids/acid fumes bases oxidizers

TLV: Not established

PEL: Not established

For Occupational Exposure Limits (OEL) for ingredients, see section 3 - Composition/Information on Ingredients.:

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance: Clear, colorless liquid

Physical State: Liquid
Molecular Weight: Not applicable
Odor: Odorless
Odor Threshold: Not applicable
pH: 6.4
Metal Corrosivity:
Corrosivity Classification: Not classified as corrosive to metals according to GHS criteria.
Steel: Not determined
Aluminum: Not determined
Specific Gravity/ Relative Density (water = 1; air = 1): Estimation: 1.15
Viscosity: Not determined
Solubility:
Water: Miscible
Acid: Miscible
Other: Not determined
Partition Coefficient (n-octanol / water): Not applicable
Coefficient of Water / Oil: Not applicable
Melting Point: Estimation: -49 °C (-56 °F)
Decomposition Temperature: Not determined
Boiling Point: Estimation: 113 °C (236 °F)
Vapor Pressure: Estimation: 16.5 mm Hg (2.2 kPa) at 20 °C (68 °F)
Vapor Density (air = 1): Estimation: 0.62
Evaporation Rate (water = 1): Estimation: 0.87
Volatile Organic Compounds Content: Not applicable
Flammable Properties: Material will not burn. Material is not classified as flammable according to GHS criteria. During a fire, irritating and highly toxic gases may be generated by thermal decomposition.
Flash Point: Not applicable
Method: Not applicable
Flammability Limits:
Lower Explosion Limits: Not applicable
Upper Explosion Limits: Not applicable
Autoignition Temperature: Not applicable
Explosive Properties:
Not classified according to GHS criteria.
Oxidizing Properties:
Not classified according to GHS criteria.
Reactivity Properties:
Not classified as self-reactive, pyrophoric, self-heating or emitting flammable gases in contact with water according to GHS criteria.
Gas under Pressure:
Not classified according to GHS criteria.

10. STABILITY AND REACTIVITY

Chemical Stability: Stable when stored under proper conditions.
Mechanical Impact: None reported
Static Discharge: None reported.
Reactivity / Incompatibility: May react violently in contact with: bromine trifluoride acids alkalies oxidizers
Hazardous Decomposition: Heating to decomposition releases toxic and/or corrosive fumes of: chlorides potassium oxide
Conditions to Avoid: Contact with acid or acid fumes Contact with oxidizers Excessive heat Extreme temperatures Incompatibles

11. TOXICOLOGICAL INFORMATION

Toxicokinetics, Metabolism and Distribution: No information available for mixture.
Toxicologically Synergistic Products: None reported
Acute Toxicity: Acute Toxicity Estimate (ATE) - Calculated from Ingredient Toxicity Data Practically Non-toxic Based on classification principles, the classification criteria are not met.
Specific Target Organ Toxicity - Single Exposure (STOT-SE): Based on classification principles, the classification criteria are not met. Summary of findings reported in the literature follow.

Potassium Chloride: Oral Man LDLo = 20 mg/kg/Cardiac arrhythmias, nausea, vomiting and change in blood clotting factors; Oral Man TDLo = 214.29 mg/kg/Hypermotility, diarrhea, nausea and vomiting.

Specific Target Organ Toxicity - Repeat Exposure (STOT-RE): Based on classification principles, the classification criteria are not met. Summary of findings reported in the literature follow.

Potassium Chloride: Oral Rat TDLo = 75.6 g/kg/6 wk/Increased urine volume; Oral Rat TDLo = 938 g/kg/78 wk/Acute renal failure; Oral Rat TDLo = 1536 g/kg/130 wk/Adrenal cortex hyperplasia.

Skin Corrosion/Irritation: Mildly irritating to skin.

Eye Damage: Irritating to eyes.

Sensitization: Based on classification principles, the classification criteria are not met. Contains a sensitizing compound. Glutaric Dialdehyde

CMR Effects/Properties (carcinogenic, mutagenic or toxic to reproduction): Based on classification principles, the classification criteria are not met. Summary of findings reported in the literature follow.

Potassium Chloride: Oral Rat Unscheduled DNA Synthesis - 1.5 mg/kg; Mouse Lymphocyte Mutation in Microorganism - 2048 mg/L; Human Leukocyte DNA Damage - 1 mmol/L/2 hr; Hamster Ovary Sister Chromatid Exchange - 180 mmol/L Glutaric Dialdehyde: Sister chromatid exchange - Hamster: ovary 110 µg/l; Teratogenic effects; Oral rat TDLo 4370 mg/kg/Reproductive effects; There are grounds for suspecting carcinogenic potential, of Glutaric Dialdehyde however, the data is not conclusive

This product does NOT contain any IARC listed chemicals.

This product does NOT contain any NTP listed chemicals.

This product does NOT contain any OSHA listed carcinogens.

Symptoms/Effects:

Ingestion: Very large doses may cause: anorexia blood pressure problems cardiac depression fever gastroenteritis

Inhalation: Very large doses may cause: allergic respiratory reaction

Skin Absorption: No effects anticipated

Chronic Effects: Chronic overexposure may cause allergic skin reactions allergic respiratory reactions

Medical Conditions Aggravated: Pre-existing: Cardiovascular diseases Kidney conditions Skin conditions Eye conditions Respiratory conditions

12. ECOLOGICAL INFORMATION

Product Ecological Information: --

No ecological data available for this product. Mobility in soil: Highly mobile Based on classification principles, not classified as hazardous to the environment. No bioaccumulation potential

Method Used for Estimation of Aquatic Toxicity of Mixture Summation Method M-factor (Multiplier) for highly toxic ingredients: 1

Ingredient Ecological Information: Glutaric dialdehyde: 96 hr Fish LC50 = 3.5 mg/L; 48 hr Crustaceans EC50 = 0.75 mg/L; Potassium Chloride: 48 hr Daphnia magna EC50 = 83 mg/L; 96 hr Pimephales promelas LC50 = 880 mg/L; 7 days Pimephales promelas NOEC = 500 mg/L

CEPA categorization for ingredients are as follows:

Glutaric Dialdehyde: Not persistent, bioaccumulative or inherently toxic to aquatic organisms

Potassium Chloride, Potassium Phosphate, Monobasic and Sodium Phosphate, Dibasic, Water: Persistent, not bioaccumulative or inherently toxic to aquatic organisms.

13. DISPOSAL CONSIDERATIONS

EPA Waste ID Number: Not applicable

Special Instructions (Disposal): Check with local municipal and state authorities and waste contractors for pertinent local information regarding the proper disposal of chemicals.

Empty Containers: Working in a well-ventilated area, Rinse three times with an appropriate solvent. Collect rinsate and dispose of according to local, state or federal regulations. In the US, rinsate from empty containers is classified as hazardous waste and should be disposed of at an E.P. A. approved facility. Rinsate from empty containers may contain sufficient product to require disposal as hazardous waste. Dispose of empty container as normal trash.

NOTICE (Disposal): These disposal guidelines are based on federal regulations and may be superseded by more stringent state or local requirements. Please consult your local environmental regulators for more information. In Europe: Chemical and analysis solutions must be disposed of in compliance with the respective national regulations. Product packaging must be disposed of in compliance with the country-specific regulations or must be passed to a packaging return system.

14. TRANSPORT INFORMATION

D.O.T.:

D.O.T. Proper Shipping Name: Not Currently Regulated

--
Hazard Class: NA
Subsidiary Risk: NA
ID Number: NA
Packing Group: NA
T.D.G.:
Proper Shipping Name: Not Currently Regulated
--
Hazard Class: NA
Subsidiary Risk: NA
UN Number/PIN: NA
Packing Group: NA
I.C.A.O.:
I.C.A.O. Proper Shipping Name: Not Currently Regulated
--
Hazard Class: NA
Subsidiary Risk: NA
ID Number: NA
Packing Group: NA
I.M.O.:
Proper Shipping Name: Not Currently Regulated
--
Hazard Class: NA
Subsidiary Risk: NA
ID Number: NA
Packing Group: NA
Additional Information: There is a possibility that this product could be contained in a reagent set or kit composed of various compatible dangerous goods. If the item is NOT in a set or kit, the classification given above applies. If the item IS part of a set or kit, the classification would change to the following: UN3316 Chemical Kit, Class 9, PG II or III. If the item is not regulated, the Chemical Kit classification does not apply.

15. REGULATORY INFORMATION

U.S. Federal Regulations:

O.S.H.A.: This product meets the criteria for a hazardous substance as defined in the Hazard Communication Standard [29 CFR 1910.1200(d)(5)(iv)].

E.P.A.:

S.A.R.A. Title III Section 311/312 Categorization (40 CFR 370): Immediate (Acute) Health Hazard

S.A.R.A. Title III Section 313 (40 CFR 372): This product does NOT contain any chemical subject to the reporting requirements of Section 313 of Title III of SARA.

--

302 (EHS) TPQ (40 CFR 355): Not applicable

304 CERCLA RQ (40 CFR 302.4): Sodium Phosphate, Dibasic 5000 lbs.

304 EHS RQ (40 CFR 355): Not applicable

Clean Water Act (40 CFR 116.4): Sodium phosphate, dibasic - RQ 5000 lbs.

RCRA: Contains no RCRA regulated substances.

State Regulations:

California Prop. 65: No Prop. 65 listed chemicals are present in this product.

Identification of Prop. 65 Ingredient(s): --

California Perchlorate Rule CCR Title 22 Chap 33: Not applicable

Trade Secret Registry: Not applicable

National Inventories:

U.S. Inventory Status: All ingredients in this product are listed on the TSCA 8(b) Inventory (40 CFR 710).

CAS Number: Not applicable

Canadian Inventory Status: All ingredients of this product are DSL Listed.

EEC Inventory Status: All ingredients used to make this product are listed on EINECS / ELINCS.

Australian Inventory (AICS) Status: All ingredients are listed.

New Zealand Inventory (NZIoC) Status: All components either listed or exempt.

Korean Inventory (KECI) Status: All components of this product are either listed, listed as the anhydrous compound or exempt.

Japan (ENCS) Inventory Status: All components either listed or exempt.

China (PRC) Inventory (MEP) Status: All components either listed or exempt.

16. OTHER INFORMATION

References: 29 CFR 1900 - 1910 (Code of Federal Regulations - Labor). CCINFO RTECS. Canadian Centre for Occupational Health and Safety. Hamilton, Ontario Canada: 30 June 1993. Technical Judgment.

Complete Text of H phrases referred to in Section 3: Not applicable

Revision Summary: . Substantial revision to comply with EU Reg 1272/2008, Reg 1907/2006 and UN GHS (ST/SG/AC.10/36/Add.3).

Date of MSDS Preparation:

Day: 08

Month: March

Year: 2014

MSDS Prepared: MSDS prepared by Product Compliance Department extension 3350

CCOHS Evaluation Note: It is offered under exemption from WHMIS labeling as specified in the Controlled Products Regulation (CPR) Section 17. It is offered under the interim policy that was established by Health Canada permitting use of GHS-formatted safety data sheets in Canada prior to revision of CPR to GHS. This product has been classified and labeled in accordance with the requirements of GHS (ST/SG/AC.10/36/Add.3). This SDS has been prepared in accordance with the requirements of GHS (ST/SG/AC.10/36/Add.3).

Legend:

NA - Not Applicable	w/w - weight/weight
ND - Not Determined	w/v - weight/volume
NV - Not Available	v/v - volume/volume

USER RESPONSIBILITY: Each user should read and understand this information and incorporate it in individual site safety programs in accordance with applicable hazard communication standards and regulations.

THE INFORMATION CONTAINED HEREIN IS BASED ON DATA CONSIDERED TO BE ACCURATE. HOWEVER, NO WARRANTY IS EXPRESSED OR IMPLIED REGARDING THE ACCURACY OF THESE DATA OR THE RESULTS TO BE OBTAINED FROM THE USE THEREOF.

HACH COMPANY ©2015

World Headquarters
Hach Company
P.O.Box 389
Loveland, CO USA 80539
(970) 669-3050

MSDS No: M00368

SAFETY DATA SHEET

1. CHEMICAL PRODUCT AND COMPANY IDENTIFICATION

Product Name: Buffer Solution pH 4.01 \pm 0.02
Catalog Number: 2283449

Hach Company
P.O.Box 389
Loveland, CO USA 80539
(970) 669-3050

Emergency Telephone Numbers:
(Medical and Transportation)
(303) 623-5716 24 Hour Service
(515)232-2533 8am - 4pm CST

MSDS Number: M00368
Chemical Name: Not applicable
CAS Number: Not applicable
Additional CAS No. (for hydrated forms): Not applicable
Chemical Formula: Not applicable
Chemical Family: Mixture
Intended Use: Laboratory Reagent Buffer

2. HAZARDS IDENTIFICATION

GHS Classification:
Hazard categories: Serious Eye Damage/Eye Irritation:Eye Irrit. 2
GHS Label Elements:
WARNING



Hazard statements: Causes serious eye irritation.
Precautionary statements: Wear protective gloves / protective clothing / eye protection / face protection. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention.

HMIS:
Health: 0
Flammability: 0
Reactivity: 0
Protective Equipment: X - See protective equipment, Section 8.
NFPA:
Health: 0
Flammability: 0
Reactivity: 0
Symbol: Not applicable
WHMIS Hazard Classification: Not applicable
WHMIS Symbols: Not applicable

3. COMPOSITION / INFORMATION ON INGREDIENTS

Hazardous Components according to GHS:
Potassium Acid Phthalate

CAS Number: 877-24-7

Chemical Formula: $C_8H_5KO_4$

GHS Classification: Acute Tox. 5, H303; Acute Tox. Derm. 4, H312; Eye Irrit. 2A, H319

Percent Range (Trade Secret): <1.5

Percent Range Units: weight / weight

PEL: 15 mg/m³ as total dust; 5 mg/m³ as respirable dust

TLV: 10 mg/m³ as inhalable dust; 3 mg/m³ as respirable dust

WHMIS Symbols: Not applicable

Formaldehyde

CAS Number: 50-00-0

Chemical Formula: CH_2O

GHS Classification: Flam. Liq. 4, H227; Acute Tox. 3 -Orl, H301; Acute Tox. 3 -Derm, H311; Skin Corr. 1B, H314; Skin Sens. 1, H317; Acute Tox. 3 -Inh, H331; Resp. Sens. 1, H334; Muta. 2, H341; Carc. 2, H351; Repr. 2, H361; STOT Single 1, H370; Aquatic Acute 2, H401

Percent Range (Trade Secret): < 0.1

Percent Range Units: weight / weight

PEL: 0.75 ppm

TLV: 0.3 ppm

WHMIS Symbols: Acute PoisonOther Toxic Effects

Methyl Alcohol

CAS Number: 67-56-1

Chemical Formula: CH_3OH

GHS Classification: Flam. Liq. 2, H225; Acute Tox. 3 -Orl, H301; Acute Tox. 3 -Derm, H311; Skin Irrit. 2, H315; Eye Irrit. 2A, H319; Acute Tox. 3 -Inh, H331; Muta. 2, H341; Repr. 2, H361; STOT SE1, H370

Percent Range (Trade Secret): < 0.1

Percent Range Units: weight / weight

PEL: 200 ppm

TLV: 200 ppm

WHMIS Symbols: Acute PoisonFlammable / CombustibleOther Toxic Effects

Hazardous Components according to GHS: No

Demineralized Water

CAS Number: 7732-18-5

Chemical Formula: H_2O

GHS Classification: Not a dangerous substance according to GHS.

Percent Range (Trade Secret): >98

Percent Range Units: weight / weight

PEL: Not established

TLV: Not established

WHMIS Symbols: Not applicable

SPADNS

CAS Number: 23647-14-5

Chemical Formula: $C_{16}H_9N_2O_{11}S_3 \cdot 3Na$

GHS Classification: Non-hazardous

Percent Range (Trade Secret): < 0.01

Percent Range Units: weight / weight

PEL: 15 mg/m³ as inhalable dust; 5 mg/m³ as respirable dust

TLV: 10 mg/m³ as inhalable dust; 3 mg/m³ as respirable dust

WHMIS Symbols: Not applicable

4. FIRST AID MEASURES

General Information: In the event of exposure, show this Material Safety Data Sheet and label (where possible) to a doctor.

Advice to doctor: Treat symptomatically.

Eye Contact: Immediately flush eyes with water for 15 minutes.

Skin Contact (First Aid): Wash skin with soap and plenty of water. Call physician if irritation develops. Remove contaminated clothing.

Inhalation: None required.

Ingestion (First Aid): Give large quantities of water. Never give anything by mouth to an unconscious person. Call physician immediately.

5. FIRE FIGHTING MEASURES

Flammable Properties: Material will not burn.

Fire Fighting Instruction: As in any fire, wear self-contained breathing apparatus pressure-demand and full protective gear.

Extinguishing Media: Use media appropriate to surrounding fire conditions

Extinguishing Media NOT To Be Used: Not applicable

Fire / Explosion Hazards: None reported

Hazardous Combustion Products: This material will not burn.

6. ACCIDENTAL RELEASE MEASURES

Spill Response Notice:

Only persons properly qualified to respond to an emergency involving hazardous substances may respond to a spill according to federal regulations (OSHA 29 CFR 1910.120(a)(v)) and per your company's emergency response plan and guidelines/procedures. See Section 13, Special Instructions for disposal assistance. Outside of the US, only persons properly qualified according to state or local regulations should respond to a spill involving chemicals.

Containment Technique: Stop spilled material from being released to the environment.

Clean-up Technique: If permitted by regulation, Cover spilled material with an alkali, such as soda ash or sodium bicarbonate. Scoop up slurry into a large beaker. Adjust to a pH between 6 and 9 with an alkali, such as soda ash or sodium bicarbonate. Flush reacted material to the drain with a large excess of water. Decontaminate the area of the spill with a soap solution. Otherwise, Pick up spill for disposal and place in a closed container. Dispose of in accordance with local, state and federal regulations or laws.

Evacuation Procedure: Evacuate as needed to perform spill clean-up. If conditions warrant, increase the size of the evacuation.

DOT Emergency Response Guide Number: Not applicable

7. HANDLING AND STORAGE

Handling: Avoid contact with eyes. Wash thoroughly after handling.

Storage: Keep container tightly closed when not in use.

Flammability Class: Not applicable

8. EXPOSURE CONTROLS / PERSONAL PROTECTION

Engineering Controls: Maintain general industrial hygiene practices when using this product.

Personal Protective Equipment:

Eye Protection: safety glasses with top and side shields

Skin Protection: disposable latex gloves. In the EU, the selected gloves must satisfy the specifications of EU Directive 89/686/EEC and standard EN 374 derived from it.

Inhalation Protection: adequate ventilation

Precautionary Measures: Avoid contact with: eyes. Wash thoroughly after handling.

TLV: Not established

PEL: Not established

For Occupational Exposure Limits (OEL) for ingredients, see section 3 - Composition/Information on Ingredients.:

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance: Clear, red liquid

Physical State: Liquid

Molecular Weight: Not applicable

Odor: None

Odor Threshold: Not applicable

pH: 4.01

Metal Corrosivity:

Corrosivity Classification: Not classified as corrosive to metals according to GHS criteria.

Steel: Not determined

Aluminum: Not determined

Specific Gravity/ Relative Density (water = 1; air =1): 1.002

Viscosity: ~ 1.0 mPa·s

Solubility:

Water: Soluble

Acid: Soluble

Other: Not determined

Partition Coefficient (n-octanol / water): Not applicable

Coefficient of Water / Oil: Not determined

Melting Point: ~ 0 °C (~ 32 °F)

Decomposition Temperature: Not applicable

Boiling Point: ~ 100 °C (~ 212 °F)

Vapor Pressure: 17 mm Hg at 20 °C (68 °F)

Vapor Density (air = 1): 0.62

Evaporation Rate (water = 1): 0.99

Volatile Organic Compounds Content: Not applicable

Flammable Properties: Material will not burn.

Flash Point: Not applicable

Method: Not applicable

Flammability Limits:

Lower Explosion Limits: Not applicable

Upper Explosion Limits: Not applicable

Autoignition Temperature: Not applicable

Explosive Properties:

Not classified according to GHS criteria.

Oxidizing Properties:

Not classified according to GHS criteria.

Reactivity Properties:

Not classified as self-reactive, pyrophoric, self-heating or emitting flammable gases in contact with water according to GHS criteria.

Gas under Pressure:

Not applicable

10. STABILITY AND REACTIVITY

Chemical Stability: Stable when stored under proper conditions.

Mechanical Impact: None reported

Static Discharge: None reported.

Reactivity / Incompatibility: None reported

Hazardous Decomposition: None reported

Conditions to Avoid: Extreme temperatures

11. TOXICOLOGICAL INFORMATION

Toxicokinetics, Metabolism and Distribution: No information available for mixture.

Toxicologically Synergistic Products: None reported

Acute Toxicity: Acute Toxicity Estimate (ATE) - Calculated from Ingredient Toxicity Data Practically Non-toxic Based on classification principles, the classification criteria are not met.

Specific Target Organ Toxicity - Single Exposure (STOT-SE): Based on classification principles, the classification criteria are not met.

Specific Target Organ Toxicity - Repeat Exposure (STOT-RE): Based on classification principles, the classification criteria are not met.

Skin Corrosion/Irritation: Based on classification principles, the classification criteria are not met.

Eye Damage: Irritating to eyes.

Sensitization: Based on classification principles, the classification criteria are not met.

CMR Effects/Properties (carcinogenic, mutagenic or toxic to reproduction): Data supporting mutagenicity was found. Data insufficient for classification

Potassium Acid Phthalate has been listed as a mutagen and reproductive toxin on the Danish Advisory List (QSAR).

There is no mammalian data to support classification.

An ingredient of this mixture is: IARC Group 1: Recognized Carcinogen

Formaldehyde

An ingredient of this mixture is: NTP Listed Group 1: Recognized Carcinogen

Formaldehyde

An ingredient of this product is an OSHA listed carcinogen.

Formaldehyde

Symptoms/Effects:

Ingestion: No Effects Anticipated

Inhalation: No effects anticipated

Skin Absorption: No effects anticipated

Chronic Effects: No effects anticipated

Medical Conditions Aggravated: None reported

12. ECOLOGICAL INFORMATION

Product Ecological Information: --

No ecological data available for this product. No bioaccumulation potential Mobility in soil: Highly mobile Based on classification principles, not classified as hazardous to the environment.

Ingredient Ecological Information: --

Ecological data for ingredients is not indicative of likely ecological harm.

13. DISPOSAL CONSIDERATIONS

EPA Waste ID Number: Not applicable

Special Instructions (Disposal): Adjust to a pH between 6 and 9 with an alkali, such as soda ash or sodium bicarbonate.

If permitted by regulation, Open cold water tap completely, slowly pour the reacted material to the drain. Otherwise, Check with local municipal and state authorities and waste contractors for pertinent local information regarding the proper disposal of chemicals.

Empty Containers: Rinse three times with an appropriate solvent. Collect rinsate and dispose of according to local, state or federal regulations. Dispose of empty container as normal trash. In the US, rinsate from empty containers is classified as hazardous waste and should be disposed of at an E.P. A. approved facility. Rinsate from empty containers may contain sufficient product to require disposal as hazardous waste.

NOTICE (Disposal): These disposal guidelines are based on federal regulations and may be superseded by more stringent state or local requirements. Please consult your local environmental regulators for more information. In Europe: Chemical and analysis solutions must be disposed of in compliance with the respective national regulations. Product packaging must be disposed of in compliance with the country-specific regulations or must be passed to a packaging return system.

14. TRANSPORT INFORMATION

D.O.T.:

D.O.T. Proper Shipping Name: Not Currently Regulated

--

Hazard Class: NA

Subsidiary Risk: NA

ID Number: NA

Packing Group: NA

T.D.G.:

Proper Shipping Name: Not Currently Regulated

--

Hazard Class: NA

Subsidiary Risk: NA

UN Number/PIN: NA

Packing Group: NA

I.C.A.O.:

I.C.A.O. Proper Shipping Name: Not Currently Regulated

--

Hazard Class: NA
Subsidiary Risk: NA
ID Number: NA
Packing Group: NA
I.M.O.:
Proper Shipping Name: Not Currently Regulated
--

Hazard Class: NA
Subsidiary Risk: NA
ID Number: NA
Packing Group: NA

Additional Information: There is a possibility that this product could be contained in a reagent set or kit composed of various compatible dangerous goods. If the item is NOT in a set or kit, the classification given above applies. If the item IS part of a set or kit, the classification would change to the following: UN3316 Chemical Kit, Class 9, PG II or III. If the item is not regulated, the Chemical Kit classification does not apply.

15. REGULATORY INFORMATION

U.S. Federal Regulations:

O.S.H.A.: This product does not meet the criteria for a hazardous substance as defined in the Hazard Communication Standard. (29 CFR 1910.1200)

E.P.A.:

S.A.R.A. Title III Section 311/312 Categorization (40 CFR 370): This product is not hazardous under 29 CFR.1910.1200 and therefore is not covered by Title III under SARA.

S.A.R.A. Title III Section 313 (40 CFR 372): This product does NOT contain any chemical subject to the reporting requirements of Section 313 of Title III of SARA.

--

302 (EHS) TPQ (40 CFR 355): Not applicable

304 CERCLA RQ (40 CFR 302.4): Not applicable

304 EHS RQ (40 CFR 355): Not applicable

Clean Water Act (40 CFR 116.4): Not applicable

RCRA: Contains no RCRA regulated substances.

State Regulations:

California Prop. 65: No Prop. 65 listed chemicals are present in this product.

Identification of Prop. 65 Ingredient(s): None

California Perchlorate Rule CCR Title 22 Chap 33: Not applicable

Trade Secret Registry: Not applicable

National Inventories:

U.S. Inventory Status: All ingredients in this product are listed on the TSCA 8(b) Inventory (40 CFR 710).

CAS Number: Not applicable

Canadian Inventory Status: All ingredients of this product are DSL Listed.

EEC Inventory Status: All ingredients used to make this product are listed on EINECS / ELINCS.

Australian Inventory (AICS) Status: All ingredients are listed.

New Zealand Inventory (NZIoC) Status: All components either listed or exempt.

Korean Inventory (KECI) Status: All components of this product are either listed, listed as the anhydrous compound or exempt.

Japan (ENCs) Inventory Status: All components either listed or exempt.

China (PRC) Inventory (MEP) Status: All components either listed or exempt.

16. OTHER INFORMATION

References: 29 CFR 1900 - 1910 (Code of Federal Regulations - Labor). Air Contaminants, Federal Register, Vol. 54, No. 12. Thursday, January 19, 1989. pp. 2332-2983. TLV's Threshold Limit Values and Biological Exposure Indices for 1992-1993. American Conference of Governmental Industrial Hygienists, 1992. Technical Judgment. In-house information. Fire Protection Guide on Hazardous Materials, 10th Ed. Quincy, MA: National Fire Protection Fire Protection Guide on Hazardous Materials, 10th Ed. Quincy, MA: National Fire Protection Association, 1991.
Complete Text of H phrases referred to in Section 3: H225 Highly flammable liquid and vapour. H227 Combustible liquid. H301 Toxic if swallowed. Not applicable H311 Toxic in contact with skin. H312 Harmful in contact with skin. H314 Causes severe skin burns and eye damage. H315 Causes skin irritation. H317 May cause an allergic skin reaction. H319 Causes serious eye irritation. H331 Toxic if inhaled. H334 May cause allergy or asthma symptoms or breathing

difficulties if inhaled. H341 Suspected of causing genetic defects. H361fd Suspected of damaging fertility. Suspected of damaging the unborn child. H370 Causes damage to organs. H401 Toxic to aquatic life.

Revision Summary: Substantial revision to comply with EU Reg 1272/2008, Reg 1907/2006 and UN GHS (ST/SG/AC.10/36/Add.3).

Date of MSDS Preparation:

Day: 13

Month: March

Year: 2014

MSDS Prepared: MSDS prepared by Product Compliance Department extension 3350

CCOHS Evaluation Note: It is offered under exemption from WHMIS labeling as specified in the Controlled Products Regulation (CPR) Section 17. It is offered under the interim policy that was established by Health Canada permitting use of GHS-formatted safety data sheets in Canada prior to revision of CPR to GHS. This product has been classified and labeled in accordance with the requirements of GHS (ST/SG/AC.10/36/Add.3). This SDS has been prepared in accordance with the requirements of GHS (ST/SG/AC.10/36/Add.3).

Legend:

NA - Not Applicable	w/w - weight/weight
ND - Not Determined	w/v - weight/volume
NV - Not Available	v/v - volume/volume

USER RESPONSIBILITY: Each user should read and understand this information and incorporate it in individual site safety programs in accordance with applicable hazard communication standards and regulations.

THE INFORMATION CONTAINED HEREIN IS BASED ON DATA CONSIDERED TO BE ACCURATE. HOWEVER, NO WARRANTY IS EXPRESSED OR IMPLIED REGARDING THE ACCURACY OF THESE DATA OR THE RESULTS TO BE OBTAINED FROM THE USE THEREOF.

HACH COMPANY ©2015

World Headquarters
Hach Company
P.O.Box 389
Loveland, CO USA 80539
(970) 669-3050

MSDS No: M00369

SAFETY DATA SHEET

1. CHEMICAL PRODUCT AND COMPANY IDENTIFICATION

Product Name: Buffer Solution pH 7.00 \pm 0.02
Catalog Number: 2283549

Hach Company
P.O.Box 389
Loveland, CO USA 80539
(970) 669-3050

Emergency Telephone Numbers:
(Medical and Transportation)
(303) 623-5716 24 Hour Service
(515)232-2533 8am - 4pm CST

MSDS Number: M00369
Chemical Name: Not applicable
CAS Number: Not applicable
Additional CAS No. (for hydrated forms): Not applicable
Chemical Formula: Not applicable
Chemical Family: Not applicable
Intended Use: Laboratory Reagent Buffer

2. HAZARDS IDENTIFICATION

This mixture is not classified as hazardous per GHS (UN publication ST/SG/AC.10/36/Add.3)

GHS Classification:
Hazard categories: Not applicable
GHS Label Elements:
Not applicable

Hazard statements: Not applicable
Precautionary statements: Not applicable

HMIS:
Health: 0
Flammability: 0
Reactivity: 0
Protective Equipment: X - See protective equipment, Section 8.

NFPA:
Health: 0
Flammability: 0
Reactivity: 0
Symbol: Not applicable
WHMIS Hazard Classification: Not applicable
WHMIS Symbols: Not applicable

3. COMPOSITION / INFORMATION ON INGREDIENTS

Hazardous Components according to GHS:

Potassium Phosphate, Monobasic

CAS Number: 7778-77-0
Chemical Formula: KH_2PO_4
GHS Classification: Acute Tox. 4 - Orl, H302; Eye Irrit. 2, H319
Percent Range (Trade Secret): < 0.5
Percent Range Units: weight / weight

PEL: 15 mg/m³ as total dust; 5 mg/m³ as respirable dust
TLV: 10 mg/m³ as inhalable fraction; 3 mg/m³ as respirable fraction

WHMIS Symbols: Not applicable

Sodium Phosphate, Dibasic

CAS Number: 7558-79-4
Chemical Formula: Na₂HPO₄
GHS Classification: Skin Irrit. 2, H315; Eye Irrit. 2A, H319
Percent Range (Trade Secret): < 0.5
Percent Range Units: weight / weight
PEL: 15 mg/m³ as inhalable dust; 5 mg/m³ as respirable dust
TLV: 10 mg/m³ as inhalable dust; 3 mg/m³ as respirable dust

WHMIS Symbols: Not applicable

Magnesium Nitrate Hexahydrate

CAS Number: 10377-60-3
Chemical Formula: Mg(NO₃)₂ · 6H₂O
GHS Classification: Ox. Sol. 3, H272; Skin Irrit. 2, H315; Eye Irrit. 2A, H319; STOT Single 3, H335
Percent Range (Trade Secret): < 0.1
Percent Range Units: weight / weight
PEL: 15 mg/m³ as inhalable dust; 5 mg/m³ as respirable dust
TLV: 10 mg/m³ as inhalable dust; 3 mg/m³ as respirable dust

WHMIS Symbols: Oxidizing

Tartrazine

CAS Number:
Chemical Formula: C₁₆H₁₂N₄O₉S₂ · 3 Na
GHS Classification:
Percent Range (Trade Secret): < 0.01
Percent Range Units: weight / weight
PEL: 15 mg/m³ as total dust; 5 mg/m³ as respirable dust
TLV: 10 mg/m³ as inhalable dust; 3 mg/m³ as respirable dust

WHMIS Symbols: Other Toxic Effects

5-Chloro-2-methyl-2H-isothiazol-3-one

CAS Number: 26172-55-4
Chemical Formula: C₄H₄ClNOS
GHS Classification: Acute Tox. 4, H302; Skin Corr. 1B, H314; Resp. Sens. 1, H317; Skin Sens. 1, H334; Aquatic Acute 1, H400
Percent Range (Trade Secret): < 0.01
Percent Range Units: weight / weight
PEL: Not established
TLV: Not established

WHMIS Symbols:

2-Methyl-2H-isothiazol-3-one

CAS Number: 2682-20-4
Chemical Formula: C₄H₅NOS
GHS Classification: Acute Tox. 4-Oral., H302; Acute Tox. 3-Inh., H331; Skin Corr. 1B, H314; Skin Sens. 1, H317; STOT SE 3, H335; Aquatic Acute 1, H400
Percent Range (Trade Secret): < 0.001
Percent Range Units: weight / weight
PEL: Not determined
TLV: Not established

WHMIS Symbols:

Hazardous Components according to GHS: No

Demineralized Water

CAS Number: 7732-18-5

Chemical Formula: H₂O

GHS Classification: Not a dangerous substance according to GHS.

Percent Range (Trade Secret): > 99.0

Percent Range Units: weight / weight

PEL: Not established

TLV: Not established

WHMIS Symbols: Not applicable

Magnesium Chloride

CAS Number: 7786-30-3

Chemical Formula: MgCl₂ · 6H₂O

GHS Classification: Not hazardous per GHS

Percent Range (Trade Secret): < 0.01

Percent Range Units: weight / weight

PEL: Not established

TLV: Not established

WHMIS Symbols: Not applicable

4. FIRST AID MEASURES

General Information: In the event of exposure, show this Material Safety Data Sheet and label (where possible) to a doctor.

Advice to doctor: Treat symptomatically.

Eye Contact: Call physician if irritation develops. Immediately flush eyes with water for 15 minutes.

Skin Contact (First Aid): Wash skin with plenty of water. Call physician if irritation develops. Remove contaminated clothing.

Inhalation: Remove to fresh air.

Ingestion (First Aid): Give large quantities of water. Never give anything by mouth to an unconscious person. Call physician immediately.

5. FIRE FIGHTING MEASURES

Flammable Properties: Material will not burn. Material is not classified as flammable according to GHS criteria. Material will not burn.

Fire Fighting Instruction: As in any fire, wear self-contained breathing apparatus pressure-demand and full protective gear.

Extinguishing Media: Use media appropriate to surrounding fire conditions

Extinguishing Media NOT To Be Used: Not applicable

Fire / Explosion Hazards: None reported

Hazardous Combustion Products: None reported

6. ACCIDENTAL RELEASE MEASURES

Spill Response Notice:

Only persons properly qualified to respond to an emergency involving hazardous substances may respond to a spill according to federal regulations (OSHA 29 CFR 1910.120(a)(v)) and per your company's emergency response plan and guidelines/procedures. See Section 13, Special Instructions for disposal assistance. Outside of the US, only persons properly qualified according to state or local regulations should respond to a spill involving chemicals.

Containment Technique: Stop spilled material from being released to the environment.

Clean-up Technique: If permitted by regulation, Absorb spilled liquid with non-reactive sorbent material. Place material in a plastic bag. Mark bag 'Non-hazardous trash', and dispose of as normal refuse. Otherwise, Pick up spill for disposal and place in a closed container. Dispose of in accordance with local, state and federal regulations or laws.

Evacuation Procedure: Evacuate as needed to perform spill clean-up. If conditions warrant, increase the size of the evacuation.

DOT Emergency Response Guide Number: Not applicable

7. HANDLING AND STORAGE

Handling: Avoid contact with eyes Wash thoroughly after handling.

Storage: Protect from: heat Keep container tightly closed when not in use.

Flammability Class: Not applicable

8. EXPOSURE CONTROLS / PERSONAL PROTECTION

Engineering Controls: Maintain general industrial hygiene practices when using this product.

Personal Protective Equipment:

Eye Protection: safety glasses with top and side shields

Skin Protection: disposable latex gloves

Inhalation Protection: adequate ventilation

Precautionary Measures: Avoid contact with: eyes Wash thoroughly after handling.

TLV: Not established

PEL: Not established

For Occupational Exposure Limits (OEL) for ingredients, see section 3 - Composition/Information on Ingredients.:

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance: Clear, yellow liquid

Physical State: Liquid

Molecular Weight: Not applicable

Odor: None

Odor Threshold: Odorless

pH: 7.0 at 20 °C (68 °F)

Metal Corrosivity:

Corrosivity Classification: Not classified as corrosive to metals according to GHS criteria.

Steel: Not determined

Aluminum: Not determined

Specific Gravity/ Relative Density (water = 1; air = 1): ~ 1.0

Viscosity: ~ 1.0 mPa·s

Solubility:

Water: Soluble

Acid: Soluble

Other: Not determined

Partition Coefficient (n-octanol / water): Not applicable

Coefficient of Water / Oil: Not determined

Melting Point: ~ 0 °C (~ 32 °F)

Decomposition Temperature: Not determined

Boiling Point: ~ 100 °C (~ 212 °F)

Vapor Pressure: ~ 18 mm Hg at 20 °C (68 °F)

Vapor Density (air = 1): 0.62

Evaporation Rate (water = 1): ~ 1

Volatile Organic Compounds Content: Not applicable

Flammable Properties: Material is not classified as flammable according to GHS criteria. Material will not burn.

Flash Point: Not applicable

Method: Not applicable

Flammability Limits:

Lower Explosion Limits: Not applicable

Upper Explosion Limits: Not applicable

Autoignition Temperature: Not applicable

Explosive Properties:

Not classified according to GHS criteria.

Oxidizing Properties:

Not classified according to GHS criteria.

Reactivity Properties:

Not classified as self-reactive, pyrophoric, self-heating or emitting flammable gases in contact with water according to GHS criteria.

Gas under Pressure:

Not classified according to GHS criteria.

10. STABILITY AND REACTIVITY

Chemical Stability: Stable when stored under proper conditions.

Mechanical Impact: None reported

Static Discharge: None reported.

Reactivity / Incompatibility: None reported

Hazardous Decomposition: None reported

Conditions to Avoid: Heat Evaporation

11. TOXICOLOGICAL INFORMATION

Toxicokinetics, Metabolism and Distribution: No information available for mixture. No health effects are anticipated in normal use.

Toxicologically Synergistic Products: None reported

Acute Toxicity: Practically Non-toxic Based on classification principles, the classification criteria are not met.

Specific Target Organ Toxicity - Single Exposure (STOT-SE): Based on classification principles, the classification criteria are not met.

Specific Target Organ Toxicity - Repeat Exposure (STOT-RE): Based on classification principles, the classification criteria are not met.

Skin Corrosion/Irritation: Based on classification principles, the classification criteria are not met.

Eye Damage: Based on classification principles, the classification criteria are not met.

Sensitization: Based on classification principles, the classification criteria are not met.

CMR Effects/Properties (carcinogenic, mutagenic or toxic to reproduction): Based on classification principles, the classification criteria are not met.

This product does NOT contain any IARC listed chemicals.

This product does NOT contain any NTP listed chemicals.

This product does NOT contain any OSHA listed carcinogens.

Symptoms/Effects:

Ingestion: No Effects Anticipated

Inhalation: No effects anticipated

Skin Absorption: No effects anticipated

Chronic Effects: No effects anticipated

Medical Conditions Aggravated: None reported

12. ECOLOGICAL INFORMATION

Product Ecological Information: --

No ecological data available for this product. Mobility in soil: Highly mobile No bioaccumulation potential Based on classification principles, not classified as hazardous to the environment.

Ingredient Ecological Information:

Ecological data for ingredients is not indicative of likely ecological harm. CEPA categorization for each and every ingredient: Persistent Not bioaccumulative and not inherently toxic to aquatic organisms.

13. DISPOSAL CONSIDERATIONS

EPA Waste ID Number: Not applicable

Special Instructions (Disposal): If permitted by regulation, Open cold water tap completely, slowly pour the material to the drain. Otherwise, Check with local municipal and state authorities and waste contractors for pertinent local information regarding the proper disposal of chemicals.

Empty Containers: Rinse three times with an appropriate solvent. Collect rinsate and dispose of according to local, state or federal regulations. Dispose of empty container as normal trash. In the US, rinsate from empty containers is classified as hazardous waste and should be disposed of at an E.P. A. approved facility. Rinsate from empty containers may contain sufficient product to require disposal as hazardous waste.

NOTICE (Disposal): These disposal guidelines are based on federal regulations and may be superseded by more stringent state or local requirements. Please consult your local environmental regulators for more information. In Europe: Chemical and analysis solutions must be disposed of in compliance with the respective national regulations. Product packaging must be disposed of in compliance with the country-specific regulations or must be passed to a packaging return system.

14. TRANSPORT INFORMATION

D.O.T.:

D.O.T. Proper Shipping Name: Not Currently Regulated

--

Hazard Class: NA

Subsidiary Risk: NA

ID Number: NA

Packing Group: NA

T.D.G.:

Proper Shipping Name: Not Currently Regulated

--

Hazard Class: NA

Subsidiary Risk: NA

UN Number/PIN: NA

Packing Group: NA

I.C.A.O.:

I.C.A.O. Proper Shipping Name: Not Currently Regulated

--

Hazard Class: NA

Subsidiary Risk: NA

ID Number: NA

Packing Group: NA

I.M.O.:

Proper Shipping Name: Not Currently Regulated

--

Hazard Class: NA

Subsidiary Risk: NA

ID Number: NA

Packing Group: NA

Additional Information: There is a possibility that this product could be contained in a reagent set or kit composed of various compatible dangerous goods. If the item is NOT in a set or kit, the classification given above applies. If the item IS part of a set or kit, the classification would change to the following: UN3316 Chemical Kit, Class 9, PG II or III. If the item is not regulated, the Chemical Kit classification does not apply.

15. REGULATORY INFORMATION

U.S. Federal Regulations:

O.S.H.A.: This product does not meet the criteria for a hazardous substance as defined in the Hazard Communication Standard. (29 CFR 1910.1200)

E.P.A.:

S.A.R.A. Title III Section 311/312 Categorization (40 CFR 370): This product is not hazardous under 29 CFR.1910.1200 and therefore is not covered by Title III under SARA.

S.A.R.A. Title III Section 313 (40 CFR 372): This product does NOT contain any chemical subject to the reporting requirements of Section 313 of Title III of SARA.

302 (EHS) TPQ (40 CFR 355): Not applicable

304 CERCLA RQ (40 CFR 302.4): Not applicable

304 EHS RQ (40 CFR 355): Not applicable

Clean Water Act (40 CFR 116.4): Not applicable

RCRA: Contains no RCRA regulated substances.

State Regulations:

California Prop. 65: No Prop. 65 listed chemicals are present in this product.

Identification of Prop. 65 Ingredient(s): None

California Perchlorate Rule CCR Title 22 Chap 33: Not applicable

Trade Secret Registry: Not applicable

National Inventories:

U.S. Inventory Status: All ingredients in this product are listed on the TSCA 8(b) Inventory (40 CFR 710).

CAS Number: Not applicable

Canadian Inventory Status: All ingredients of this product are DSL Listed.

EEC Inventory Status: All ingredients used to make this product are listed on EINECS / ELINCS.

Australian Inventory (AICS) Status: All ingredients are listed.

New Zealand Inventory (NZIoC) Status: All components either listed or exempt.

Korean Inventory (KECI) Status: All components of this product are either listed, listed as the anhydrous compound or exempt.

Japan (ENCS) Inventory Status: All components either listed or exempt.

China (PRC) Inventory (MEP) Status: All components either listed or exempt.

16. OTHER INFORMATION

References: 29 CFR 1900 - 1910 (Code of Federal Regulations - Labor). Air Contaminants, Federal Register, Vol. 54, No. 12. Thursday, January 19, 1989. pp. 2332-2983. TLV's Threshold Limit Values and Biological Exposure Indices for 1992-1993. American Conference of Governmental Industrial Hygienists, 1992. Technical Judgment. In-house information.

Complete Text of H phrases referred to in Section 3: H301 Toxic if swallowed. H302 Harmful if swallowed. H311 Toxic in contact with skin. H314 Causes severe skin burns and eye damage. H315 Causes skin irritation. H317 May cause an allergic skin reaction. H319 Causes serious eye irritation. H331 Toxic if inhaled. H351 Suspected of causing cancer.

Revision Summary: Substantial revision to comply with EU Reg 1272/2008, Reg 1907/2006 and UN GHS (ST/SG/AC.10/36/Add.3).

Date of MSDS Preparation:

Day: 01

Month: June

Year: 2015

MSDS Prepared: MSDS prepared by Product Compliance Department extension 3350

CCOHS Evaluation Note: It is offered under exemption from WHMIS labeling as specified in the Controlled Products Regulation (CPR) Section 17. It is offered under the interim policy that was established by Health Canada permitting use of GHS-formatted safety data sheets in Canada prior to revision of CPR to GHS. This product has been classified and labeled in accordance with the requirements of GHS (ST/SG/AC.10/36/Add.3). This SDS has been prepared in accordance with the requirements of GHS (ST/SG/AC.10/36/Add.3).

Legend:

NA - Not Applicable

ND - Not Determined

NV - Not Available

w/w - weight/weight

w/v - weight/volume

v/v - volume/volume

USER RESPONSIBILITY: Each user should read and understand this information and incorporate it in individual site safety programs in accordance with applicable hazard communication standards and regulations.

THE INFORMATION CONTAINED HEREIN IS BASED ON DATA CONSIDERED TO BE ACCURATE. HOWEVER, NO WARRANTY IS EXPRESSED OR IMPLIED REGARDING THE ACCURACY OF THESE DATA OR THE RESULTS TO BE OBTAINED FROM THE USE THEREOF.

HACH COMPANY ©2015

World Headquarters
Hach Company
P.O.Box 389
Loveland, CO USA 80539
(970) 669-3050

MSDS No: M00370

SAFETY DATA SHEET

1. CHEMICAL PRODUCT AND COMPANY IDENTIFICATION

Product Name: Buffer Solution pH 10.01 \pm 0.02
Catalog Number: 2283649

Hach Company
P.O.Box 389
Loveland, CO USA 80539
(970) 669-3050

Emergency Telephone Numbers:
(Medical and Transportation)
(303) 623-5716 24 Hour Service
(515)232-2533 8am - 4pm CST

MSDS Number: M00370
Chemical Name: Not applicable
CAS Number: Not applicable
Additional CAS No. (for hydrated forms): Not applicable
Chemical Formula: Not applicable
Chemical Family: Mixture
Intended Use: Buffer

2. HAZARDS IDENTIFICATION

This mixture is not classified as hazardous per GHS (UN publication ST/SG/AC.10/36/Add.3)

GHS Classification:
Hazard categories: Not applicable
GHS Label Elements:
Not applicable

Hazard statements: Not applicable
Precautionary statements: Not applicable

HMIS:
Health: 0
Flammability: 0
Reactivity: 0
Protective Equipment: X - See protective equipment, Section 8.

NFPA:
Health: 0
Flammability: 0
Reactivity: 0
Symbol: Not applicable
WHMIS Hazard Classification: Not applicable
WHMIS Symbols: Not applicable

3. COMPOSITION / INFORMATION ON INGREDIENTS

Hazardous Components according to GHS:

Sodium Bicarbonate

CAS Number: 144-55-8
Chemical Formula: NaHCO₃
GHS Classification: Acute Tox. 5-Orl, H303; Skin irrit. 3, H316
Percent Range (Trade Secret): < 1.0
Percent Range Units: weight / weight

PEL: 15 mg/m³ as inhalable dust; 5 mg/m³ as respirable dust
TLV: 10 mg/m³ as inhalable dust; 3 mg/m³ as respirable dust

WHMIS Symbols: Not applicable

Formaldehyde

CAS Number: 50-00-0

Chemical Formula: CH₂O

GHS Classification: Flam. Liq. 4, H227; Acute Tox. 3 -Orl, H301; Acute Tox. 3 -Derm, H311; Skin Corr. 1B, H314; Skin Sens. 1, H317; Acute Tox. 3-Inh, H331; Resp. Sens. 1, H334; Muta. 2, H341; Carc. 2, H351; Repr. 2, H361; STOT Single 1, H370; Aquatic Acute 2, H401

Percent Range (Trade Secret): < 0.1

Percent Range Units: weight / weight

PEL: 0.75 ppm

TLV: 0.3 ppm

WHMIS Symbols: Acute PoisonOther Toxic Effects

Methyl Alcohol

CAS Number:

Chemical Formula: CH₃OH

GHS Classification:

Percent Range (Trade Secret): < 0.1

Percent Range Units: weight / weight

PEL: 200 ppm

TLV: 200 ppm

WHMIS Symbols: Acute PoisonFlammable / CombustibleOther Toxic Effects

Sodium Carbonate

CAS Number: 497-19-8

Chemical Formula: Na₂CO₃

GHS Classification: Eye Irrit. 2, H319; Acute Tox. Inh. 4, H332; Acute Tox. Orl. 5, H303

Percent Range (Trade Secret): < 0.1

Percent Range Units: weight / weight

PEL: Not established

TLV: Not established

WHMIS Symbols: Other Toxic Effects

Dye, Turquoise Blue Pyralor S-400

CAS Number: 1330-38-7

Chemical Formula: C₃₂H₁₄CuN₈O₆S₂.2Na

GHS Classification: Acute Tox. 5-Orl, BH303; Skin Irrit. 2, H315; Eye Irrit. 2, H319; STOT SE 3, H335; Aquatic Chronic 3, H412

Percent Range (Trade Secret): < 0.01

Percent Range Units: weight / weight

PEL: 15 mg/m³ as inhalable dust; 5 mg/m³ as respirable dust

TLV: 10 mg/m³ as inhalable dust; 3 mg/m³ as respirable dust

WHMIS Symbols: Not applicable

Hazardous Components according to GHS: No

Demineralized Water

CAS Number: 7732-18-5

Chemical Formula: H₂O

GHS Classification: Not a dangerous substance according to GHS.

Percent Range (Trade Secret): > 99.0

Percent Range Units: weight / weight

PEL: Not established

TLV: Not established

WHMIS Symbols: Not applicable

4. FIRST AID MEASURES

General Information: In the event of exposure, show this Material Safety Data Sheet and label (where possible) to a doctor.

Advice to doctor: Treat symptomatically.

Eye Contact: Flush eyes with water. Call physician if irritation develops.

Skin Contact (First Aid): Wash skin with plenty of water. Call physician if irritation develops. Remove contaminated clothing.

Inhalation: Remove to fresh air.

Ingestion (First Aid): Never give anything by mouth to an unconscious person. Rinse mouth with plenty of water. Give large quantities of water. Call physician immediately.

5. FIRE FIGHTING MEASURES

Flammable Properties: Material will not burn. Material is not classified as flammable according to GHS criteria.

Fire Fighting Instruction: As in any fire, wear self-contained breathing apparatus pressure-demand and full protective gear.

Extinguishing Media: Use media appropriate to surrounding fire conditions

Extinguishing Media NOT To Be Used: Not applicable

Fire / Explosion Hazards: This product will not burn or explode. May react violently with: strong acids strong bases

Hazardous Combustion Products: carbon monoxide, carbon dioxide.

6. ACCIDENTAL RELEASE MEASURES

Spill Response Notice:

Only persons properly qualified to respond to an emergency involving hazardous substances may respond to a spill according to federal regulations (OSHA 29 CFR 1910.120(a)(v)) and per your company's emergency response plan and guidelines/procedures. See Section 13, Special Instructions for disposal assistance.

Containment Technique: Stop spilled material from being released to the environment. Absorb spilled liquid with non-reactive sorbent material.

Clean-up Technique: If permitted by regulation, Cover spilled material with a dry acid, such as citric or boric. Scoop up slurry into a large beaker. Adjust to a pH between 6 and 9 with an acid, such as sulfuric or citric. Flush reacted material to the drain with a large excess of water. Otherwise, Pick up spill for disposal and place in a closed container. Dispose of in accordance with local, state and federal regulations or laws.

Evacuation Procedure: Evacuate as needed to perform spill clean-up. If conditions warrant, increase the size of the evacuation.

DOT Emergency Response Guide Number: Not applicable

7. HANDLING AND STORAGE

Handling: Avoid contact with eyes Wash thoroughly after handling. Maintain general industrial hygiene practices when using this product.

Storage: Protect from: heat Keep container tightly closed when not in use. Keep away from: acids bases oxidizers

Flammability Class: Not applicable

8. EXPOSURE CONTROLS / PERSONAL PROTECTION

Engineering Controls: Maintain general industrial hygiene practices when using this product.

Personal Protective Equipment:

Eye Protection: safety glasses with top and side shields

Skin Protection: disposable latex gloves

Inhalation Protection: adequate ventilation

Precautionary Measures: Avoid contact with: eyes Wash thoroughly after handling. Protect from: heat Keep away from: acids/acid fumes bases oxidizers

TLV: Not established

PEL: Not established

For Occupational Exposure Limits (OEL) for ingredients, see section 3 - Composition/Information on Ingredients.:

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance: Clear, blue

Physical State: Liquid

Molecular Weight: Not applicable

Odor: Odorless

Odor Threshold: Not applicable

pH: 10.0

Metal Corrosivity:

Corrosivity Classification: Not classified as corrosive to metals according to GHS criteria.

Steel: Not determined

Aluminum: Not determined

Specific Gravity/ Relative Density (water = 1; air =1): 0.990

Viscosity: ~ 1.0 mPa·s

Solubility:

Water: Soluble

Acid: Soluble

Other: Not determined

Partition Coefficient (n-octanol / water): Not applicable

Coefficient of Water / Oil: Not applicable

Melting Point: ~ 0 °C (~ 32 °F)

Decomposition Temperature: Not applicable

Boiling Point: ~ 100 °C (~ 212 °F)

Vapor Pressure: 17.5 mm Hg (2.27 kPa) at 20 °C (68 °F)

Vapor Density (air = 1): 0.62

Evaporation Rate (water = 1): 0.76

Volatile Organic Compounds Content: < 0.1%

Flammable Properties: Material will not burn. Material is not classified as flammable according to GHS criteria.

Flash Point: Not applicable

Method: Not applicable

Flammability Limits:

Lower Explosion Limits: Not applicable

Upper Explosion Limits: Not applicable

Autoignition Temperature: Not applicable

Explosive Properties:

Not classified according to GHS criteria.

Oxidizing Properties:

Not classified according to GHS criteria.

Reactivity Properties:

Not classified as self-reactive, pyrophoric, self-heating or emitting flammable gases in contact with water according to GHS criteria.

Gas under Pressure:

Not classified according to GHS criteria.

10. STABILITY AND REACTIVITY

Chemical Stability: Stable when stored under proper conditions.

Mechanical Impact: No

Static Discharge: No

Reactivity / Incompatibility: None reported

Hazardous Decomposition: Heating to decomposition releases toxic and/or corrosive fumes of: carbon dioxide carbon monoxide formaldehyde sodium oxides nitrogen oxides

Conditions to Avoid: Heat Evaporation Extreme temperatures Contact with acid or acid fumes Contact with heat, sparks, open flames or other ignition sources.

11. TOXICOLOGICAL INFORMATION

Toxicokinetics, Metabolism and Distribution: No information available for mixture. No health effects are anticipated in normal use.

Toxicologically Synergistic Products: None reported

Acute Toxicity: Acute Toxicity Estimate (ATE) - Calculated from Ingredient Toxicity Data Practically Non-toxic Based on classification principles, the classification criteria are not met.

Specific Target Organ Toxicity - Single Exposure (STOT-SE): Based on classification principles, the classification criteria are not met.

Specific Target Organ Toxicity - Repeat Exposure (STOT-RE): Based on classification principles, the classification criteria are not met.

Skin Corrosion/Irritation: Based on classification principles, the classification criteria are not met.

Eye Damage: Based on classification principles, the classification criteria are not met.

Sensitization: Based on classification principles, the classification criteria are not met.

CMR Effects/Properties (carcinogenic, mutagenic or toxic to reproduction): Based on classification principles, the classification criteria are not met.

An ingredient of this mixture is: IARC Group 1: Recognized Carcinogen
Formaldehyde

An ingredient of this mixture is: NTP Listed Group 1: Recognized Carcinogen
Formaldehyde

An ingredient of this product is an OSHA listed carcinogen.
Formaldehyde

Symptoms/Effects:

Ingestion: Practically non-toxic Large doses may cause: gastrointestinal tract irritation nausea vomiting

Inhalation: No effects anticipated

Skin Absorption: No effects anticipated

Chronic Effects: None reported

Medical Conditions Aggravated: None reported

12. ECOLOGICAL INFORMATION

Product Ecological Information: --

Based on classification principles, not classified as hazardous to the environment. No bioaccumulation potential Mobility in soil: Highly mobile

Method Used for Estimation of Aquatic Toxicity of Mixture Summation Method M-factor (Multiplier) for highly toxic ingredients: 1

Ingredient Ecological Information: Dye, Turquoise Blue: 48 hr Daphnia pulex LC50 = 100 mg/L; Formaldehyde: 96 hr Morone saxatilis LC50 = 6.7 mg/L; 96 hr Fish LC50 = 52.5 mg/L; 48 hr Daphnia pulex EC50 = 5.8 mg/L; 48 hr Daphnia magna EC50 = 29 mg/L; 48 hr Crustacea EC50 = 14 mg/L

CEPA categorization for ingredients are as follows:

Sodium Carbonate, Sodium Bicarbonate, Dye, Turquoise Blue, Water: Persistent, not bioaccumulative or inherently toxic to aquatic organisms.

Formaldehyde: Not persistent or bioaccumulative. Inherently toxic to aquatic organisms; Methanol: Not persistent, bioaccumulative or inherently toxic to aquatic organisms.

13. DISPOSAL CONSIDERATIONS

EPA Waste ID Number: Not applicable

Special Instructions (Disposal): Adjust to a pH between 6 and 9 with an acid, such as sulfuric or citric. If permitted by regulation, Open cold water tap completely, slowly pour the reacted material to the drain. Otherwise, Check with local municipal and state authorities and waste contractors for pertinent local information regarding the proper disposal of chemicals.

Empty Containers: Working in a well-ventilated area, Rinse three times with an appropriate solvent. Collect rinsate and dispose of according to local, state or federal regulations. Dispose of empty container as normal trash. In the US, rinsate from empty containers is classified as hazardous waste and should be disposed of at an E.P. A. approved facility. Rinsate from empty containers may contain sufficient product to require disposal as hazardous waste.

NOTICE (Disposal): These disposal guidelines are based on federal regulations and may be superseded by more stringent state or local requirements. Please consult your local environmental regulators for more information. In Europe: Chemical and analysis solutions must be disposed of in compliance with the respective national regulations. Product packaging must be disposed of in compliance with the country-specific regulations or must be passed to a packaging return system.

14. TRANSPORT INFORMATION

D.O.T.:

D.O.T. Proper Shipping Name: Not Currently Regulated

--

Hazard Class: NA
Subsidiary Risk: NA
ID Number: NA
Packing Group: NA

T.D.G.:

Proper Shipping Name: Not Currently Regulated

--

Hazard Class: NA
Subsidiary Risk: NA
UN Number/PIN: NA
Packing Group: NA

I.C.A.O.:

I.C.A.O. Proper Shipping Name: Not Currently Regulated

--

Hazard Class: NA
Subsidiary Risk: NA
ID Number: NA
Packing Group: NA

I.M.O.:

Proper Shipping Name: Not Currently Regulated

--

Hazard Class: NA
Subsidiary Risk: NA
ID Number: NA
Packing Group: NA
Marine Pollutant:

Additional Information: There is a possibility that this product could be contained in a reagent set or kit composed of various compatible dangerous goods. If the item is NOT in a set or kit, the classification given above applies. If the item IS part of a set or kit, the classification would change to the following: UN3316 Chemical Kit, Class 9, PG II or III. If the item is not regulated, the Chemical Kit classification does not apply.

15. REGULATORY INFORMATION

U.S. Federal Regulations:

O.S.H.A.: This product does not meet the criteria for a hazardous substance as defined in the Hazard Communication Standard. (29 CFR 1910.1200)

E.P.A.:

S.A.R.A. Title III Section 311/312 Categorization (40 CFR 370): This product is not hazardous under 29 CFR.1910.1200 and therefore is not covered by Title III under SARA.

S.A.R.A. Title III Section 313 (40 CFR 372): This product contains a chemical(s) subject to the reporting requirements of Section 313 of Title III of SARA.

Formaldehyde and Methanol

302 (EHS) TPQ (40 CFR 355): Formaldehyde 500 lbs.

304 CERCLA RQ (40 CFR 302.4): 100 lbs. Formaldehyde 5000 lbs. Methanol

304 EHS RQ (40 CFR 355): Formaldehyde - RQ 100 lbs

Clean Water Act (40 CFR 116.4): Formaldehyde - RQ 100 lbs.

RCRA: Contains RCRA regulated substances. See Section 13, EPA Waste ID Number.

State Regulations:

California Prop. 65: WARNING - This product contains a chemical known to the State of California to cause cancer.
WARNING - This product contains a chemical known to the State of California to cause birth defects or other reproductive harm.

Identification of Prop. 65 Ingredient(s): Formaldehyde and Methanol

California Perchlorate Rule CCR Title 22 Chap 33: Not applicable

Trade Secret Registry: Not applicable

National Inventories:

U.S. Inventory Status: All ingredients in this product are listed on the TSCA 8(b) Inventory (40 CFR 710).

CAS Number: Not applicable

Canadian Inventory Status: All ingredients of this product are DSL Listed.

EEC Inventory Status: All ingredients used to make this product are listed on EINECS / ELINCS.

Australian Inventory (AICS) Status: All ingredients are listed.

New Zealand Inventory (NZIoC) Status: All components either listed or exempt.

Korean Inventory (KECI) Status: All components of this product are either listed, listed as the anhydrous compound or exempt.

Japan (ENCS) Inventory Status: All components either listed or exempt.

China (PRC) Inventory (MEP) Status: All components either listed or exempt.

16. OTHER INFORMATION

References: 29 CFR 1900 - 1910 (Code of Federal Regulations - Labor). Air Contaminants, Federal Register, Vol. 54, No. 12. Thursday, January 19, 1989. pp. 2332-2983. TLV's Threshold Limit Values and Biological Exposure Indices for 1992-1993. American Conference of Governmental Industrial Hygienists, 1992. Technical Judgment. In-house information.

Complete Text of H phrases referred to in Section 3: H225 Highly flammable liquid and vapour. H301 Toxic if swallowed. Not applicable H311 Toxic in contact with skin. H314 Causes severe skin burns and eye damage. H315 Causes skin irritation. H316 Causes mild skin irritation. H317 May cause an allergic skin reaction. H319 Causes serious eye irritation. H331 Toxic if inhaled. H332 Harmful if inhaled. H334 May cause allergy or asthma symptoms or breathing difficulties if inhaled. H335 May cause respiratory irritation. H341 Suspected of causing genetic defects. H351 Suspected of causing cancer. H361d Suspected of damaging the unborn child. H370 Causes damage to organs. H412 Harmful to aquatic life with long lasting effects.

Revision Summary: . Updates in Section(s) 9,

Date of MSDS Preparation:

Day: 15

Month: December

Year: 2015

MSDS Prepared: MSDS prepared by Product Compliance Department extension 3350

CCOHS Evaluation Note: It is offered under exemption from WHMIS labeling as specified in the Controlled Products Regulation (CPR) Section 17. It is offered under the interim policy that was established by Health Canada permitting use of GHS-formatted safety data sheets in Canada prior to revision of CPR to GHS. This product has been classified and labeled in accordance with the requirements of GHS (ST/SG/AC.10/36/Add.3). This SDS has been prepared in accordance with the requirements of GHS (ST/SG/AC.10/36/Add.3).

Legend:

NA - Not Applicable	w/w - weight/weight
ND - Not Determined	w/v - weight/volume
NV - Not Available	v/v - volume/volume

USER RESPONSIBILITY: Each user should read and understand this information and incorporate it in individual site safety programs in accordance with applicable hazard communication standards and regulations.

THE INFORMATION CONTAINED HEREIN IS BASED ON DATA CONSIDERED TO BE ACCURATE. HOWEVER, NO WARRANTY IS EXPRESSED OR IMPLIED REGARDING THE ACCURACY OF THESE DATA OR THE RESULTS TO BE OBTAINED FROM THE USE THEREOF.

HACH COMPANY ©2015