

qTOWER³ and qTOWER³ touch

qPCR Thermal Cyclers

All qTOWER³ models include qTOWER³ system, Analytik Jena's PC-based qPCRsoft software to conduct qPCR and perform analysis, 115V power cable, USB cable, and user manual. The starter kit includes a color module, laptop (for non-touch models), ten 96-well plates, sealing films, and lab notebook.

- Extremely homogenous ($\pm 0.15^{\circ}\text{C}$) 96-well qPCR for up to six targets
- Rapid heating/cooling rates of $8.0^{\circ}\text{C}/\text{sec}$ and $6.0^{\circ}\text{C}/\text{sec}$, respectively
- Industry-leading 10-year optical system warranty; 2-year warranty on remainder of system (free upgrade to 3 years upon registration)

Description	Cat. No.	Unit
qTOWER ³ Real-Time PCR	75830-232	Each
qTOWER ³ G Real-Time PCR	75830-234	Each
qTOWER ³ touch Real-Time	75830-236	Each
qTOWER ³ G touch Real-Time PCR	75830-238	Each

Starter kit is bundled with the qTOWER³. No redemption necessary. This promotion is not available in Canada. Please contact your VWR Sales Representative to learn about easy access to similar options available in your region.

Promo Code | **4713**

Buy qTOWER³, Receive
Starter Kit FREE!

analytikjena
An Endress+Hauser Company

