

1. PRODUCT AND COMPANY IDENTIFICATION

Product code 95053-100, 95053-500, 95053-02K.
Product name AccuQuant™ SYBR® Green SuperMix for iQ™

Manufacturer
 QUANTA BIOSCIENCES INC.
 202 PERRY PKWY., SUITE 1
 GAITHERSBURG, MD 20877
 301-956-1670

2. COMPOSITION/INFORMATION ON INGREDIENTS

Hazardous/Non-hazardous Components

Chemical Name	CAS-No	Weight %
Dimethyl Sulfoxide	67-68-5	<10

3. HAZARDS IDENTIFICATION

Emergency Overview

Irritating to eyes. Irritating to skin. Components of the product may be absorbed into the body through the skin.

Form
 Liquid

Potential Health Hazards

Eyes Irritating to eyes.
Skin Irritating to skin. Components of the product may be absorbed into the body through the skin
Inhalation May cause irritation of respiratory tract.
Ingestion May be harmful if swallowed.

Specific Effects

Carcinogenic effects No information available
Mutagenic effects No information available
Reproductive toxicity No information available
Sensitization No information available

Target Organ Effects Eyes. Skin.

4. FIRST AID MEASURES

Skin contact Wash off immediately with plenty of soap and water. If irritation develops, get medical attention.
Eye contact Rinse thoroughly with plenty of water, also under the eyelids. If irritation persists, get medical attention.
Ingestion Never give anything by mouth to an unconscious person
Inhalation Move to fresh air
Notes to physician Treat symptomatically

5. FIRE-FIGHTING MEASURES

Suitable extinguishing media Dry chemical
Special protective equipment for firefighters Wear self-contained breathing apparatus and protective suit

6. ACCIDENTAL RELEASE MEASURES

Personal precautions Use personal protective equipment
Methods for cleaning up Soak up with inert absorbent material

7. HANDLING AND STORAGE

Handling Avoid contact with skin and eyes.
Storage Keep in properly labeled containers

8. EXPOSURE CONTROLS / PERSONAL PROTECTION

Occupational exposure controls

Exposure limits

Chemical Name	OSHA PEL (TWA)	OSHA PEL (Ceiling)	ACGIH OEL (TWA)	ACGIH OEL (STEL)
Dimethyl Sulfoxide	-	-	-	-

Engineering measures Ensure adequate ventilation, especially in confined areas

Personal protective equipment

Respiratory protection In case of insufficient ventilation, wear suitable respiratory equipment
Hand protection Protective gloves
Eye protection Safety glasses with side-shields
Skin and body protection Lightweight protective clothing
Hygiene measures Handle in accordance with good industrial hygiene and safety practice
Environmental exposure controls Prevent product from entering drains

9. PHYSICAL AND CHEMICAL PROPERTIES

General Information

Form Liquid

Important Health Safety and Environmental Information

Boiling point/range	°C 189	°F No data available
Melting point/range	°C 18.4	°F No data available
Flash point	°C 94	°F No data available
Autoignition temperature	°C No data available	°F No data available
Oxidizing properties	No information available	
Water solubility	Soluble	

10. STABILITY AND REACTIVITY

Stability Stable under normal conditions
Materials to avoid No information available
Hazardous decomposition products No information available
Polymerization Hazardous polymerization does not occur

11. TOXICOLOGICAL INFORMATION

Acute Toxicity

Chemical Name	LD50 (oral, rat/mouse)	LD50 (dermal, rat/rabbit)	LC50 (inhalation, rat/mouse)
Dimethyl Sulfoxide	14500 mg/kg (Rat)	No data available	No data available

Principle Routes of Exposure/ Potential Health Effects

Eyes	Irritating to eyes.
Skin	Irritating to skin. Components of the product may be absorbed into the body through the skin.
Inhalation	May cause irritation of respiratory tract.
Ingestion	May be harmful if swallowed.

Specific Effects

Carcinogenic effects	No information available
Mutagenic effects	No information available
Reproductive toxicity	No information available
Sensitization	No information available

Target Organ Effects	Eyes. Skin.
----------------------	-------------

12. ECOLOGICAL INFORMATION

Ecotoxicity effects	No information available
Mobility	No information available
Biodegradation	Inherently biodegradable
Bioaccumulation	Does not bioaccumulate.

13. DISPOSAL CONSIDERATIONS

Dispose of in accordance with local regulations

14. TRANSPORT INFORMATION

IATA

Proper shipping name	Not classified as dangerous in the meaning of transport regulations
Hazard Class	No information available
Subsidiary Class	No information available
Packing group	No information available
UN-No	No information available

15. REGULATORY INFORMATION

International Inventories

Chemical Name	TSCA	PICCS	ENCS	DSL	NDSL	AICS
Dimethyl sulfoxide	Listed	Listed	Listed	Listed	-	Listed

U.S. Federal Regulations

SARA 313

Not regulated

Clean Air Act, Section 112 Hazardous Air Pollutants (HAPs) (see 40 CFR 61)

This product contains the following HAPs:

U.S. State Regulations

Chemical Name	Massachusetts - RTK	New Jersey - RTK	Pennsylvania - RTK	Illinois - RTK	Rhode Island - RTK
Dimethyl sulfoxide	-	-	-	-	-

California Proposition 65

This product contains the following Proposition 65 chemicals:

WHMIS Hazard class

D2B Toxic materials

This product has been classified according to the hazard criteria of the CPR and the MSDS contains all of the information required by the CPR

16. OTHER INFORMATION

This material is sold for research and development purposes only. It is not for any human or animal therapeutic or clinical diagnostic use. It is not intended for food, drug, household, agricultural, or cosmetic use. An individual technically qualified to handle potentially hazardous chemicals must supervise the use of this material.

The above information was acquired by diligent search and/or investigation and the recommendations are based on prudent application of professional judgment. The information shall not be taken as being all inclusive and is to be used only as a guide. All materials and mixtures may be present unknown hazards and should be used with caution. Since Quanta BioSciences cannot control the actual methods, volumes, or conditions of use, the Company shall not be held liable for any damages or losses resulting from the handling or from contact with the product as described herein. THE INFORMATION IN THIS MSDS DOES NOT CONSTITUTE A WARRANTY, EXPRESS OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE.

AccuQuant is a trademark of Quanta BioSciences, Inc.

SYBR is a registered trademark of Molecular Probes, Inc.

iQ is a trademark of Bio-Rad Laboratories.

END of Safety Data Sheet