

Quality Certificate

VWR® SPINBAR® MAGNETIC STIR BAR

As per the manufacturer, the below product meets the following criteria:

North American Catalog No:	58948-986
Lot Number:	J231391
Description:	VWR BAR PTFE 2-1/2X1/2
Date of Manufacture:	January 15, 2018
Product Shelf Life:	Not Applicable
Country of Origin:	United States of America


Quality Conformance

The product is in conformance with the Purchase Order Contract, drawings, specifications, and other applicable documentation.

Product Specifications

- Material :** Virgin Dupont White Teflon #807N-X, Magnet, ALNICO 5
- BSE/TSE Statement:** BSE/TSE protein free
- Other Compliances:**
- Materials Meet FDA Regulations 21CFR177.1550, 21CFR177.1520
 - Teflon used to completely encapsulate the magnetic bar meets US Pharmacopeias Class VI
 - All materials used do not contain or come in contact with any animal bi products or organisms
 - All Items are carefully inspected to conform with ISO 9001:2008.

Signed:



Linda K. Shine
Manager - Supplier Qualification
Global Compliance
VWR International LLC

Date: May 03, 2018