


PDS No. 5410xx	PRODUCT DATA SHEET					Page 1 of 1
Revision 03	Cell Scrapers					 greiner bio-one
	Item-No. 5410xx					
Valid for Item-No.:	541070 (sterile)	541080 (sterile)				

1.	Description / Specification	
1.1	Description	Cell scraper, pivoted scraper blade, sterile 541070: short handle 541080: long handle
1.2	Dimensions	541070, -080: blade: width: 18,0 mm 541070: handle: length: 280,0 mm 541080: handle: length: 400,0 mm
1.3	Volume	-
1.4	Material / Resin	Handle: PS (Polystyrene), free of heavy metal Blade: HDPE (High Density Polyethylene), free of heavy metal
1.5	Colour	Handle: blue Blade: white
1.6	Sterilization	SAL 10 ⁻³
1.7	Quality Control	Raw Material-Control: physical testing Product-Control: testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	-
2.2	Temperature range	Handle: -20°C to +60°C Blade: -50°C to +100°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	-
2.5	Chemical Resistance	See homepage: www.gbo.com/bioscience →Products →Literature →Technical Information →Chemical Resistance of Resins
2.6	Shelf life	5 years after month of production
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	1
3.2	Pieces / Box	100
3.3	Lot-No.	F YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	Certificate of Quality

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 02	Date 15 July 2013	Date 15 July 2013	Date 15 July 2013	
Date 11.12.2009	Name S. Kaelberer	Name Dr. L. Breth	Name A. Schulz	