



VWR® Disposable Diamond-Shaped Weighing Boats

- Wide, flat bottom avoids tipping
- Temperature resistant up to +70 °C

These polystyrene, diamond shaped flexible boats feature smooth surfaces to provide accurate pour-out with minimal sample loss. They are ideal for use with microwaves and moisture analyzers.

VWR® Disposable Diamond-Shaped Weighing Boats

Capacity	W x D x H	Color	Cat. No.	Pack of
Standard				
5 mL	55 x 35 x 6 mm (2.1 x 1.3 x 0.2")	Black	10803-110	500
25 mL	80 x 50 x 14 mm (3.1 x 1.9 x 0.5")	Black	10803-128	500
30 mL	80 x 60 x 14 mm (3.1 x 2.3 x 0.5")	Black	10803-130	500
100 mL	125 x 100 x 19 mm (4.9 x 3.9 x 0.7")	Black	10803-158	500
5 mL	55 x 35 x 6 mm (2.1 x 1.3 x 0.2")	White	10803-114	500
25 mL	80 x 50 x 14 mm (3.1 x 1.9 x 0.5")	White	10803-124	500
30 mL	80 x 60 x 14 mm (3.1 x 2.3 x 0.5")	White	10803-132	500
100 mL	125 x 100 x 19 mm (4.9 x 3.9 x 0.7")	White	10803-156	500
Antistatic				
5 mL	55 x 35 x 6 mm (2.1 x 1.3 x 0.2")	Black	10803-112	500
25 mL	80 x 50 x 14 mm (3.1 x 1.9 x 0.5")	Black	10803-120	500
30 mL	80 x 60 x 14 mm (3.1 x 2.3 x 0.5")	Black	10803-122	500
100 mL	125 x 100 x 19 mm (4.9 x 3.9 x 0.7")	Black	10803-118	500
5 mL	55 x 35 x 6 mm (2.1 x 1.3 x 0.2")	White	10803-134	500
25 mL	80 x 50 x 14 mm (3.1 x 1.9 x 0.5")	White	10803-142	500
30 mL	80 x 60 x 14 mm (3.1 x 2.3 x 0.5")	White	10803-144	500
100 mL	125 x 100 x 19 mm (4.9 x 3.9 x 0.7")	White	10803-126	500