

ROTARY EVAPORATOR EVAPORATION FLASKS

Interested in increasing the efficiency and safety of your rotary evaporator? Consider pear-shaped evaporating flasks. The sloped sides of a pear-shaped flask provide an increased surface area over similar capacity spherical flasks, thereby increasing evaporation rates and efficiencies in your process. Available with and without plastic coating.

Pear Shaped Recovery Flasks

- Replacement drying flask for rotary evaporators
- Available with 30% heavier wall weight
- With or without side indents
- 24/40, 29/32, 29/42

Splash Guard Adapter

- Medium length joints that accept full length joints
- Overall length is kept to a minimum to positively effect distillation
- Available with or without fritted disc

Contact us for...

Custom glass and repair services.



ROTARY EVAPORATOR RECEIVING FLASKS

For the monitoring of collection rates, consider our graduated receiving vessel options. For smaller volumes from 0 to 100mL, our tapered bottom receiving tubes are the best choice. For less granular graduations and larger volumes from 250 to 2000mL, consider our graduated receiving bottles. Available with and without plastic coating.

Spherical Receiving Flask

- Replacement receiving flask for rotary evaporators
- Available with 30% heavier wall weight

Graduated Receiving Tube

- Tapered bottom with 35/20 spherical outer joint top
- Graduated in 0.5mL increments from 0-10mL
- Graduated in 5mL increments from 10-100mL

Graduated Receiving Bottle

- Round bottom with 35/20 spherical outer joint top
- 250mL & 500mL graduated in 10mL increments
- 1000mL & 2000mL graduated in 100mL increments
- Improved design for use with Dewar

Also Available...

Contact us for all your rotavap glassware replacement part needs, for both large and bench scale applications.



1.800.932.5000 | VWR.COM

Prices and product details are current when published and subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2019 Avantor, Inc. All rights reserved.