

Product Name:	HCUCBI0420
Description:	<p> VWR Catalog No. 10819-666 VWR Undercounter Freezer (Built-In) 4.2 Cu. Ft. Capacity One solid swing door, right hinged* Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms Remote alarm contacts 2 probes; 1 in air and 1 in glycol bottle Adjustable operating temperature range: -15°C – -25°C Two fixed shelves Keyed door lock Recessed door handle Magnetic door gasket for positive seal Leveling legs Probe access port (3/8") Powder coat exterior High density urethane foam cabinet and door insulation HFC-free refrigerant (R600a), foam insulation and packaging Two year parts and labor warranty, plus an additional three year compressor parts only warranty Exterior dimensions: 23 ¾ "W x 24 ½" D x 33 3/8 "H Manual defrost 115V, 60 Hz, 1.5 Amps 1/8 HP UL/C-UL Shipping weight: 132 lbs. *This product is available left hinged, product number HCUCBI0420LH </p>