



VWR® MAXIMUM PROTECTION NON-SKID HIGH-TOP BOOT COVERS

As per the manufacturer, the below product meets the following criteria:

VWR Catalog Number:	414004-415
Lot Number:	01192093
Description:	VWR BOOTCOVER HITOP XL CS200
Date of Manufacture:	November 25, 2019
Country of Origin:	CHINA

Quality System Compliance

Products are manufactured under the ISO 9001:2015 & ISO13485:2003 standard. Products are Inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP. Inspection records are reviewed and signed off by qualified personnel for product release.

QC Testing

Inspected according to AQL level II 4.0 sampling plan per MIL-STD-105E

-Break Tensile properties(N):	ASTM D5034standard,	MD= 108.79	CD= 111.53
-Elongation rate (%):	ASTM D5034standard,	MD= 71.34	CD= 86.85
-Trapezoid tearing strength(N):	FZ/T 01085-1999 standard,	MD= 436.00	CD= 443.00

Product Specifications

Material:	SF 63g/m ² White (breathable fibre) (Spunbonded Polypropylene fabric coated with breathable Polyethylene film)
BSE/TSE Statement:	The above product contains no ingredients of animal origin and no material derived from, or exposed to animals affected by or under quarantine for Transmitting, Animal Spongiform Encephalopathy (TSE)/Bovine Spongiform Encephalopathy(BSE)
Non-Pyrogenic:	NA
DNase & RNase Free:	NA
Cytotoxicity:	NA
Latex Free Statement:	The above product is latex free.

Signed:



Jamie Ethier
VP Global Quality
VWR, Part of Avantor

Date: January 04, 2020