

# SAFETY DATA SHEET

Issuing Date 10-Jun-2010 Revision Date 08-Jan-2016 Revision Number 2

## 1. PRODUCT AND COMPANY IDENTIFICATION

**Product Identifier** 

Product Name EDTA 0.5M Solution

Other means of identification

Product Code(s) E177

Synonyms No information available

Recommended use of the chemical and restrictions on use

Recommended Use For Further Manufacturing Use Only.
Uses advised against Not for Human or Animal Drug Use

Details of the supplier of the safety data sheet

Company AddressManufacturer AddressVWR International, LLCVWR International, LLCRadnor Corporate Center28600 Fountain Parkway100 Matsonford RoadSolon, Ohio 44139

Radnor, PA 19087-8660

Company Phone Number 1-800-448-4442 E-mail Address info@amresco-inc.com

**Emergency Telephone Number** 

Emergency Telephone Number Chemtrec 1-800-424-9300

### 2. HAZARDS IDENTIFICATION

### Classification

### **OSHA Regulatory Status**

This chemical is considered hazardous by the 2012 OSHA Hazard Communication Standard (29 CFR 1910.1200)

Skin Corrosion/Irritation	Category 1 Sub-category
Serious Eye Damage/Eye Irritation	Category 1

### Label elements

## **Emergency Overview**

### Danger

#### Hazard statements

H314 - Causes severe skin burns and eye damage



haze-free

Appearance Colorless to slightly yellow

Physical State Liquid

Odor No information available

**Precautionary Statements - Prevention** 

\_\_\_\_\_

Do not breathe dusts or mists

Wash face, hands and any exposed skin thoroughly after handling

Wear protective gloves/protective clothing/eye protection/face protection

### **Precautionary Statements - Response**

Immediately call a POISON CENTER or doctor

Specific treatment (see .? on this label)

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing Immediately call a POISON CENTER or doctor

IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/ shower

Wash contaminated clothing before reuse

IF INHALED: Remove person to fresh air and keep comfortable for breathing

Immediately call a POISON CENTER or doctor

IF SWALLOWED: Rinse mouth. DO NOT induce vomiting

### **Precautionary Statements - Storage**

Store locked up

## **Precautionary Statements - Disposal**

Dispose of contents/container to an approved waste disposal plant

### Hazards not otherwise classified (HNOC)

Not regulated

## Other Information

- May be harmful if swallowed
- · May be harmful in contact with skin
- · Harmful to aquatic life with long lasting effects

## 3. COMPOSITION/INFORMATION ON INGREDIENTS

Chemical Name	CAS-No	EC No.	Weight %	Trade Secret
EDTA, disodium dihydrate	6381-92-6	-	10-20	Not applicable
Sodium hydroxide	1310-73-2	215-185-5	5-10	Not applicable

### 4. FIRST AID MEASURES

**First Aid Measures** 

**Eye Contact**Rinse thoroughly with plenty of water for at least 15 minutes, lifting lower and upper eyelids.

Consult a physician.

**Skin Contact** Wash off immediately with soap and plenty of water removing all contaminated clothes and

shoes.

**Inhalation** Move to fresh air. If not breathing, give artificial respiration.

**Ingestion** Never give anything by mouth to an unconscious person. Rinse mouth.

Most important symptoms and effects, both acute and delayed

Most Important Symptoms/Effects No information available.

Indication of any immediate medical attention and special treatment needed

Notes to Physician Treat symptomatically.

## 5. FIRE-FIGHTING MEASURES

### **Suitable Extinguishing Media**

\_\_\_\_\_

Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

Unsuitable Extinguishing Media No information available.

### **Specific Hazards Arising from the Chemical**

No information available.

**Hazardous Combustion** 

Carbon oxides. Nitrogen oxides (NOx).

**Products** 

**Explosion Data** 

Sensitivity to Mechanical Impact None. Sensitivity to Static Discharge None.

### **Protective Equipment and Precautions for Firefighters**

As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear.

## 6. ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment and emergency procedures

**Personal protection** Ensure adequate ventilation, especially in confined areas.

**Environmental Precautions** 

**Environmental Precautions** See Section 12 for additional Ecological information.

Methods and material for containment and cleaning up

**Methods for Containment** Prevent further leakage or spillage if safe to do so.

Methods for Cleaning Up Absorb on sand or vermiculite and place in closed containers for disposal. Ventilate area

and wash spill site after material pickup is complete.

### 7. HANDLING AND STORAGE

**Precautions for Safe Handling** 

Handling Handle in accordance with good industrial hygiene and safety practice.

Conditions for safe storage, including any incompatibilities

Storage Keep containers tightly closed in a dry, cool and well-ventilated place.

**Incompatible Products** Strong oxidizing agents.

### 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

### Control parameters

Exposure Guidelines

Chemical Name	ACGIH TLV	OSHA PEL	NIOSH IDLH
EDTA, disodium dihydrate 6381-92-6	-	-	-
Sodium hydroxide 1310-73-2	-	TWA: 2 mg/m <sup>3</sup>	IDLH: 10 mg/m <sup>3</sup> Ceiling: 2 mg/m <sup>3</sup>

### **Appropriate engineering controls**

Engineering Measures Showers

Eyewash stations Ventilation systems.

### Individual protection measures, such as personal protective equipment

**Eye/Face Protection** Avoid contact with eyes. Safety glasses with side-shields.

**Skin and Body Protection** Wear protective gloves/clothing.

Respiratory Protection If exposure limits are exceeded or irritation is experienced, NIOSH/MSHA approved

respiratory protection should be worn. Positive-pressure supplied air respirators may be required for high airborne contaminant concentrations. Respiratory protection must be

provided in accordance with current local regulations.

Hygiene Measures Handle in accordance with good industrial hygiene and safety practice.

## 9. PHYSICAL AND CHEMICAL PROPERTIES

### Information on basic physical and chemical properties

Physical State Liquid

AppearanceColorless to slightly yellow haze-freeOdorNo information availableColorNo information availableOdor ThresholdNo information available

Property Values Remarks • Method

pH No information available

Melting point/freezing point

Boiling Point/Range No information available
Flash Point (High in °C) No information available
Evaporation Rate No information available
Flammability (solid, gas) No information available

Flammability Limit in Air

Upper flammability limit: No information available Lower flammability limit: No information available Vapor pressure No information available Vapor Density No information available **Specific Gravity** No information available Water Solubility No information available Solubility in other solvents No information available Partition coefficient No information available **Autoignition Temperature** No information available **Decomposition Temperature** No information available No information available Kinematic viscosity Dynamic viscosity No information available **Explosive Properties** No information available **Oxidizing Properties** No information available

### **Other Information**

Softening PointNo information availableMolecular WeightNo information availableVOC ContentNo information availableDensityNo information availableBulk DensityNo information available

## 10. STABILITY AND REACTIVITY

### Reactivity

No data available

### **Chemical Stability**

Stable under recommended storage conditions.

**Possibility of Hazardous Reactions** 

None under normal processing.

**Conditions to Avoid** 

Extremes of temperature and direct sunlight.

**Incompatible Materials** 

Strong oxidizing agents.

**Hazardous Decomposition Products** 

Carbon oxides. Nitrogen oxides (NOx).

## 11. TOXICOLOGICAL INFORMATION

### Information on likely routes of exposure

### **Product Information**

**Inhalation** There is no data available for this product.

**Eye Contact** May cause irreversible damage to eyes.

Skin Contact Causes burns.

**Ingestion** There is no data available for this product.

Chemical Name	LD50 Oral	LD50 Dermal	LC50 Inhalation
EDTA, disodium dihydrate	-	-	-
6381-92-6			
Sodium hydroxide	-	= 1350 mg/kg ( Rabbit )	-
1310-73-2		,	

### Information on toxicological effects

**Symptoms** No information available.

Delayed and immediate effects as well as chronic effects from short and long-term exposure

SensitizationNo information available.Mutagenic EffectsNo information available.CarcinogenicityNo information available.

Chemical Name	ACGIH	IARC	NTP	OSHA
EDTA, disodium dihydrate	-	-	-	-
6381-92-6				
Sodium hydroxide	-	-	-	-
1310-73-2				

Reproductive Toxicity
STOT - single exposure
STOT - repeated exposure
Aspiration hazard
No information available.
No information available.
No information available.

### Numerical measures of toxicity - Product Information

The following values are calculated based on chapter 3.1 of the GHS document .

## 12. ECOLOGICAL INFORMATION

### **Ecotoxicity**

0% of the mixture consists of components(s) of unknown hazards to the aquatic environment

Chemical Name	Algae/aquatic plants	Toxicity to Fish	Crustacea
EDTA, disodium dihydrate	-	-	-
6381-92-6			
Sodium hydroxide	-	45.4: 96 h Oncorhynchus mykiss	-
1310-73-2		mg/L LC50 static	

### Persistence and Degradability

No information available.

### Bioaccumulation/Accumulation

No information available.

Chemical Name	Partition coefficient
EDTA,disodium dihydrate 6381-92-6	-
Sodium hydroxide 1310-73-2	-

Other Adverse Effects No information available

## 13. DISPOSAL CONSIDERATIONS

### Waste treatment methods

Waste Disposal Method Dispose of material in accordance with all federal, state, and local regulations.

**Contaminated Packaging** Do not re-use empty containers.

Chemical Name	RCRA	RCRA - Basis for Listing	RCRA - D Series Wastes	RCRA - U Series Wastes
EDTA, disodium dihydrate	-	-	-	-
6381-92-6				

Chemical Name	RCRA - Halogenated Organic Compounds	RCRA - P Series Wastes	RCRA - F Series Wastes	RCRA - K Series Wastes
EDTA, disodium dihydrate 6381-92-6	-	-	-	-
Sodium hydroxide 1310-73-2	-	-	-	-

Chemical Name	California Hazardous Waste Status
EDTA,disodium dihydrate 6381-92-6	-
Sodium hydroxide 1310-73-2	-

## 14. TRANSPORT INFORMATION

<u>DOT</u>Not regulated<u>IATA</u>Not regulated

## 15. REGULATORY INFORMATION

**International Inventories** 

**TSCA** Does not comply **DSL/NDSL** Complies **EINECS/ELINCS** Does not Comply **ENCS** Does not Comply **IECSC** Complies **KECL** Does not Comply Complies **PICCS AICS** Complies

<u>Legend:</u>

TSCA - United States Toxic Substances Control Act Section 8(b) Inventory DSL/NDSL - Canadian Domestic Substances List/Non-Domestic Substances List

EINECS/ELINCS - European Inventory of Existing Commercial Chemical Substances/EU List of Notified Chemical Substances

**ENCS** - Japan Existing and New Chemical Substances

IECSC - China Inventory of Existing Chemical Substances

**KECL** - Korean Existing and Evaluated Chemical Substances

PICCS - Philippines Inventory of Chemicals and Chemical Substances

AICS - Australian Inventory of Chemical Substances

### U.S. Federal Regulations

### **SARA 313**

Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 (SARA). This product does not contain any chemicals which are subject to the reporting requirements of the Act and Title 40 of the Code of Federal Regulations, Part 372.

Component	SARA 313 - Threshold Values %
EDTA, disodium dihydrate	-
6381-92-6 ( 10-20 )	
Sodium hydroxide	-
1310-73-2 ( 5-10 )	

### SARA 311/312 Hazard Categories

Acute Health HazardYesChronic Health HazardNoFire HazardNoSudden Release of Pressure HazardNoReactive HazardNo

## **Clean Water Act**

This product does not contain any substances regulated as pollutants pursuant to the Clean Water Act (40 CFR 122.21 and 40 CFR 122.42)

Chemical Name	CWA - Reportable Quantities	CWA - Toxic Pollutants	CWA - Priority Pollutants	CWA - Hazardous Substances
EDTA, disodium dihydrate	-	-	-	-
6381-92-6				
Sodium hydroxide	-	-	-	-
1310-73-2				ļ

### CERCLA

This material, as supplied, contains one or more substances regulated as a hazardous substance under the Comprehensive Environmental Response Compensation and Liability Act (CERCLA) (40 CFR 302)

Chemical Name	Hazardous Substances RQs	Extremely Hazardous Substances RQs	RQ
EDTA, disodium dihydrate 6381-92-6	-	-	-
Sodium hydroxide	1000 lb	-	RQ 1000 lb final RQ
1310-73-2			RQ 454 kg final RQ

### U.S. State Regulations

## **California Proposition 65**

This product does not contain any Proposition 65 chemicals.

Chemical Name	California Prop. 65	
EDTA, disodium dihydrate - 6381-92-6	-	
Sodium hydroxide - 1310-73-2	-	

### U.S. State Right-to-Know Regulations

Chemical Name	New Jersey	Massachusetts	Pennsylvania
EDTA, disodium dihydrate	=	-	-
6381-92-6			
Sodium hydroxide	X	X	X
1310-73-2			

## U.S. EPA Label Information

EPA Pesticide Registration Number Not regulated

40 OTHER DIFFERENCE	
16. OTHER INFORMATION	

Issuing Date 10-Jun-2010 Revision Date 08-Jan-2016 Revision Note

No information available

**Disclaimer** 

The above information is believed to be correct but does not purport to be all-inclusive and shall be used only as a guide. The information in this document is based on the present state of our knowledge and is applicable to the product with regard to appropriate safety precautions. It does not represent any guarantee of the properties of the product. VWR International and its Affiliates shall not be held liable for any damage resulting from handling.

**End of Safety Data Sheet**