

# Certificate of Quality

<b>VWR Catalogue Number</b>	See List
<b>Description</b>	VWR® Microscope Slide Holders/Mailers
<b>Country of Origin</b>	Manufactured in China
<b>Date of Issue (yyyy-mm-dd)</b>	2022-03-24

Single-Slide Holder/Mailer

VWR Catalogue Number	Length x Width	Colour
95042-208	82 x 39 mm (3 <sup>1</sup> / <sub>4</sub> x 1 <sup>1</sup> / <sub>2</sub> "	Blue

Double-Slide Holders/Mailers

VWR Catalogue Number	Length x Width	Colour
95042-668	91 x 66 mm (3 <sup>9</sup> / <sub>16</sub> x 2 <sup>9</sup> / <sub>16</sub> "	Clear
95042-206	91 x 66 mm (3 <sup>9</sup> / <sub>16</sub> x 2 <sup>9</sup> / <sub>16</sub> "	Blue

## Quality System Compliance

Products are manufactured under the **ISO 9001:2015 & ISO 13485:2003** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

## QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

Dimensional Analysis: ASQ Z1.4-2003, General Insp. Level- II, AQL- 0.65. Result: Within tolerance  
Product is intended to hold a standard 25mm x 75mm glass slide.  
Overall dimensions are 82mm x 39.5mm x 6mm.

## Product Specifications

<b>Material</b>	100% Polypropylene, FDA 21 CFR Compliance, USP Class VI passed.
<b>Sterilization</b>	Not applicable.
<b>ATP Assay</b>	Not applicable.
<b>DNase &amp; RNase Free</b>	Not applicable.

VWR International LLC, Radnor Corporate Center, Building One, Suite 200, 100 Matson ford Road Radnor, PA 19087 - <http://www.vwr.com>

<b>BSE/TSE</b>	All materials used, and contact materials are of fully synthetic origin.
<b>Cytotoxicity</b>	Not applicable.
<b>Non-Pyrogenic Statement</b>	Not applicable.
<b>Latex Statement</b>	Not applicable.
<b>BPA Statement</b>	Not applicable.
<b>DEHP Statement</b>	Not applicable.
<b>RoHS</b>	Not applicable.
<b>REACH Statement</b>	Not applicable.
<b>Storage Conditions</b>	Climate controlled storage.
<b>Shelf Life</b>	No Expiration date.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.  
This document has been produced electronically and is valid without a signature.