



**KIMTECH SCIENCE®**  
**KIMWIPES® Ex-L Delicate Task Wipers**  
**Certificate of Conformance**

This document certifies that KIMTECH SCIENCE® KIMWIPES® Ex-L Delicate Task Wipers conform to Kimberly-Clark Corporation's internal specifications for this product as stated below. We use a system of in-process and audit inspections to assure that only product which meets specification is released for shipment to customers.

Certified by:

Comments: Lot: FL82542A

Code	Type	Color/ Colour	Material	Size	Count
34120	Pop-Up* Box	White	1 ply tissue	4.5" x 8.4" 11.4 x 21.6 cm	280 wipers/box 30 boxes/case
34133	Pop-Up* Box	White	1 ply tissue	12" x 12" 30.5 x 30.5 cm	196 wipers/box 15 boxes/ case
34155	Pop-Up* Box	White	1 ply tissue	4.5" x 8.4" 11.4 x 21.6 cm	280 wipers/box 60 boxes/case
34256	Pop-Up* Box	White	1 ply tissue	15" x 17" 38.1 x 43.2 cm	140 wipers/box 15 boxes/case

This product is 100% virgin cellulose fiber.

Physical property data is as produced, subsequent storage conditions and the age of the product, which are beyond Kimberly-Clark's control, may affect the properties of the product as received.

	Average	Test Method
<b>CD Dry Tensile Strength</b>		
g/3" strip	1,550	STP-578
<b>CD Wet Tensile Strength</b>		
g/3" strip	450	STP-578
<b>MD Dry Tensile Strength</b>		
g/3" strip	2100	STP-578
<b>Water Capacity (1 ply)</b>		
g/4" x 4"	1.3	RTM-6430
<b>Oil Capacity</b>		
g/4" x 4"	0.40	RTM-6430
<b>Gelbo Lint</b>		
#of particulates >50 micrometers	10	RTM-6176
<b>Silicone</b>		
ppm	< 5	BF3-GC

Should you require additional information, please contact our Technical Help Line at 888/346-4652.