



Quality Certificate

VWR® SHARPS CONTAINER SYSTEMS, RED

As per the manufacturer, the below product meets the following criteria:

North American Catalog No:	19001-010
Lot Number:	227513
Description:	VWR SHARP CONT XL RED WITH LIDS 8GL
Manufacture Date:	06/27/2017
Country of Origin:	USA

The product has been inspected and released in accordance with Plant Standard Production Methods and Process Quality Assurance Specifications. All resins used in the processing of these products have been approved and released. The product was manufactured in compliance with current cGMP's.

Materials of Construction: Polypropylene

Non-Pyrogenic: Product has not been treated for pyrogens.

DNase & RNase Free: Product has not been treated for DNase or RNase agents.

Cytotoxicity: N/A

Latex Free: N/A.

BSE/TSE: Resins used to mold the sharps container may contain additives from North American cattle. Any tallow used has been treated by transesterification or hydrolysis at 200 degrees C for 20 minutes at 40 bar or greater.

Signed:

Linda K. Shine
Manager - Supplier Qualification
Global Compliance
VWR International LLC

Date: 05/03/2018