

Boreal2 HM Stereo		
	470014-536	470016-036
Eyepiece	Paired 10x, 20mm	Paired 10x, 20mm
Head	Fixed 45° inclined head	Fixed 45° inclined head
Magnification	2X 4X	2X 4X
Stage	Frosted glass and reversible black & white	Frosted glass and reversible black & white
Working Distance	80mm	80mm
Focusing	Coarse focus with tension adjustment and slip clutch	Coarse focus with tension adjustment and slip clutch
Top Illumination	Cordless LED	Corded LED
Bottom Illumination	Cordless LED	Corded LED
Replacement Bulb	(T)Ward's Item #470330-976 (B) Ward's #470330-978	(T)Ward's Item #470330-976 (B) Ward's #470330-978
Product Includes	Removable soft eyecups and vinyl dust cover	Removable soft eyecups and vinyl dust cover
Dimensions	11 1/4" (H) x 7" (W) x 9" (D)	11 1/4" (H) x 7" (W) x 9" (D)
Warranty	Boreal Lifetime (includes shipping, excludes bulb & fuse)	Boreal Lifetime (includes shipping, excludes bulb & fuse)