

Amersham Anti-human IgG, Biotinylated Species- specific Whole Antibody (from sheep)

Product Specification Sheet

Code: RPN1003

Warning

For research use only.

Not recommended or intended for diagnosis of disease in humans or animals.

Do not use internally or externally in humans and animals.

Storage

Store at 2–8°C. Avoid repeated freeze-thaw cycles. Under these conditions the product is stable for six months.

Expiry

See outer packaging.

Safety warnings and precautions

All chemicals should be considered as potentially hazardous. We therefore recommend that this product is handled only by those persons who have been trained in laboratory techniques and that it is used in accordance with the principles of good laboratory practice. Wear suitable protective clothing such as laboratory overalls, safety glasses and gloves. Care should be taken to avoid contact with skin or eyes. In the case of contact with skin or eyes wash immediately with water. See material safety data sheet(s) and/or safety statement(s) for specific advice.

Component

Biotinylated antibody is supplied in Phosphate Buffered Saline (Sodium Phosphate 0.1 M, NaCl 0.1 M) pH 7.5, containing 1% (w/v) Bovine Serum Albumin and 0.05% (w/v) Sodium Azide.

Description

Purification to ensure species-specificity

The antibody is prepared by hyper-immunizing sheep with purified immunoglobulin fractions from normal mouse serum to produce high affinity antibodies. The pooled antiserum is used to produce an immunoglobulin preparation which is then affinity adsorbed to remove cross-reacting antibodies towards rat, mouse and rabbit immunoglobulins. These activities are thoroughly depleted to ensure species-specificity.

Finally, to select for specific binding, the antibodies are purified on affinity columns of human IgG. After washing to remove

any remaining non-specific serum components and low affinity antibodies, the human IgG specific antibodies are eluted using carefully selected, mild conditions which minimize aggregation and preserve immunological activity, yet which elute high affinity antibodies.

Preparation of labeled antibody

Biotin is attached to free amino groups on the purified antibodies via a spacer arm. Excess labelling reagent is removed by extensive dialysis.

The degree of biotinylation has been carefully selected to maximize binding to streptavidin, while retaining essentially complete antigen binding ability. Finally, the performance of the biotinylated antibody is evaluated and the concentration is adjusted to ensure that batch to batch variation is minimized and that the product will work in a consistently reproducible manner.

Quality control

The cross-reactivity with other species' IgG and antibody titre is determined in an ELISA system, with antigen bound to the solid phase and using a streptavidin-biotinylated peroxidase complex (RPN1051) to detect bound biotinylated antibody. These tests are carried out on each batch of the biotinylated antibody produced.

Applications

Species-specific antibodies are designed to detect only immunoglobulins from the species against which they are raised. This is desirable for any application where IgG of a second species is present, for example in the immunocytochemical staining of tissue sections and in two-site immunoassays.

Biotinylated antibodies are now used in many ways, including quantitative immunoassay, screening procedures for the detection of hybridoma antibodies, immunocytochemistry and for antigen identification in protein blotting. The biotin-streptavidin system offers flexibility in the choice of detection system in that any of the following may be used- [¹²⁵I]-streptavidin, streptavidin conjugated with fluorochromes or enzymes, streptavidin 'bridge' methods and streptavidin biotinylated enzyme complexes.

The working dilution of the biotinylated antibody will vary according to the application and this reagent has been evaluated for a number of uses. Guidelines for working titres are given below.

1. ELISA

If this reagent is to be used as a link antibody between human IgG and a streptavidin-biotinylated peroxidase complex detection system, we have found in our laboratories that a dilution of 1:40 000 is suitable for the detection of 1 µg of IgG. For greater sensitivity (for example down to 600 pg) the reagent should be diluted rather less (for example 1:5000). Thus 1.0 ml of stock reagent will be sufficient for up to 400 000 wells at the higher dilution, if used at 0.1 ml per well in standard microplates. A suitable diluent is Phosphate Buffered Saline containing 0.05% (v/v) Tween™20.



2. Immunocytochemistry

The use of this anti-human reagent for immunocytochemistry may be different from the use of other biotinylated anti-species antibodies. Usually such anti-species antibodies are used to recognize primary antibodies, for example from rats, rabbits and mice. However, human antibodies as primary reagents are not widely available, but there are applications such as the study of auto-antibody production in which this biotinylated anti-human reagent may be used. It must be accepted that background staining will arise as a result of the antibody recognizing ubiquitous IgG in human tissue sections. This must be considered by the user in making an interpretation of the results, and the reagent should be titrated to give a signal-to-noise ratio that is considered acceptable, and suitable controls should be included.

3. Estimation of final dye/protein (D/P) ratio

If this reagent is to be used to localize a detection system to antigen-antibody complex on a nitrocellulose or a nylon blot it would be expected to be sufficient for up to 30 blots of 96 cm², assuming that the membrane is incubated with a volume of 25 ml of biotinylated antibody diluted 1:400.

Related products

Mouse IgG, biotinylated whole antibody (from sheep)	RPN1001
Rabbit IgG, biotinylated whole antibody (from donkey)	RPN1004
Streptavidin alkaline phosphatase conjugate	RPN1234
Streptavidin biotinylated Horseradish Peroxidase complex	RPN1051
Streptavidin Horseradish Peroxidase conjugate	RPN1231
Streptavidin fluorescein	RPN1232
Streptavidin Texas Red™	RPN1233
ECL™ Western blotting detection reagents, 1000 cm ²	RPN2109
ECL Western blotting detection reagents, 2000 cm ²	RPN2209
ECL Western blotting detection reagents, 4000 cm ²	RPN2106
ECL Western blotting detection reagents, 6000 cm ²	RPN2134
ECL Plus Western blotting detection reagents, 1000 cm ²	RPN2132
ECL Plus Western blotting detection reagents, 3000 cm ²	RPN2133
ECL streptavidin-HRP and blocking reagent	RPN2195
ECL blocking agent	RPN2125
Mouse IgG, Horseradish Peroxidase linked whole antibody (from sheep)	NA931
Rabbit IgG, Horseradish Peroxidase linked whole antibody (from donkey)	NA934
Human IgG, Horseradish Peroxidase linked whole antibody (from sheep)	NA933
Mouse IgG, Horseradish Peroxidase linked F(Ab') ₂ fragment (from sheep)	NA9310
Rabbit IgG, Horseradish Peroxidase linked F(Ab') ₂ fragment (from donkey)	NA9340
ECL protein molecular weight markers	RPN2107
Rainbow™ colored protein molecular weight markers (molecular weight range 2350–46000)	RPN755E
Rainbow colored protein molecular weight markers (molecular weight range 14 300–200 000)	RPN756E

Full Range Rainbow recombinant protein molecular weight markers (molecular weight range 10 000–250 000)

RPN800E

Hybond™-ECL nitrocellulose membrane

RPN2020D

Legal

GE, imagination at work and GE Monogram are trademarks of General Electric Company.

Amersham, ECL, Rainbow and Hybond are trademarks of GE Healthcare companies

All third party trademarks are the property of their respective owners.

© 2002–2009 General Electric Company – All rights reserved. Previously published 2002.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.

GE Healthcare offices:

GE Healthcare Bio-Sciences AB
Björkgatan 30, 751 84 Uppsala,
Sweden

GE Healthcare Europe GmbH
Munzinger Strasse 5, D-79111 Freiburg,
Germany

GE Healthcare Bio-Sciences Corp.
800 Centennial Avenue, P.O. Box 1327,
Piscataway, NJ 08855-1327,
USA

GE Healthcare Bio-Sciences KK
Sanken Bldg. 3-25-1, Hyakunincho,
Shinjuku-ku, Tokyo 169-0073,
Japan

For your local office contact information, visit
www.gelifesciences.com/contact

GE Healthcare UK Limited
Amersham Place
Little Chalfont, Buckinghamshire,
HP7 9NA, UK

<http://www.gelifesciences.com>



imagination at work