


AC-1242xx	CUSTOMER DRAWING	Page 1 of 1
Revision 03	Cryo.s™ 5 ml, round bottom	 greiner bio-one
	Item-No. 124 2xx	



NOT TO SCALE



Cat.-No.	Description	A	B	C	D	E	F
124 261	Cryo.s™ 5 ml, natural	Ø 10.7	Ø 12.5	Ø 12.6	76.0	85.85	17.0
124 263	Cryo.s™ 5 ml, natural	Ø 10.7	Ø 12.5	Ø 12.6	76.0	85.85	17.0
124 273	Cryo.s™ 5 ml, green	Ø 10.7	Ø 12.5	Ø 12.6	76.0	85.85	17.0
124 274	Cryo.s™ 5 ml, yellow	Ø 10.7	Ø 12.5	Ø 12.6	76.0	85.85	17.0
124 275	Cryo.s™ 5 ml, blue	Ø 10.7	Ø 12.5	Ø 12.6	76.0	85.85	17.0
124 276	Cryo.s™ 5 ml, red	Ø 10.7	Ø 12.5	Ø 12.6	76.0	85.85	17.0

All dimensions in mm

Customer drawing subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 02	Date November 20, 2008	Date November 20, 2008	Date November 21, 2000	
Date November 17, 2006	Name F. Falter	Name F. Sonntag	Name F. Sonntag	