

Certificate Of Analysis

Whatman™

Product Information

Product Number: 1851-025

Product Name: QMA 2.5CM 100/PK

Lot Number/Serial Number: 17238127

Date of Manufacture: 14-Dec-2020

Furnish Number: GS20193

Quality Control

Test	Range	Results
Grammage,gsm	73-93	85
Thickness, um@53kPa	400-550	496
Porosity, S/300ml/1.0sq.in	14-24	22
Air Retention,% 0.3µm@5.33cm/sec	≥99.95	99.99

Conformance & Quality systems statement

This is to certify that this product conforms to GE Healthcare Life Sciences specifications.

All products are manufactured via a Quality System certified to ISO9001:2015 and tested in accordance with documented quality procedures and approved as a result of meeting the required specification.

Electronic signature

This document has electronically produced and is valid without a signature.

Version AH



www.gelifesciences.com

Global Life Sciences Solutions Operations UK Ltd.

Amersham Place Little Chalfont

Buckinghamshire HP7 9NA UK