

The following list contains the Material Safety Data Sheets you requested. Please scroll down to view the requested MSDS(s).

Product	MSDS	Distributor	Format	Language	Quantity
148002	70532	Hach Company	OSHA	English	1
148002	105599	Hach Company	OSHA	English	1
148002	220399	Hach Company	OSHA	English	1

---

Total Enclosures: 3

World Headquarters  
Hach Company  
P.O.Box 389  
Loveland, CO USA 80539  
(970) 669-3050

MSDS No: M00740

# MATERIAL SAFETY DATA SHEET

---

## 1. CHEMICAL PRODUCT AND COMPANY IDENTIFICATION

**Product Name:** Sulfite 3 Reagent  
**Catalog Number:** 70532

Hach Company  
P.O.Box 389  
Loveland, CO USA 80539  
(970) 669-3050

Emergency Telephone Numbers:  
(Medical and Transportation)  
(303) 623-5716 24 Hour Service  
(515)232-2533 8am - 4pm CST

**MSDS Number:** M00740  
**Chemical Name:** Not applicable  
**CAS No.:** Not applicable  
**Chemical Formula:** Not applicable  
**Chemical Family:** Not applicable  
**Hazard:** Practically non-toxic.  
**Date of MSDS Preparation:**  
**Day:** 06  
**Month:** May  
**Year:** 2013

---

## 2. COMPOSITION / INFORMATION ON INGREDIENTS

### Demineralized Water

**CAS No.:** 7732-18-5  
**TSCA CAS Number:** 7732-18-5  
**Percent Range:** > 99.0  
**Percent Range Units:** volume / volume  
**LD50:** LD50 oral rat = >89,800 mg/kg  
**LC50:** None reported  
**TLV:** Not established  
**PEL:** Not established  
**Hazard:** No effects anticipated.

### Potassium Iodate

**CAS No.:** 7758-05-6  
**TSCA CAS Number:** 7758-05-6  
**Percent Range:** < 1.0  
**Percent Range Units:** weight / volume  
**LD50:** Oral mouse LDLo = 531 mg/kg  
**LC50:** None reported  
**TLV:** Not established  
**PEL:** Not established  
**Hazard:** Oxidizer.May cause irritation.

---

## 3. HAZARDS IDENTIFICATION

**Emergency Overview:**  
**Appearance:** Clear, colorless liquid  
**Odor:** None  
**HMIS:**

**Health:** 0  
**Flammability:** 0  
**Reactivity:** 0  
**Protective Equipment:** X - See protective equipment, Section 8.

**NFPA:**

**Health:** 0  
**Flammability:** 0  
**Reactivity:** 0  
**Symbol:** Not applicable

**Potential Health Effects:**

**Eye Contact:** No effects are anticipated  
**Skin Contact:** No effects are anticipated  
**Skin Absorption:** No effects anticipated  
**Target Organs:** Not applicable  
**Ingestion:** No Effects Anticipated  
**Target Organs:** Not applicable  
**Inhalation:** No effects anticipated  
**Target Organs:** Not applicable  
**Medical Conditions Aggravated:** None reported  
**Chronic Effects:** No effects anticipated  
**Cancer / Reproductive Toxicity Information:**  
This product does NOT contain any OSHA listed carcinogens.  
This product does NOT contain any IARC listed chemicals.  
This product does NOT contain any NTP listed chemicals.  
**Additional Cancer / Reproductive Toxicity Information:** None reported  
**Toxicologically Synergistic Products:** None reported

---

## 4. FIRST AID

**Eye Contact:** Flush eyes with water.  
**Skin Contact (First Aid):** Wash skin with plenty of water.  
**Ingestion (First Aid):** Give large quantities of water. Call physician immediately.  
**Inhalation:** None required.

---

## 5. FIRE FIGHTING MEASURES

**Flammable Properties:** Material will not burn.  
**Flash Point:** Not applicable  
**Method:** Not applicable  
**Flammability Limits:**  
**Lower Explosion Limits:** Not applicable  
**Upper Explosion Limits:** Not applicable  
**Autoignition Temperature:** Not applicable  
**Hazardous Combustion Products:** None  
**Fire / Explosion Hazards:** Not applicable  
**Static Discharge:** None reported.  
**Mechanical Impact:** None reported  
**Extinguishing Media:** Use media appropriate to surrounding fire conditions  
**Fire Fighting Instruction:** As in any fire, wear self-contained breathing apparatus pressure-demand and full protective gear.

---

## 6. ACCIDENTAL RELEASE MEASURES

**Spill Response Notice:**  
Only persons properly qualified to respond to an emergency involving hazardous substances may respond to a spill according to federal regulations (OSHA 29 CFR 1910.120(a)(v)) and per your company's emergency response plan and guidelines/procedures. See Section 13, Special Instructions for disposal assistance.  
**Containment Technique:** Stop spilled material from being released to the environment.  
**Clean-up Technique:** Absorb spilled liquid with non-reactive sorbent material. Place material in a plastic bag. Mark bag 'Non-hazardous trash', and dispose of as normal refuse.

**Evacuation Procedure:** Evacuate as needed to perform spill clean-up. If conditions warrant, increase the size of the evacuation.

**Special Instructions (for accidental release):** Not applicable

**304 EHS RQ (40 CFR 355):** Not applicable

**D.O.T. Emergency Response Guide Number:** Not applicable

---

## 7. HANDLING / STORAGE

**Handling:** Avoid contact with eyes Wash thoroughly after handling.

**Storage:** Keep container tightly closed when not in use.

**Flammability Class:** Not applicable

---

## 8. EXPOSURE CONTROLS / PROTECTIVE EQUIPMENT

**Engineering Controls:** Maintain general industrial hygiene practices when using this product.

**Personal Protective Equipment:**

**Eye Protection:** safety glasses with top and side shields

**Skin Protection:** lab coat

**Inhalation Protection:** adequate ventilation

**Precautionary Measures:** Avoid contact with: eyes Wash thoroughly after handling.

**TLV:** Not established

**PEL:** Not established

---

## 9. PHYSICAL / CHEMICAL PROPERTIES

**Appearance:** Clear, colorless liquid

**Physical State:** Liquid

**Molecular Weight:** Not applicable

**Odor:** None

**pH:** 5.8

**Vapor Pressure:** Not determined

**Vapor Density (air = 1):** Not determined

**Boiling Point:** ~100°C

**Melting Point:** Not applicable

**Specific Gravity/ Relative Density (water = 1; air =1):** 0.983

**Evaporation Rate (water = 1):** 0.84

**Volatile Organic Compounds Content:** Not applicable

**Partition Coefficient (n-octanol / water):** Not applicable

**Solubility:**

**Water:** Miscible

**Acid:** Miscible

**Other:** Not determined

**Metal Corrosivity:**

**Steel:** None

**Aluminum:** None

---

## 10. STABILITY / REACTIVITY

**Chemical Stability:** Stable when stored under proper conditions.

**Conditions to Avoid:** Heat Evaporation

**Reactivity / Incompatibility:** None reported

**Hazardous Decomposition:** None reported

**Hazardous Polymerization:** Will not occur.

---

## 11. TOXICOLOGICAL INFORMATION

**Product Toxicological Data:**

**LD50:** None reported

**LC50:** None reported

**Dermal Toxicity Data:** None reported

*Skin and Eye Irritation Data:* None reported  
*Mutation Data:* None reported  
*Reproductive Effects Data:* None reported  
*Ingredient Toxicological Data:* Potassium Iodate: Oral mouse LDLo = 531 mg/kg

---

## 12. ECOLOGICAL INFORMATION

*Product Ecological Information:* --  
No ecological data available for this product.  
*Ingredient Ecological Information:* --  
No ecological data available for the ingredients of this product.  
CEPA Statement: Potassium Iodate: Persistent, not bioaccumulative or inherently toxic to aquatic organisms.

---

## 13. DISPOSAL CONSIDERATIONS

*EPA Waste ID Number:* Not applicable  
*Special Instructions (Disposal):* Open cold water tap completely, slowly pour the material to the drain. Flush system with plenty of water.  
*Empty Containers:* Rinse three times with an appropriate solvent. Dispose of empty container as normal trash.  
*NOTICE (Disposal):* These disposal guidelines are based on federal regulations and may be superseded by more stringent state or local requirements. Please consult your local environmental regulators for more information.

---

## 14. TRANSPORT INFORMATION

*D.O.T.:*  
*D.O.T. Proper Shipping Name:* Not Currently Regulated  
--  
*DOT Hazard Class:* NA  
*DOT Subsidiary Risk:* NA  
*DOT ID Number:* NA  
*DOT Packing Group:* NA  
*I.C.A.O.:*  
*I.C.A.O. Proper Shipping Name:* Not Currently Regulated  
--  
*ICAO Hazard Class:* NA  
*ICAO Subsidiary Risk:* NA  
*ICAO ID Number:* NA  
*ICAO Packing Group:* NA  
*I.M.O.:*  
*I.M.O. Proper Shipping Name:* Not Currently Regulated  
--  
*I.M.O. Hazard Class:* NA  
*I.M.O. Subsidiary Risk:* NA  
*I.M.O. ID Number:* NA  
*I.M.O. Packing Group:* NA

*Additional Information:* There is a possibility that this product could be contained in a reagent set or kit composed of various compatible dangerous goods. If the item is NOT in a set or kit, the classification given above applies. If the item IS part of a set or kit, the classification would change to the following: UN3316 Chemical Kit, Class 9, PG II or III. If the item is not regulated, the Chemical Kit classification does not apply.

---

## 15. REGULATORY INFORMATION

*U.S. Federal Regulations:*  
*O.S.H.A.:* This product does not meet the criteria for a hazardous substance as defined in the Hazard Communication Standard. (29 CFR 1910.1200)  
*E.P.A.:*  
*S.A.R.A. Title III Section 311/312 Categorization (40 CFR 370):* This product is not hazardous under 29 CFR.1910.1200 and therefore is not covered by Title III under SARA.  
*S.A.R.A. Title III Section 313 (40 CFR 372):* This product does NOT contain any chemical subject to the reporting requirements of Section 313 of Title III of SARA.

*302 (EHS) TPQ (40 CFR 355):* Not applicable  
*304 CERCLA RQ (40 CFR 302.4):* Not applicable  
*304 EHS RQ (40 CFR 355):* Not applicable  
*Clean Water Act (40 CFR 116.4):* Not applicable  
*RCRA:* Contains no RCRA regulated substances.

*C.P.S.C.:* Not applicable

**State Regulations:**

*California Prop. 65:* No Prop. 65 listed chemicals are present in this product.

*Identification of Prop. 65 Ingredient(s):* Not applicable

*California Perchlorate Rule CCR Title 22 Chap 33:* No

*Trade Secret Registry:* Not applicable

**National Inventories:**

*U.S. Inventory Status:* All ingredients in this product are listed on the TSCA 8(b) Inventory (40 CFR 710).

*TSCA CAS Number:* Not applicable

---

## 16. OTHER INFORMATION

**Intended Use:** Sulfite determination

**References:** 29 CFR 1900 - 1910 (Code of Federal Regulations - Labor). Air Contaminants, Federal Register, Vol. 54, No. 12. Thursday, January 19, 1989. pp. 2332-2983. TLV's Threshold Limit Values and Biological Exposure Indices for 1992-1993. American Conference of Governmental Industrial Hygienists, 1992. Technical Judgment. In-house information. Fire Protection Guide on Hazardous Materials, 10th Ed. Quincy, MA: National Fire Protection Fire Protection Guide on Hazardous Materials, 10th Ed. Quincy, MA: National Fire Protection Association, 1991.

**Revision Summary:** Substantial revision to comply with EU Reg 1272/2008, Reg 1907/2006 and UN GHS (ST/SG/AC.10/36/Add.3).

---

**Legend:**

NA - Not Applicable	w/w - weight/weight
ND - Not Determined	w/v - weight/volume
NV - Not Available	v/v - volume/volume

**USER RESPONSIBILITY:** Each user should read and understand this information and incorporate it in individual site safety programs in accordance with applicable hazard communication standards and regulations.

**THE INFORMATION CONTAINED HEREIN IS BASED ON DATA CONSIDERED TO BE ACCURATE. HOWEVER, NO WARRANTY IS EXPRESSED OR IMPLIED REGARDING THE ACCURACY OF THESE DATA OR THE RESULTS TO BE OBTAINED FROM THE USE THEREOF.**

**HACH COMPANY ©2013**

World Headquarters  
Hach Company  
P.O.Box 389  
Loveland, CO USA 80539  
(970) 669-3050

MSDS No: M00007

# MATERIAL SAFETY DATA SHEET

---

## 1. CHEMICAL PRODUCT AND COMPANY IDENTIFICATION

**Product Name:** Sulfamic Acid  
**Catalog Number:** 105599

Hach Company  
P.O.Box 389  
Loveland, CO USA 80539  
(970) 669-3050

Emergency Telephone Numbers:  
(Medical and Transportation)  
(303) 623-5716 24 Hour Service  
(515)232-2533 8am - 4pm CST

**MSDS Number:** M00007  
**Chemical Name:** Sulfamic Acid  
**CAS No.:** 5329-14-6  
**Chemical Formula:**  $\text{H}_3\text{NO}_3\text{S}$   
**Chemical Family:** Inorganic Acid  
**Hazard:** Causes eye burns.  
**Date of MSDS Preparation:**  
**Day:** 25  
**Month:** 06  
**Year:** 2010

---

## 2. COMPOSITION / INFORMATION ON INGREDIENTS

### Other component

**CAS No.:** Not applicable  
**TSCA CAS Number:** Not applicable  
**Percent Range:** < 1.0  
**Percent Range Units:** weight / weight  
**LD50:** Not applicable  
**LC50:** Not applicable  
**TLV:** Not established  
**PEL:** Not established  
**Hazard:** Any ingredient(s) of this product listed as "Other component(s)" is not considered a health hazard to the user of this product.

### Sulfamic Acid

**CAS No.:** 5329-14-6  
**TSCA CAS Number:** 5329-14-6  
**Percent Range:** > 99.0  
**Percent Range Units:** weight / weight  
**LD50:** Oral rat LD50 = 3160 mg/kg.  
**LC50:** None reported  
**TLV:** Not established  
**PEL:** Not established  
**Hazard:** Causes eye burns.

---

## 3. HAZARDS IDENTIFICATION

**Emergency Overview:**  
**Appearance:** White crystals  
**Odor:** None

CAUSES EYE BURNS CAUSES SKIN AND RESPIRATORY TRACT IRRITATION

**HMIS:**

**Health:** 2

**Flammability:** 1

**Reactivity:** 1

**Protective Equipment:** X - See protective equipment, Section 8.

**NFPA:**

**Health:** 2

**Flammability:** 1

**Reactivity:** 1

**Symbol:** Not applicable

**Potential Health Effects:**

**Eye Contact:** Causes eye burns.

**Skin Contact:** Causes severe irritation

**Skin Absorption:** None reported

**Target Organs:** None reported

**Ingestion:** May cause: irritation of the mouth and esophagus gastrointestinal irritation

**Target Organs:** None reported

**Inhalation:** May cause: irritation of nose and throat

**Target Organs:** None reported

**Medical Conditions Aggravated:** Pre-existing: Eye conditions Skin conditions Respiratory conditions

**Chronic Effects:** None reported

**Cancer / Reproductive Toxicity Information:**

O.S.H.A. Listed: No

IARC Listed: No

NTP Listed: No

**Additional Cancer / Reproductive Toxicity Information:** Not applicable

**Toxicologically Synergistic Products:** None reported

---

## 4. FIRST AID

**Eye Contact:** Immediately flush eyes with water for 15 minutes. Call physician.

**Skin Contact (First Aid):** Wash skin with plenty of water for 15 minutes. Call physician immediately.

**Ingestion (First Aid):** Do not induce vomiting. Give 1-2 glasses of water. Call physician immediately. Never give anything by mouth to an unconscious person.

**Inhalation:** Remove to fresh air.

---

## 5. FIRE FIGHTING MEASURES

**Flammable Properties:** During a fire, irritating and highly toxic gases may be generated by thermal decomposition.

**Flash Point:** Not applicable

**Method:** Not applicable

**Flammability Limits:**

**Lower Explosion Limits:** Not applicable

**Upper Explosion Limits:** Not applicable

**Autoignition Temperature:** Not applicable

**Hazardous Combustion Products:** Toxic fumes of: ammonia nitrogen oxides. sulfur oxides.

**Fire / Explosion Hazards:** May react violently with: chlorine / chlorine compounds metal nitrates metal nitrites nitric acid

**Static Discharge:** None reported.

**Mechanical Impact:** None reported

**Extinguishing Media:** Dry chemical. Water.

**Fire Fighting Instruction:** As in any fire, wear self-contained breathing apparatus pressure-demand and full protective gear.

---

## 6. ACCIDENTAL RELEASE MEASURES

**Spill Response Notice:**

Only persons properly qualified to respond to an emergency involving hazardous substances may respond to a spill according to federal regulations (OSHA 29 CFR 1910.120(a)(v)) and per your company's emergency response plan and guidelines/procedures. See Section 13, Special Instructions for disposal assistance.



**Containment Technique:** Stop spilled material from being released to the environment. Cover spilled solid material with sand or other inert material.

**Clean-up Technique:** Scoop up spilled material into a large beaker and dissolve with water. Adjust to a pH between 6 and 9 with an alkali, such as soda ash or sodium bicarbonate. Flush reacted material to the drain with a large excess of water. Decontaminate the area of the spill with a soap solution.

**Evacuation Procedure:** Evacuate local area (15 foot radius or as directed by your facility's emergency response plan) when: a pound or more of loose powder is spilled. If conditions warrant, increase the size of the evacuation.

**Special Instructions (for accidental release):** Product is regulated as RCRA hazardous waste in the U.S.

**304 EHS RQ (40 CFR 355):** Not applicable

**D.O.T. Emergency Response Guide Number:** 154

---

## 7. HANDLING / STORAGE

**Handling:** Avoid contact with eyes skin Do not breathe dust. Maintain general industrial hygiene practices when using this product.

**Storage:** Store away from: oxidizers alkalies chlorine/chlorinated metals Protect from: heat moisture

**Flammability Class:** Not applicable

---

## 8. EXPOSURE CONTROLS / PROTECTIVE EQUIPMENT

**Engineering Controls:** Have an eyewash station nearby. Maintain general industrial hygiene practices when using this product.

**Personal Protective Equipment:**

**Eye Protection:** safety glasses with top and side shields

**Skin Protection:** disposable latex gloves lab coat

**Inhalation Protection:** adequate ventilation

**Precautionary Measures:** Avoid contact with: eyes skin Do not breathe: dust Wash thoroughly after handling. Keep away from: alkalies metals Protect from: heat moisture

**TLV:** Not established

**PEL:** Not established

---

## 9. PHYSICAL / CHEMICAL PROPERTIES

**Appearance:** White crystals

**Physical State:** Solid

**Molecular Weight:** 97.10

**Odor:** None

**pH:** 1% soln = 1.18

**Vapor Pressure:** Not applicable

**Vapor Density (air = 1):** Not applicable

**Boiling Point:** Not applicable

**Melting Point:** Product decomposes at 205 °C; 401 °F

**Specific Gravity/ Relative Density (water = 1; air =1):** 2.15

**Evaporation Rate (water = 1):** Not applicable

**Volatile Organic Compounds Content:** Not applicable

**Partition Coefficient (n-octanol / water):** None reported

**Solubility:**

**Water:** 1:2 ratio @ 80 ° C (176 °F)

**Acid:** Soluble

**Other:** Slightly soluble in alcohol, methanol.

**Metal Corrosivity:**

**Steel:** 0.814 in/yr

**Aluminum:** 0.212 in/yr

---

## 10. STABILITY / REACTIVITY

**Chemical Stability:** Stable when stored under proper conditions.

**Conditions to Avoid:** Heating to decomposition. Excess moisture

**Reactivity / Incompatibility:** May react violently in contact with: chlorates metal nitrates metal nitrites nitric acid

Incompatible with: alkalies oxidizers

**Hazardous Decomposition:** Heating to decomposition releases toxic and/or corrosive fumes of: ammonia nitrogen oxides sulfur oxides

**Hazardous Polymerization:** Will not occur.

---

## 11. TOXICOLOGICAL INFORMATION

**Product Toxicological Data:**

**LD50:** Oral rat LD50 = 3160 mg/kg.

**LC50:** None reported

**Dermal Toxicity Data:** None reported

**Skin and Eye Irritation Data:** Skin Human 4%/5 days intermittent MILD, Skin rabbit 500 mg/24H SEVERE, Eye rabbit 20mg MODERATE, Eye rabbit 250µg/24H SEVERE.

**Mutation Data:** None reported

**Reproductive Effects Data:** None reported

**Ingredient Toxicological Data:** --

Not applicable

---

## 12. ECOLOGICAL INFORMATION

**Product Ecological Information:** --

No ecological data available for this product.

**Ingredient Ecological Information:** --

Not applicable

---

## 13. DISPOSAL CONSIDERATIONS

**EPA Waste ID Number:** None

**Special Instructions (Disposal):** Dilute to 3 to 5 times the volume with cold water. Adjust to a pH between 6 and 9 with an alkali, such as soda ash or sodium bicarbonate. Open cold water tap completely, slowly pour the reacted material to the drain.

**Empty Containers:** Rinse three times with an appropriate solvent. Dispose of empty container as normal trash.

**NOTICE (Disposal):** These disposal guidelines are based on federal regulations and may be superseded by more stringent state or local requirements. Please consult your local environmental regulators for more information.

---

## 14. TRANSPORT INFORMATION

**D.O.T.:**

**D.O.T. Proper Shipping Name:** Sulphamic Acid

--

**DOT Hazard Class:** 8

**DOT Subsidiary Risk:** NA

**DOT ID Number:** UN2967

**DOT Packing Group:** III

**I.C.A.O.:**

**I.C.A.O. Proper Shipping Name:** Sulphamic Acid

--

**ICAO Hazard Class:** 8

**ICAO Subsidiary Risk:** NA

**ICAO ID Number:** UN2967

**ICAO Packing Group:** III

**I.M.O.:**

**I.M.O. Proper Shipping Name:** Sulphamic Acid

--

**I.M.O. Hazard Class:** 8

**I.M.O. Subsidiary Risk:** NA

**I.M.O. ID Number:** UN2967

**I.M.O. Packing Group:** III

**Additional Information:** There is a possibility that this product could be contained in a reagent set or kit composed of various compatible dangerous goods. If the item is NOT in a set or kit, the classification given above applies. If the item IS

part of a set or kit, the classification would change to the following: UN3316 Chemical Kit, Class 9, PG II or III. If the item is not regulated, the Chemical Kit classification does not apply.

---

## 15. REGULATORY INFORMATION

### *U.S. Federal Regulations:*

**O.S.H.A.:** This product meets the criteria for a hazardous substance as defined in the Hazard Communication Standard. (29 CFR 1910.1200)

### *E.P.A.:*

**S.A.R.A. Title III Section 311/312 Categorization (40 CFR 370):** Immediate (Acute) Health Hazard

**S.A.R.A. Title III Section 313 (40 CFR 372):** This product does NOT contain any chemical subject to the reporting requirements of Section 313 of Title III of SARA.

--

**302 (EHS) TPQ (40 CFR 355):** Not applicable

**304 CERCLA RQ (40 CFR 302.4):** Not applicable

**304 EHS RQ (40 CFR 355):** Not applicable

**Clean Water Act (40 CFR 116.4):** Not applicable

**RCRA:** Contains no RCRA regulated substances.

**C.P.S.C.:** Not applicable

### *State Regulations:*

**California Prop. 65:** No Prop. 65 listed chemicals are present in this product.

**Identification of Prop. 65 Ingredient(s):** Not applicable

**California Perchlorate Rule CCR Title 22 Chap 33:**

**Trade Secret Registry:** Not applicable

### *National Inventories:*

**U.S. Inventory Status:** TSCA Listed: Yes

**TSCA CAS Number:** 5329-14-6

---

## 16. OTHER INFORMATION

**Intended Use:** Laboratory Use

**References:** Vendor Information. NIOSH Registry of Toxic Effects of Chemical Substances, 1985-86. Cincinnati: U.S. Department of Health and Human Services, April, 1987. Gosselin, R. E. et al. Clinical Toxicology of Commercial Products, 5th Ed. Baltimore: The Williams and Wilkins Co., 1984. Fire Protection Guide on Hazardous Materials, 10th Ed. Quincy, MA: National Fire Protection Fire Protection Guide on Hazardous Materials, 10th Ed. Quincy, MA: National Fire Protection Association, 1991. Outside Testing. Technical Judgment. The Merck Index, 11th Ed. Rahway, New Jersey: Merck and Co., Inc., 1989. Sax, N. Irving. Dangerous Properties of Industrial Materials, 7th Ed. New York: Van Nostrand Reinhold Co., 1989. Air Contaminants, Federal Register, Vol. 54, No. 12. Thursday, January 19, 1989. pp. 2332-2983. TLV's Threshold Limit Values and Biological Exposure Indices for 1992-1993. American Conference of Governmental Industrial Hygienists, 1992.

**Revision Summary:** Updates in Section(s) 14,

---

### **Legend:**

NA - Not Applicable	w/w - weight/weight
ND - Not Determined	w/v - weight/volume
NV - Not Available	v/v - volume/volume

**USER RESPONSIBILITY:** Each user should read and understand this information and incorporate it in individual site safety programs in accordance with applicable hazard communication standards and regulations.

**THE INFORMATION CONTAINED HEREIN IS BASED ON DATA CONSIDERED TO BE ACCURATE. HOWEVER, NO WARRANTY IS EXPRESSED OR IMPLIED REGARDING THE ACCURACY OF THESE DATA OR THE RESULTS TO BE OBTAINED FROM THE USE THEREOF.**

**HACH COMPANY ©2010**

World Headquarters  
Hach Company  
P.O.Box 389  
Loveland, CO USA 80539  
(970) 669-3050

MSDS No: M00011

# MATERIAL SAFETY DATA SHEET

---

## 1. CHEMICAL PRODUCT AND COMPANY IDENTIFICATION

**Product Name:** Sulfite 1 Reagent  
**Catalog Number:** 220399

Hach Company  
P.O.Box 389  
Loveland, CO USA 80539  
(970) 669-3050

Emergency Telephone Numbers:  
(Medical and Transportation)  
(303) 623-5716 24 Hour Service  
(515)232-2533 8am - 4pm CST

**MSDS Number:** M00011  
**Chemical Name:** Not applicable  
**CAS No.:** Not applicable  
**Chemical Formula:** Not applicable  
**Chemical Family:** Not applicable  
**Hazard:** Causes irritation.  
**Date of MSDS Preparation:**  
**Day:** 20  
**Month:** April  
**Year:** 2012

---

## 2. COMPOSITION / INFORMATION ON INGREDIENTS

### Potassium Iodide

**CAS No.:** 7681-11-0  
**TSCA CAS Number:** 7681-11-0  
**Percent Range:** >95.0  
**Percent Range Units:** weight / weight  
**LD50:** Oral Mouse LD50 = 1862 mg/kg  
**LC50:** None reported  
**TLV:** Not established  
**PEL:** Not established  
**Hazard:** Causes irritation.

### Other component

**CAS No.:** Not applicable  
**TSCA CAS Number:** Not applicable  
**Percent Range:** < 1.0  
**Percent Range Units:** weight / weight  
**LD50:** Not applicable  
**LC50:** Not applicable  
**TLV:** Not established  
**PEL:** Not established  
**Hazard:** Any ingredient(s) of this product listed as "Other component(s)" is not considered a health hazard to the user of this product.

### Potassium Carbonate

**CAS No.:** 584-08-7  
**TSCA CAS Number:** 584-08-7  
**Percent Range:** < 5.0  
**Percent Range Units:** weight / weight  
**LD50:** Oral rat LD<sub>50</sub> = 1870 mg/kg

**LC50:** None reported  
**TLV:** Not established  
**PEL:** Not established  
**Hazard:** Causes irritation.

#### **Starch**

**CAS No.:** 9005-84-9  
**TSCA CAS Number:** 9005-84-9  
**Percent Range:** < 5.0  
**Percent Range Units:** weight / weight  
**LD50:** None reported  
**LC50:** None reported  
**TLV:** 10 mg/m<sup>3</sup> total dust  
**PEL:** 15 mg/m<sup>3</sup> total dust  
**Hazard:** May cause allergic reaction. May cause irritation.

---

### **3. HAZARDS IDENTIFICATION**

#### ***Emergency Overview:***

**Appearance:** White to tan powder or crystals

**Odor:** None

CAUSES EYE IRRITATION MAY CAUSE SKIN AND RESPIRATORY TRACT IRRITATION

#### ***HMIS:***

**Health:** 2

**Flammability:** 0

**Reactivity:** 0

**Protective Equipment:** X - See protective equipment, Section 8.

#### ***NFPA:***

**Health:** 1

**Flammability:** 0

**Reactivity:** 0

**Symbol:** Not applicable

#### ***Potential Health Effects:***

**Eye Contact:** Causes irritation

**Skin Contact:** May cause irritation

**Skin Absorption:** None reported

**Target Organs:** None reported

**Ingestion:** May cause: abdominal pain May cause iodism, which symptoms include skin rash, conjunctivitis, runny nose, sneezing, bronchitis, headache, fever and irritation of mucous membranes.

**Target Organs:** None reported

**Inhalation:** May cause: respiratory tract irritation

**Target Organs:** None reported

**Medical Conditions Aggravated:** Pre-existing: Skin conditions

**Chronic Effects:** None reported

#### ***Cancer / Reproductive Toxicity Information:***

This product does NOT contain any OSHA listed carcinogens.

This product does NOT contain any IARC listed chemicals.

This product does NOT contain any NTP listed chemicals.

**Additional Cancer / Reproductive Toxicity Information:** None reported

**Toxicologically Synergistic Products:** None reported

---

### **4. FIRST AID**

**Eye Contact:** Immediately flush eyes with water for 15 minutes. Call physician.

**Skin Contact (First Aid):** Wash skin with plenty of water.

**Ingestion (First Aid):** Give large quantities of water. Call physician immediately.

**Inhalation:** Remove to fresh air. Give artificial respiration if necessary. Call physician.

---

### **5. FIRE FIGHTING MEASURES**

**Flammable Properties:** Does not burn, but may melt in a fire, releasing toxic fumes.

**Flash Point:** Not applicable

**Method:** Not applicable

**Flammability Limits:**

**Lower Explosion Limits:** Not applicable

**Upper Explosion Limits:** Not applicable

**Autoignition Temperature:** Not applicable

**Hazardous Combustion Products:** Toxic fumes of: potassium oxides iodine iodine compounds

**Fire / Explosion Hazards:** May react violently with: bromine trifluoride

**Static Discharge:** None reported.

**Mechanical Impact:** None reported

**Extinguishing Media:** Dry chemical. Carbon dioxide Alcohol foam.

**Fire Fighting Instruction:** As in any fire, wear self-contained breathing apparatus pressure-demand and full protective gear.

---

## 6. ACCIDENTAL RELEASE MEASURES

**Spill Response Notice:**

Only persons properly qualified to respond to an emergency involving hazardous substances may respond to a spill according to federal regulations (OSHA 29 CFR 1910.120(a)(v)) and per your company's emergency response plan and guidelines/procedures. See Section 13, Special Instructions for disposal assistance. Outside of the US, only persons properly qualified according to state or local regulations should respond to a spill involving chemicals.

**Containment Technique:** Stop spilled material from being released to the environment.

**Clean-up Technique:** Scoop up spilled material into a large beaker and dissolve with water. Flush the spilled material to the drain with a large excess of water. Decontaminate the area of the spill with a weak acid solution.

**Evacuation Procedure:** Evacuate as needed to perform spill clean-up. If conditions warrant, increase the size of the evacuation.

**Special Instructions (for accidental release):** Not applicable

**304 EHS RQ (40 CFR 355):** Not applicable

**D.O.T. Emergency Response Guide Number:** None

---

## 7. HANDLING / STORAGE

**Handling:** Avoid contact with eyes skin Do not breathe dust. Wash thoroughly after handling. Maintain general industrial hygiene practices when using this product.

**Storage:** Protect from: moisture Keep container tightly closed when not in use.

**Flammability Class:** Not applicable

---

## 8. EXPOSURE CONTROLS / PROTECTIVE EQUIPMENT

**Engineering Controls:** Maintain general industrial hygiene practices when using this product.

**Personal Protective Equipment:**

**Eye Protection:** safety glasses with top and side shields

**Skin Protection:** disposable latex gloves In the EU, the selected gloves must satisfy the specifications of EU Directive 89/686/EEC and standard EN 374 derived from it. lab coat

**Inhalation Protection:** adequate ventilation

**Precautionary Measures:** Avoid contact with: eyes skin Do not breathe: dust Wash thoroughly after handling. Protect from: moisture

**TLV:** Not applicable

**PEL:** Not applicable

---

## 9. PHYSICAL / CHEMICAL PROPERTIES

**Appearance:** White to tan powder or crystals

**Physical State:** Solid

**Molecular Weight:** Not applicable

**Odor:** None

**pH:** 5% soln = 9.4

**Vapor Pressure:** Not applicable

**Vapor Density (air = 1):** Not applicable

**Boiling Point:** Not applicable  
**Melting Point:** Decomposes at 130°C (266°F)  
**Specific Gravity/ Relative Density (water = 1; air =1):** 2.97  
**Evaporation Rate (water = 1):** Not applicable  
**Volatile Organic Compounds Content:** Not applicable  
**Partition Coefficient (n-octanol / water):** Not applicable  
**Solubility:**  
    **Water:** Soluble  
    **Acid:** Soluble  
    **Other:** Not determined  
**Metal Corrosivity:**  
    **Steel:** 0.007 in/yr  
    **Aluminum:** 0.002 in/yr

---

## 10. STABILITY / REACTIVITY

**Chemical Stability:** Stable when stored under proper conditions.  
**Conditions to Avoid:** Excess moisture Exposure to air. Heating to decomposition.  
**Reactivity / Incompatibility:** Incompatible with: alkaloidal salts metallic salts bromine trifluoride chlorine trifluoride  
**Hazardous Decomposition:** Toxic fumes of: potassium oxide iodine iodine compounds  
**Hazardous Polymerization:** Will not occur.

---

## 11. TOXICOLOGICAL INFORMATION

**Product Toxicological Data:**  
    **LD50:** None reported  
    **LC50:** None reported  
    **Dermal Toxicity Data:** None reported  
    **Skin and Eye Irritation Data:** None reported  
    **Mutation Data:** None reported  
    **Reproductive Effects Data:** None reported  
**Ingredient Toxicological Data:** Potassium Iodide: Oral mouse LD<sub>Lo</sub> = 1862 mg/kg; Potassium Carbonate: Oral rat LD<sub>50</sub> = 1870 mg/kg  
No toxicological data available for the ingredients of this product.

---

## 12. ECOLOGICAL INFORMATION

**Product Ecological Information:** --  
No ecological data available for this product. Mobility in soil: No data available  
**Ingredient Ecological Information:** --  
No ecological data available for the ingredients of this product.

---

## 13. DISPOSAL CONSIDERATIONS

**EPA Waste ID Number:** None  
**Special Instructions (Disposal):** Dilute material with excess water making a weaker than 5% solution. Open cold water tap completely, slowly pour the material to the drain. Allow cold water to run for 5 minutes to completely flush the system.  
**Empty Containers:** Rinse three times with an appropriate solvent. Dispose of empty container as normal trash.  
**NOTICE (Disposal):** These disposal guidelines are based on federal regulations and may be superseded by more stringent state or local requirements. Please consult your local environmental regulators for more information.

---

## 14. TRANSPORT INFORMATION

**D.O.T.:**  
    **D.O.T. Proper Shipping Name:** Not Currently Regulated  
    --  
    **DOT Hazard Class:** NA  
    **DOT Subsidiary Risk:** NA  
    **DOT ID Number:** NA  
    **DOT Packing Group:** NA

**I.C.A.O.:**

**I.C.A.O. Proper Shipping Name:** Not Currently Regulated

--

**ICAO Hazard Class:** NA

**ICAO Subsidiary Risk:** NA

**ICAO ID Number:** NA

**ICAO Packing Group:** NA

**I.M.O.:**

**I.M.O. Proper Shipping Name:** Not Currently Regulated

--

**I.M.O. Hazard Class:** NA

**I.M.O. Subsidiary Risk:** NA

**I.M.O. ID Number:** NA

**I.M.O. Packing Group:** NA

**Additional Information:** There is a possibility that this product could be contained in a reagent set or kit composed of various compatible dangerous goods. If the item is NOT in a set or kit, the classification given above applies. If the item IS part of a set or kit, the classification would change to the following: UN3316 Chemical Kit, Class 9, PG II or III. If the item is not regulated, the Chemical Kit classification does not apply.

---

## 15. REGULATORY INFORMATION

**U.S. Federal Regulations:**

**O.S.H.A.:** This product meets the criteria for a hazardous substance as defined in the Hazard Communication Standard. (29 CFR 1910.1200)

**E.P.A.:**

**S.A.R.A. Title III Section 311/312 Categorization (40 CFR 370):** Immediate (Acute) Health Hazard

**S.A.R.A. Title III Section 313 (40 CFR 372):** This product does NOT contain any chemical subject to the reporting requirements of Section 313 of Title III of SARA.

--

**302 (EHS) TPQ (40 CFR 355):** Not applicable

**304 CERCLA RQ (40 CFR 302.4):** Not applicable

**304 EHS RQ (40 CFR 355):** Not applicable

**Clean Water Act (40 CFR 116.4):** Not applicable

**RCRA:** Contains no RCRA regulated substances.

**C.P.S.C.:** Not applicable

**State Regulations:**

**California Prop. 65:** No Prop. 65 listed chemicals are present in this product.

**Identification of Prop. 65 Ingredient(s):** None

**California Perchlorate Rule CCR Title 22 Chap 33:** Not applicable

**Trade Secret Registry:** Not applicable

**National Inventories:**

**U.S. Inventory Status:** All ingredients in this product are listed on the TSCA 8(b) Inventory (40 CFR 710).

**TSCA CAS Number:** Not applicable

---

## 16. OTHER INFORMATION

**Intended Use:** Laboratory Reagent Sulfite determination

**References:** 29 CFR 1900 - 1910 (Code of Federal Regulations - Labor). Air Contaminants, Federal Register, Vol. 54, No. 12. Thursday, January 19, 1989. pp. 2332-2983. TLV's Threshold Limit Values and Biological Exposure Indices for 1992-1993. American Conference of Governmental Industrial Hygienists, 1992. Technical Judgment. In-house information. The Merck Index, 11th Ed. Rahway, New Jersey: Merck and Co., Inc., 1989.

**Revision Summary:** Substantial revision to comply with EU Reg 1272/2008, Reg 1907/2006 and UN GHS (ST/SG/AC.10/36/Add.3).

---

**Legend:**

NA - Not Applicable

ND - Not Determined

NV - Not Available

w/w - weight/weight

w/v - weight/volume

v/v - volume/volume



**USER RESPONSIBILITY:** Each user should read and understand this information and incorporate it in individual site safety programs in accordance with applicable hazard communication standards and regulations.

**THE INFORMATION CONTAINED HEREIN IS BASED ON DATA CONSIDERED TO BE ACCURATE.  
HOWEVER, NO WARRANTY IS EXPRESSED OR IMPLIED REGARDING THE ACCURACY OF THESE DATA  
OR THE RESULTS TO BE OBTAINED FROM THE USE THEREOF.**

**HACH COMPANY ©2012**