



VWR® Cell Lifter, 19mm Flat Blade with 9mm J-Hook end

Description: Cell Lifter, 19mm Flat Blade with 9mm J-Hook end, Sterilized ,Non-pyrogenic


Purpose: Flat Blade Cell Lifter is manufactured with high-grade polyethylene(PE) which is tensile to ensure minimal cell damage during manual harvesting

Materials

19mm Flat Blade: HDPE/High-density Polyethylene Color: Purple


9mm J-Hook end: HDPE/High-density Polyethylene Color: Purple

Dimensions (Unit: mm)



Features

- Convenient operation, special blade design minimizes cell damage
- Individually wrapped
- Non-pyrogenic
- Sterilation: sterilized by gamma irradiation SAL10⁻⁶(ISO11137)
- Shelf life: 3 years after month of production



VWR® Cell Lifter, 19mm Flat Blade with 9mm J-Hook end

VWR NA ENA. No.	VWR EU Cat. No.	Description	Sterile	Qty/Tray	Qty/Pk	Qty/Case
76036-004	734-2978	19mm Flat Blade with 9mm J-Hook end, Sterilized ,Non-pyrogenic	Yes	1	100	100



VWR® Cell Lifter, 19mm Flat Blade with 2.5mm Narrow Blade end

Description: Cell Lifter, 19mm Flat Blade with 2.5mm Narrow Blade end, Sterilized ,Non-pyrogenic


Purpose: Flat Blade Cell Lifter is manufactured with high-grade polyethylene(PE) which is tensile to ensure minimal cell damage during manual harvesting

Materials

19mm Flat Blade: HDPE/High-density Polyethylene Color: Purple

2.5mm Narrow Blade end: HDPE/High-density Polyethylene Color: Purple

Dimensions (Unit: mm)



- Features

 - Convenient operation,special blade design minimizes cell damage
 - Individually wrapped
 - Non-pyrogenic
 - Sterilation: sterilized by gamma irradiation SAL10⁻⁶(ISO11137)
 - Shelf life: 3 years after month of production



VWR® Cell Lifter, 19mm Flat Blade with 2.5mm Narrow Blade end

VWR NA ENA. No.	VWR EU Cat. No.	Description	Sterile	Qty/Tray	Qty/Pk	Qty/Case
76036-006	734-2979	19mm Flat Blade with 2.5mm Narrow Blade end, Sterilized ,Non-pyrogenic	Yes	1	100	100