



Anachemia
A VWR Company

Product Specification

Material

71426-320

Material description

Platinum-Cobalt colour standard

Characteristics

Specifications

Absorbance (430 nm)

0,110 --> 0,120

Absorbance (455 nm)

0,130 --> 0,145

Absorbance (480 nm)

0,105 --> 0,120

Absorbance (510 nm)

0,055 --> 0,065

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Alicia Bourdeau, Head of the Laboratory - Rouses Point
VWR Chemicals, LLC, 3 Lincoln Blvd, Rouses Point, NY
12979-1045, USA