

VWR®

vwr™
part of avantor

Wire Shelving with Round Posts

01. DURABLE PLATED
SHELVES AND POSTS

02. ADJUSTABLE SHELVES

03. ADJUSTABLE FEET





Starter Unit



Add-On Unit

SHELVING STARTER AND ADD-ON UNITS

- Durable plated shelves and posts
- 800 lb. shelf capacity for shelves up to 48" wide; shelves wider than 48" have 600 lbs. capacity (evenly distributed)
- Shelves are adjustable in 1" increments
- Adjustable feet assure level positioning
- Shipped unassembled for compact storage and freight savings

Starter Units: 54" Posts

Cat. No.	D	W	H
89409-900	18	36	54
89409-902	18	48	54
89409-904	18	60	54
89409-906	18	72	54
89409-908	24	36	54
89409-910	24	48	54
89409-912	24	60	54
89409-914	24	72	54

Starter Units: 63" Posts

Cat. No.	D	W	H
89409-916	18	36	63
89409-918	18	48	63
89409-920	18	60	63
89409-922	18	72	63
89409-924	24	36	63
89409-926	24	48	63
89409-928	24	60	63
89409-930	24	72	63

Starter Units: 74" Posts

Cat. No.	D	W	H
89409-932	18	36	74
89409-934	18	48	74
89409-936	18	60	74
89409-938	18	72	74
89409-940	24	36	74
89409-942	24	48	74
89409-944	24	60	74
89409-946	24	72	74

Starter Units: 86" Posts

Cat. No.	D	W	H
89409-948	18	36	86
89409-950	18	48	86
89409-952	18	60	86
89409-954	18	72	86
89409-956	24	36	86
89409-958	24	48	86
89409-960	24	60	86
89409-962	24	72	86

Add-On Units: 54" Posts

Cat. No.	D	W	H
89409-964	18	36	54
89409-966	18	48	54
89409-968	18	60	54
89409-970	18	72	54
89409-972	24	36	54
89409-974	24	48	54
89409-976	24	60	54
89409-978	24	72	54

Add-On Units: 63" Posts

Cat. No.	D	W	H
89409-980	18	36	63
89409-982	18	48	63
89409-984	18	60	63
89409-986	18	72	63
89409-988	24	36	63
89409-990	24	48	63
89409-992	24	60	63
89409-994	24	72	63

Add-On Units: 74" Posts

Cat. No.	D	W	H
89409-996	18	36	74
89409-998	18	48	74
89410-000	18	60	74
89410-002	18	72	74
89410-004	24	36	74
89410-006	24	48	74
89410-008	24	60	74
89410-010	24	72	74

Add-On Units: 86" Posts

Cat. No.	D	W	H
89410-012	18	36	86
89410-014	18	48	86
89410-016	18	60	86
89410-018	18	72	86
89410-020	24	36	86
89410-022	24	48	86
89410-023	24	60	86
89410-024	24	72	86



WIRE CARTS

- 1200 lb. total capacity
- 5" non-marking swivel casters
- Shelves are adjustable in 1" increments
- Two- and three-shelf utility cart models feature integrated push handles

Tall Carts 60" Height

Cat. No.	D	W	H
89410-026	18	36	60
89410-028	18	48	60
89410-030	18	60	60
89410-032	18	72	60
89410-034	24	36	60
89410-036	24	48	60
89410-038	24	60	60
89410-040	24	72	60

Tall Carts 68" Height

Cat. No.	D	W	H
89410-042	18	36	68
89410-044	18	48	68
89410-046	18	60	68
89410-048	18	72	68
89410-050	24	36	68
89410-052	24	48	68
89410-054	24	60	68
89410-056	24	72	68

Tall Carts 80" Height

Cat. No.	D	W	H
89410-058	18	36	80
89410-060	18	48	80
89410-062	18	60	80
89410-064	18	72	80
89410-066	24	36	80
89410-068	24	48	80
89410-070	24	60	80
89410-072	24	72	80

2-Shelf Utility Carts

Cat. No.	D	W
89410-086	18	36
89410-088	18	48
89410-090	18	60
89410-092	24	36
89410-094	24	48
89410-096	24	60

3-Shelf Utility Carts

Cat. No.	D	W
89410-074	18	36
89410-076	18	48
89410-078	18	60
89410-080	24	36
89410-082	24	48
89410-084	24	60



Ready when you are!

VWR® REDISHIP

If you need essential laboratory furniture at a moment's notice, VWR Rediship is your source.

Our fast-track inventory program ensures delivery of top-quality furniture items in a matter of days, helping you operate without interruption.

In addition to the products listed below, VWR Rediship includes carts, tables, shelving and much more.

VWR Contour™ Casework:

This U.S.-made collection of steel laboratory furniture features smooth, easy-to-maintain surfaces. Heavy-gauge steel and full-frame construction are specially engineered to withstand a full ton of pressure.

VWR Lab Chairs and Stools:

This complete line of ergonomic, comfortable seating features three adjustable-height ranges. You can select upholstered black vinyl for maximum comfort and support, or self-skinned urethane, which resists moisture, stains and damage. These products include a five-year warranty and ship to your location within 48 hours.

Labconco® Hoods and Enclosures:

Protector® Premier Laboratory Hoods, Xpert™ Balance Enclosures, and Purifier Delta Series® Class II Safety Cabinets are all available through VWR Rediship.

To learn more about VWR Rediship, contact your furniture specialist.

Most Rediship items ship within 72 hours of order. VWR Lab Chairs ship within 48 hours. All items are subject to availability at time of order. Large-quantity orders may require longer lead times. 72-hour shipping time is based on immediate credit approval.



VWR.COM

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by country, federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. Offers valid in countries listed above, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2020 Avantor, Inc. All rights reserved.