

Biology

From the highest quality living and preserved specimens to the widest selection of media, microscope slides, and more, we have more options for your biology classroom than anywhere else.



	SAM 3G®	SAM II	SAM Basic™	SAM Basic™ without wi-fi	PAT™ the Pediatric Auscultation Trainer	PAT Basic™	PAT Basic™ without wi-fi
Item No.	<u>470317-952</u>	<u>470317-920</u>	<u>470317-928</u>	<u>470317-924</u>	<u>470317-922</u>	<u>470317-926</u>	<u>470317-930</u>
Hardware							
Laptop with pre-installed software	✓	✓	-	-	✓	-	-
Stethoscope used	any stethoscope	any stethoscope	SimScope™ wi-fi	SimScope™	any stethoscope	SimScope™ wi-fi	SimScope™
Carotid pulse	✓	✓	-	-	-	-	-
Can be connected to the speaker system 470317-948 for group instruction	✓	✓	-	-	✓	-	-
Auscultation sites							
Heart	4	4	4	4	4	4	4
Lung	8	8	8	8	8	8	8
Bowel	2	2	2	2	2	2	2
Bruit	1	1	1	1	1	1	1
Auscultation sounds							
Heart	49	35	18	18	76	25	25
Lung	34	24	19	19	17	12	12
Bowel	11	11	5	5	4	4	4
Bruit	6	6	0	0	0	0	0
Software							
Phonocardiogram (recording of all the sounds made by the heart during a cardiac cycle)	✓	✓	-	-	✓	-	-
Sounds in the Real Sound Library	24	0	0	0	0	0	0
Echo Cardiogram Videos	12	0	0	0	0	0	0
Case Videos	24	0	0	0	0	0	0
Users can create and save own videos	✓	-	-	-	-	-	-
Users can setup password protected lectures	✓	✓	-	-	✓	-	-