

PRODUCT DATA SHEET

PART NUMBER:414004-268

PART DESCRIPTION:Aspirating Pipet-25ml

MATERIAL

Polystyrene

PRODUCT REQUIREMENTS

Non-pyrogenic,Gamma Irradiated

Dimensions(mm)

Tip:

Length=40
Interior Dia.=2.4
Outside Dia.=4.75

PIPET SIZE:25ml

INCREMENTS:N/A

GRADUATIONS:N/A

Mouthpiece:

Length=25
Interior Dia.=5.7
Outside Dia.=8

OVERFILL CAPACITY:N/A


COLOR CODE:N/A

Extrusion

Interior Dia.=13.15
Outside Dia.=15.1
Length=245

GENERAL DIMENSIONS

A=15.1
B=8
C=4.75
D=310



PACKAGING REQUIREMENTS

Individually wrapped packaging,50/bag,200/case

INSPECTION REQUIREMENTS

Representative production samples are collected and inspected in accordance with current applicable product specifications.