

## Quality Certificate

### VWR® SPINBAR® MAGNETIC STIR BARS, POLYGON WITHOUT PIVOT RING

As per the manufacturer, the below product meets the following criteria:

<b>North American Catalog No:</b>	58948-981
<b>Lot Number:</b>	J255389
<b>Description:</b>	VWR Bar PTFE 3x1/2 Polygon.
<b>Date of Manufacture:</b>	January 23, 2020
<b>Product Shelf Life:</b>	Not Applicable
<b>Country of Origin:</b>	United States of America

### Quality Conformance

The product is in conformance with the Purchase Order Contract, drawings, specifications, and other applicable documentation.

### Product Specifications

**Material:** Dupont White Teflon #807N-X; Magnet, ALNICO 5

**Measurement:** 3" x 1/2" (76mm x 12.7mm), +/- .020"

**Other Compliances:**

- Autoclavable at 121° C
- Materials Meet FDA Regulations 21CFR177.1550, 21CFR177.1520.
- Teflon used to totally encapsulate the magnetic bar meets US Pharmacopeias Class VI.
- Based on the information provided the product does not contain or come in contact with lead, latex or other contaminants.
- Based on information provided the product does not contain or come in contact with any animal origin.
- All items are carefully inspected to conform with our ISO 9001:2008

**Signed:**



Jamie Ethier  
VP Global Quality  
VWR, Part of Avantor

**Date:** October 15, 2020