



## VWR® Anti-Static Weigh Boats

These flexible weighing boats with smooth surfaces enable accurate pour-out with minimal sample loss. They are biologically inert and resistant to dilute and weak acids, aqueous solutions, alcohols, and bleaches. Ideal for use with microwaves, balances, or moisture analysers.

Available in white, black, or blue for best color contrast to ensure sample is completely removed.

Biologically inert and resistant to dilute and weak acids, aqueous solutions, alcohols, and bases.

Temperature range: -10 to 70°C (14 to 158°F)

- Can be used with microwave moisture analyzers
- Helps reduce the risk of building up static on samples
- FDA/food-grade compliant
- Disposable

### Technical Data

Test Item	Limit	Unit	MDL	001
FDA21.CFR 177 Total Residual Styrene	0.5	% (w/w)	0.1	ND
Result	Pass			

## VWR® Square Weigh Boats

- Shallow form with wide, flat bottom avoids tipping standard
- Rounded corners simplify transfer

Size	Color	Dimensions, mm (in)	Weight, kg (lbs)	Cat. No.	Pack of
Small	White	46 x 46 x 8 (6.2 x 4 x 2.4)	0.34 (0.75)	10770-442	500
Medium	White	85 x 85 x 24 (9.6 x 7.4 x 4.1)	1.43 (3.15)	10770-448	500
Medium	Blue	85 x 85 x 24 (9.5 x 7.5 x 4)	1.15 (2.54)	10770-406	500
Large	White	140 x 140 x 22 (11.8 x 9.2 x 6.4)	3.06 (6.75)	10770-450	500

## VWR® Hexagonal Weigh Boats

- Flexible boats bend at points to form pour spouts
- Flat bottoms and sloping sides provide accurate transfer

Size	Color	Dimensions, mm (in)	Weight, kg (lbs)	Cat. No.	Pack of
Small	White	38 x 25 x 9 (7.5 x 4.6 x 2.4)	0.32 (0.7)	10770-452	500
Medium	White	64 x 51 x 19 (8.4 x 7.6 x 3.8)	0.84 (1.85)	10770-454	500
Large	White	102 x 76 x 22 (11.8 x 8.1 x 5.7)	2.06 (4.55)	10770-456	500
X-Large	White	127 x 89 x 31 (12.75 x 8.4 x 6.5)	2.97 (6.55)	10770-458	500

## VWR® Pour Boat Weighing Dishes

- Pour-spout design easily dispenses samples without flexing the dish
- Flat-bottom dishes are made from 12-mil thick polystyrene for better strength, stability, and handling

Size	Color	Dimensions, mm (in)	Weight, kg (lbs)	Cat. No.	Pack of
Small	White	43 x 58 x 13 (9.5 x 4 x 2.85)	0.34 (0.76)	10770-408	500
Medium	White	83 x 132 x 26 (9.75 x 8 x 5.75)	1.47 (3.24)	10770-438	500
Large	White	108 x 183 x 26 (10.2 x 9.1 x 8.25)	2.54 (5.6)	10770-440	500

