

VWR® PLASTICWARE

Wire Drying Rack

Laboratory Notebooks

Weighboats



TUBE RACKS


Rota-Rack® Duo

Dilution and concentration studies made easy

- made from polypropylene
- transfer solutions between 0.2 mL PCR® tubes and 0.5 and 1.5 mL micro centrifuge tubes
- rotate and customize to suit your application needs. 4 combinations possible
- blue section holds 48 x 0.5 mL and 48 x 1.5 mL microfuge tubes on the reverse side
- purple section holds 96 x 0.2 mL tubes in a standard array pattern one and 48 x 1.5 mL tubes on the reverse side
- tapered wells securely suspend tubes and makes removal easy
- alphanumeric well identification
- **autoclavable**



| cat. no. | tubes | wells | well dia. | l x w x h | pkg. | color |
|-----------|------------|-------|-----------|---|------|-------------|
| 89140-272 | 0.2 mL | 96 | 5.0 mm | 11.6 x 4.6 x 1.33 in (29.6 x 11.8 x 3.4 cm) | 1/ea | blue/purple |
| | 0.5 mL | 48 | 8.0 mm | | | |
| | 1.5/2.0 mL | 96 | 11 mm | | | |



PCR® Racks Reversible

Removable hinged lid fits both sides

- PCR® side of the rack holds 8 or 12 x 0.2 mL tube strips or individual tubes
- reverse side holds 0.5 mL tubes
- both sides of the rack have 12 wells for 1.5 mL tubes
- **autoclavable**

| cat. no. | tubes | wells | well dia. | l x w x h | pkg. | color |
|-----------|------------|-------|-----------|--|------|-------|
| 82024-454 | 0.2 mL | 168 | 5.5 mm | 8.11 x 5.16 x 2.13 in (20.6 x 13.1 x 5.4 cm) | 5/pk | ■ ■ |
| | 0.5 mL | 40 | 7.6 mm | 8.11 x 5.16 x 2.13 in (20.6 x 13.1 x 5.4 cm) | 5/pk | ■ ■ |
| | 1.5/2.0 mL | 12 | 11.4 mm | 8.11 x 5.16 x 2.13 in (20.6 x 13.1 x 5.4 cm) | 5/pk | ■* |

*assorted colors (blue, green, pink, yellow, orange)




Storage Racks for 80-Well Racks

Strong, clear acrylic construction

these racks help you keep large numbers of microtubes organized on your benchtop.

- holds up to eight 80-well microtube racks filled with either 1.5 or 2.0 mL microtubes - see page 20
- plenty of room to slide the racks in and out

| cat. no. | description | l x w x d |
|-----------|-----------------|--------------------|
| 82024-482 | vertical rack | 157 x 236 x 234 mm |
| 82024-480 | horizontal rack | 310 x 236 x 122 mm |



Cube Rack

Interlocking multi-tube racks

- Made of polypropylene
- Allows you to work with different sized tubes all in the same plane
- **Autoclavable**



| cat. no. | tubes | wells | well dia. | l x w x h | pkg. | color |
|------------------|------------|-------|-----------|--|-------|-----------|
| 82010-746 | 50 mL | 4 | 30 mm | 4.25 x 4.45 x 4.45 in (10.8 x 11.3 x 11.3 cm) | 5/pkg | ■ ■ ■ ■ * |
| | 15 mL | 10 | 17 mm | | | |
| | 0.5 mL | 12 | 11.5 mm | | | |
| | 1.5/2.0 mL | 16 | 12.7 mm | | | |

*Assorted colors (Blue, Green, Pink, Yellow, Orange)



Cryogenic Vial Holders

Safely store and transport cryogenic vials

- made from polycarbonate—can withstand cryogenic freezing down to -196°C
- well tabs allow one-handed opening and closing of vials that have a slot feature
- embossed index numbers and squared, upper left corner provide easy orientation
- ribbed grips make handling with gloves easier
- **autoclavable**

| cat. no. | tubes | wells | wells dia. | l x w x h | pkg. | color |
|------------------|--------|-------|------------|--|-------|--------|
| 82024-334 | 2-5 mL | 50 | 12.9 mm | 4.1 x .79 x 0.9 in (10.5 x 2 x 2.3 cm) | 4/pkg | blue ■ |



Reversible Racks 96-Well

Easy-to-grasp tabs for sample handling and storage

reverse 96 wells for 0.5 mL tubes on one side and 96 wells for 1.5/2.0 mL tubes on the reverse side

- made of polypropylene
- wells are arranged in standard 8 x 12 array
- tabs on each end of rack snap lid into place
- 0.5 mL wells separated by molded-in grid lines, making them easier to see
- raised rings around the 1.5 to 2.0 mL wells make tubes easier to remove
- **autoclavable**

| cat. no. | tubes | wells | well dia. | l x w x h | pkg. | color |
|------------------|-------------------|-------|-----------|---------------------------------------|-------|-----------|
| 82024-494 | 1.5/2.0 mL 0.5 mL | 96 | 11 / 7.6 | 9.6 x 4.7 x 2 in (24.6 x 12.1 x 5 cm) | 5/pkg | ■ ■ ■ ■ * |
| 82024-484 | 1.5/2.0 mL 0.5 mL | 96 | 11 / 7.6 | 9.6 x 4.7 x 2 in (24.6 x 12.1 x 5 cm) | 5/pkg | natural □ |
| 82024-488 | 1.5/2.0 mL 0.5 mL | 96 | 11 / 7.6 | 9.6 x 4.7 x 2 in (24.6 x 12.1 x 5 cm) | 5/pkg | blue ■ |
| 82024-490 | 1.5/2.0 mL 0.5 mL | 96 | 11 / 7.6 | 9.6 x 4.7 x 2 in (24.6 x 12.1 x 5 cm) | 5/pkg | green ■ |
| 82024-492 | 1.5/2.0 mL 0.5 mL | 96 | 11 / 7.6 | 9.6 x 4.7 x 2 in (24.6 x 12.1 x 5 cm) | 5/pkg | pink ■ |
| 82024-486 | 1.5/2.0 mL 0.5 mL | 96 | 11 / 7.6 | 9.6 x 4.7 x 2 in (24.6 x 12.1 x 5 cm) | 5/pkg | yellow ■ |

*assorted colors (blue, green, pink, yellow, orange)

TUBE RACKS

4-Way Rack

Link together to cater for your application needs

create as large a rack as you require for your bench top applications

- made from polypropylene
- **autoclavable**



| cat. no. | tubes | wells | well dia. | l x w x h | pkg. | color | |
|------------------|------------|-------|-----------|--|------|-----------------|-------------|
| 82010-752 | 50 mL | 4 | 30 mm | 6.85 x 3.74 x 2.04 in (17.4 x 9.5 x 5.2 cm) | 5/pk | Natural | □ |
| | 15 mL | 12 | 17.5 mm | | | | |
| | 0.5 mL | 32 | 8.2 mm | | | | |
| | 1.5/2.0 mL | 32 | 11.6 mm | | | | |
| 82024-460 | 50 mL | 4 | 30 mm | 6.85 x 3.74 x 2.04 in (17.4 x 9.5 x 5.2 cm) | 5/pk | Green | ■ |
| | 15 mL | 12 | 17.5 mm | | | | |
| | 0.5 mL | 32 | 8.2 mm | | | | |
| | 1.5/2.0 mL | 32 | 11.6 mm | | | | |
| 82024-462 | 50 mL | 4 | 30 mm | 6.85 x 3.74 x 2.04 in (17.4 x 9.5 x 5.2 cm) | 5/pk | Pink | ■ |
| | 15 mL | 12 | 17.5 mm | | | | |
| | 0.5 mL | 32 | 8.2 mm | | | | |
| | 1.5/2.0 mL | 32 | 11.6 mm | | | | |
| 82024-464 | 50 mL | 4 | 30 mm | 6.85 x 3.74 x 2.04 in (17.4 x 9.5 x 5.2 cm) | 5/pk | Yellow | ■ |
| | 15 mL | 12 | 17.5 mm | | | | |
| | 0.5 mL | 32 | 8.2 mm | | | | |
| | 1.5/2.0 mL | 32 | 11.6 mm | | | | |
| 82024-458 | 50 mL | 4 | 30 mm | 6.85 x 3.74 x 2.04 in (17.4 x 9.5 x 5.2 cm) | 5/pk | Orange | ■ |
| | 15 mL | 12 | 17.5 mm | | | | |
| | 0.5 mL | 32 | 8.2 mm | | | | |
| | 1.5/2.0 mL | 32 | 11.6 mm | | | | |
| 82010-750 | 50 mL | 4 | 30 mm | 6.85 x 3.74 x 2.04 in (17.4 x 9.5 x 5.2 cm) | 5/pk | Assorted colors | ■ ■ ■ ■ ■ * |
| | 15 mL | 12 | 17.5 mm | | | | |
| | 0.5 mL | 32 | 8.2 mm | | | | |
| | 1.5/2.0 mL | 32 | 11.6 mm | | | | |

*Assorted colors (Blue, Green, Pink, Yellow, Orange)

Microtube Racks 80-Well

Popular fraction-collector-plate style rack

- made of polypropylene
- easily fill and transfer 80 tubes in compact arrangement compact arrangement (5 x 16 array)
- molded-in coordinates to aid organization
- **autoclavable**



| cat. no. | tubes | wells | well dia. | l x w x h | pkg. | color | |
|------------------|------------|-------|-----------|--|------|-----------------|-------------|
| 82024-467 | 1.5/2.0 mL | 80 | 11 mm | 8.9 x 2.6 x 1 in (22.5 x 6.7 x 2.8 cm) | 5/pk | natural | □ |
| 82024-473 | 1.5/2.0 mL | 80 | 11 mm | 8.9 x 2.6 x 1 in (22.5 x 6.7 x 2.8 cm) | 5/pk | blue | ■ |
| 82024-475 | 1.5/2.0 mL | 80 | 11 mm | 8.9 x 2.6 x 1 in (22.5 x 6.7 x 2.8 cm) | 5/pk | green | ■ |
| 82024-477 | 1.5/2.0 mL | 80 | 11 mm | 8.9 x 2.6 x 1 in (22.5 x 6.7 x 2.8 cm) | 5/pk | pink | ■ |
| 82024-471 | 1.5/2.0 mL | 80 | 11 mm | 8.9 x 2.6 x 1 in (22.5 x 6.7 x 2.8 cm) | 5/pk | yellow | ■ |
| 82024-469 | 1.5/2.0 mL | 80 | 11 mm | 8.9 x 2.6 x 1 in (22.5 x 6.7 x 2.8 cm) | 5/pk | orange | ■ |
| 82024-479 | 1.5/2.0 mL | 80 | 11 mm | 8.9 x 2.6 x 1 in (22.5 x 6.7 x 2.8 cm) | 5/pk | Assorted colors | ■ ■ ■ ■ ■ * |

*assorted colors (blue, green, pink, yellow, orange)



Foam Tube Racks

High storage, minimal footprint

- two-piece racks are constructed of expanded polystyrene (EPS) foam
- lightweight microtube racks are ideal for long-term studies or tube storage
- extra spacing between wells allows for easy gripping with fingertips
- racks stack securely with nesting features on lid and base

| cat. no. | tubes | wells | wells dia. | l x w x h | pkg. | color |
|------------------|------------|-------|------------|--|-------|--------------------------------|
| 82017-624 | 1.5/2.0 mL | 100 | 11.5 mm | 13.22 x 3.74 x 2.87 in (33.6 x 9.5 x 7.3 cm) | 10/pk | white <input type="checkbox"/> |
| 89194-766 | 1.5/2.0 mL | 50 | 11.5 mm | 8.26 x 4.33 x 2.7 in (21.0 x 11.0 x 7.1 cm) | 10/pk | white <input type="checkbox"/> |

Floating Microtube Rack

- made from polypropylene
- use on the bench, with water baths or ice containers
- keeps microtubes afloat in a water bath or can be pushed into crushed ice — perfect for controlling sample temperature.
- detachable legs provide two options:
 - attach legs so rack stands upright on bench top
 - detach legs so samples can be quickly popped out by pushing down onto any surface
- *not recommended for autoclaving*

| cat. no. | tubes | wells | wells dia. | dimension | color |
|------------------|------------|-------|------------|----------------------------|--------------------------------|
| 82024-504 | 1.5/2.0 mL | 16 | 10.9 mm | square 4 in (10.2 cm) | white <input type="checkbox"/> |
| 82024-506 | 1.5/2.0 mL | 16 | 10.9 mm | square 4 in (10.2 cm) | white <input type="checkbox"/> |
| 82024-508 | 1.5/2.0 mL | 20 | 10.9 mm | round 3.8 in (9.8 cm dia.) | white <input type="checkbox"/> |
| 82024-510 | 1.5/2.0 mL | 20 | 10.9 mm | round 3.8 in (9.8 cm dia.) | white <input type="checkbox"/> |



detachable legs

82017-632

82017-628

detachable handle

82017-630

82017-634

Floating Foam Tube Racks

- made from high density polyethylene
- reusable and flexible for a variety of applications
- brightly colored, high density polyethylene foam racks are ideal for floating tubes in water baths, tubs, or beakers
- closed-cell foam resists water absorption so racks can be washed and reused
- racks include detachable carrying handles—easily lift samples out of baths
- the blue round rack features unique “x” slits to hold a variety of tube sizes
- *not recommended for autoclaving*

| cat. no. | tubes | wells | wells dia. | l x w x h | color |
|------------------|------------------------|-------|------------|---|--|
| 82017-628 | 50 mL | 4 | 29 mm | 6.4 x 5.3 x 2.5 in (16.5 x 13.7 x 6.5 cm) | green <input checked="" type="checkbox"/> |
| 82017-630 | 15 mL | 8 | 16 mm | 6.7 x 2.2 x 2.5 in (17.2 x 5.7 x 6.5 cm) | blue <input checked="" type="checkbox"/> |
| 82017-632 | 1.5-2.0 mL | 24 | 9.3 mm | 6 x 3.2 x 2.7 in (15 x 8.3 x 7 cm) | yellow <input checked="" type="checkbox"/> |
| 82017-634 | 0.2/0.5/ 1.5-2.0 mL | 18 | 9.3 | 2.7 in (7 cm dia.) | blue <input checked="" type="checkbox"/> |



TUBE RACKS

Peg Racks

Strong, tough, versatile racks for bench or tube storage

fiberglass reinforced polypropylene racks are ideal for bench work or tube storage

- made from polypropylene
- spacious peg design allows easy viewing of tube contents
- use as a drying rack for glass plates or upside-down tubes
- pegs remain straight even after repeated autoclaving
- embossed indexing, ribbed grips, and drainage holes
- **autoclavable**



| cat. no. | tubes | pegs | wells dia. | l x w x h | pkg. | color | |
|------------------|-------|------|------------|--|------|-------|---|
| 82024-498 | 10 mL | 50 | 14–17 mm | 7.4 x 4.1 x 2.71 in (18.8 x 10.4 x 6.9 cm) | 2/pk | white | □ |
| 82024-502 | 10 mL | 50 | 14–17 mm | 7.4 x 4.1 x 2.71 in (18.8 x 10.4 x 6.9 cm) | 2/pk | blue | ■ |
| 82024-496 | 7 mL | 96 | 10–13 mm | 3.2 x 4.96 x 2.71 in (8.2 x 12.6 x 6.9 cm) | 2/pk | white | □ |
| 82024-500 | 7 mL | 96 | 10–13 mm | 3.2 x 4.96 x 2.71 in (8.2 x 12.6 x 6.9 cm) | 2/pk | blue | ■ |

Transformable Cube Rack

For storage of large samples and 15 and 50 ml tubes

- polypropylene construction complete with two snap-in grids
- imprinted grid and first tube orientation and keyed lid
- use without grids for larger sized product storage
- fill with ice and use as a temporary low temperature workstation
- clear lid, translucent sides
- **autoclavable**

| cat. no. | tubes | wells | wells dia. | l x w x h | pkg. | color | |
|------------------|-------|-------|------------|--|------|--------|---|
| 89140-266 | 15 mL | 16 | 18 mm | 4.9 x 4.9 x 5 in (12.5 x 12.5 x 12.9 cm) | 5/pk | purple | ■ |
| | 50 mL | 9 | 31 mm | | | | |



upright cube freezer rack




chest cube freezer rack

fill with ice and use
as a temporary low
temperature workstation



holds either:

- 15 mL or
- 50 mL samples




89031-458

89031-460



Rota-Rack® Tube Racks

Get organized with this fun, interlocked rack

- made from polypropylene
- rotate each of the four modules independently to configure for the tubes that you want to use.
- suitable for a massive variety of culture, centrifuge, micro centrifuge, round, flat or conical bottomed tubes
- two versions to accommodate small and large tubes
- **autoclavable**



twist and turn to suit your application

| cat. no. | tubes | wells | wells dia. | l x w x h | | pkg. | color | |
|-----------|------------|-------|------------|---|------|-----------------------------|---|--|
| 89031-458 | 50 mL | 2 | 29 mm | 8.81 x 4.33 x 4.33 in (22.4 x 11 x 11 cm) | 1/pk | purple, blue, green, yellow |  | |
| | 10-15 mL | 8 | 20 mm | | | | | |
| | 15-20 mL | 6 | 16.0 mm | | | | | |
| | 10-13 mL | 10 | 13.5 mm | | | | | |
| 89031-460 | 15 mL | 6 | 16.4 mm | 8.2 x 3.9 x 3.9 in (21 x 10 x 10 cm) | 1/pk | purple, blue, green, yellow |  | |
| | 1.5/2.0 mL | 9 | 11.2 mm | | | | | |
| | .5 mL | 12 | 8.0 mm | | | | | |
| | .2 mL | 32 | 6.5 mm | | | | | |

twist and turn to suit your application




adjustable slot adapts to different diameter tubes

Universal Tube Rack

Ideal for holding a wide range of tubes

- made from polypropylene
- holds up to 12 similar size tubes
- imprinted references for easy identification
- side handles for easy transport
- **not recommended for autoclaving**

| cat. no. | tubes | wells | well dia. | l x w x h | pkg. | color |
|-----------|----------|-------|-----------|---------------------------------------|------|---|
| 89234-744 | 15-50 mL | 12 | 16-30 mm | 7.6 x 5.5 x 2 in (19.4 x 14 x 5.1 cm) | 1/ea | blue/purple  |



Pop-Up™ Racks

Collapse for compact storage and “pop-up” for quick use

- made from polypropylene
- easy to open — easy to close
- designed to help organize 15 and 50 mL tubes
- imprinted grid aids quick retrieval
- stack racks in open or closed position
- saves space when racks are not in use
- *not recommended for cleaning in autoclave*



easy-to-open



easy-to-close

| cat. no. | tubes | wells | wells dia. | l x w x h (expanded) | l x w x h (collapsed) | pkg. | color |
|-----------|------------|---------|------------|--|--|------|-------|
| 89140-268 | 15 & 50 mL | 21 / 12 | 17 / 30 mm | 10 x 5.39 x 2.83 in (25.5 x 13.7 x 7.2 cm) | 10 x 5.39 x 0.83 in (25.5 x 13.7 x 2.1 cm) | 2/pk | blue |
| 89140-270 | 15 & 50 mL | 21 / 12 | 17 / 30 mm | 10 x 5.39 x 2.83 in (25.5 x 13.7 x 7.2 cm) | 10 x 5.39 x 0.83 in (25.5 x 13.7 x 2.1 cm) | 2/pk | pink |



racks link together

Puzzle Rack

Fun and functional racks slide to connect

- made from polypropylene
- two-sided puzzle-shaped pieces hold an assortment of 1.5/2.0, 15 and 50 mL test tubes
- pieces can be used separately for small sample applications or connect in a square or row configurations
- imprinted references aid with easy retrieval of tubes
- **autoclavable**

| cat. no. | tubes | wells | wells dia. | l x w x h (expanded) | pkg. | color |
|-----------|-------------|-------|------------|---|------|-----------|
| 89234-740 | 15 mL | 9 | 16.9 mm | 5.9 x 5.31 x 2.59 in (15.1 x 13.5 x 6.6 cm) | 4/pk | ■ ■ ■ ■ * |
| | 50 mL | 4 | 29.5 mm | | | |
| | 1.5/ 2.0 mL | 8/4 | 11.7 mm | | | |

*assorted colors (yellow, blue, green, orange)

Wire Drying Rack

Saves glassware from chips and cracks

drain and dry test tubes, beakers, flasks, and other labware

- place self-standing rack on a table or near a sink, or use hooks to hang rack on the wall
- holds up to 90 pieces of labware on 18 x 127 mm pins and 72 x 76 mm pins
- constructed of spot-welded, vinyl-coated steel wire for durability
- easily accessible and removable bottom drainage pan

| cat. no. | l x w x h |
|------------------|--------------------|
| 82024-450 | 462 x 182 x 525 mm |



Clamp-ette® Pipette Holders

Create an instant pipette holding point within your lab

secures to tabletops, benches and shelves with a twist

- chemical-resistant polypropylene and fits most major pipette brands
- available in a single or triple format
- **autoclavable**

| cat. no. | description | l x w x h | color |
|---------------------------------------|----------------|-------------------|--------|
| 82007-226 | tri clamp-ette | 157 x 120 x 95 mm | blue |
| 82007-228 | tri clamp-ette | 157 x 120 x 95 mm | red |
| 82007-230 | tri clamp-ette | 157 x 120 x 95 mm | purple |
| * assorted colors (blue, red, purple) | | | |



PIPETTES

Pipette Racks

Protects and organizes glass and plastic pipettes

angled four shelved compartments keep your 0.1 to 50 mL pipettes or cans within easy reach while preventing damage and clutter.

- **heavy duty acrylic** version provides crystal-clear viewing from any angle
- **abs versions** have cut-out viewing window and come with two removable attachments that hold pipette bulbs, fillers, and pumps
- **HS20612M magnetized ABS** version has strong magnets on both sides and mounts securely to most metal surfaces

| cat. no. | material | l x w x h | pkg. | color |
|-----------|----------|-------------------|------|--|
| 82024-536 | acrylic | 291 x 88 x 401 mm | 1/ea | clear <input type="checkbox"/> |
| 82007-220 | ABS | 291 x 88 x 401 mm | 1/ea | blue <input checked="" type="checkbox"/> |




Pipette Tip Disposal Box

Quick flip of port cover makes tip ejection easy

convenient for one-handed disposal of used pipette tips. eliminates any need to handle pipette tips used with hazardous fluids

- the 3/8" (10 mm) thick acrylic provides protection from beta radiation
- the large, hinged lid allows easy disposal of contents
- suggest using in conjunction with small biohazard disposal bags (page 81), if hazardous fluids are not acrylic friendly

| cat. no. | l x w x h |
|-----------|--------------------|
| 82010-754 | 154 x 171 x 206 mm |



Cargo Pipette Rack

Customize your storage of pipettes

Organize the cargo pipette rack to meet your laboratory or personal needs. click together horizontally or vertically to build the desired configuration, or use as individual holders, utilizing the three lid sections provided.

- Magnets on both sides of units connect securely and safely to most metal surfaces
- Angled shelves store 0.1 to 50 mL pipettes
- Protective lid encloses the unit to prevent dust build up and provides additional shelf space

| Cat. no. | Lx W x H | Pkg. | Color |
|-----------|----------------------|------|--|
| 89140-276 | 301 x 100 x 127.5 mm | 3/pk | Blue <input checked="" type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> |

Units lock securely together horizontally or vertically

Fillers and Bulbs

| latex dropper bulbs (natural rubber) | cat. no. | capacity | length | OD | hole dia | pkg. | color | |
|---|------------------|----------|--------|-------|----------|-------|---------|---|
| | 82024-550 | 1 mL | 43 mm | 13 mm | 5 mm | 72/pk | natural | ■ |
| | 82024-554 | 2 mL | 52 mm | 13 mm | 6 mm | 72/pk | natural | ■ |

- designed for use with Pasteur pipettes
- small pipettes, and medicine droppers

| rubber bulbs (natural rubber) | cat. no. | capacity | length | OD | hole dia | pkg. | color | |
|----------------------------------|------------------|----------|--------|-------|----------|------|-------|---|
| | 82024-564 | 12 mL | 51 mm | 34 mm | 4 mm | 1/ea | red | ■ |
| | 82024-566 | 15 mL | 56 mm | 34 mm | 4 mm | 1/ea | red | ■ |
| | 82024-568 | 30 mL | 66 mm | 46 mm | 6 mm | 1/ea | red | ■ |
| | 82024-570 | 60 mL | 80 mm | 54 mm | 6 mm | 1/ea | red | ■ |
| | 82024-572 | 150 mL | 103 mm | 73 mm | 6 mm | 1/ea | red | ■ |


- contoured shape and heavy-wall neck increase grip for easier use with pipettes and syringes

| vacuum/pressure bulb (synthetic rubber) | cat. no. | capacity | length | OD | hole dia | pkg. | color | |
|--|------------------|----------|--------|-------|----------|------|-------|---|
| | 82024-578 | 85 mL | 137 mm | 53 mm | 7.4 mm | 1/ea | black | ■ |

- ideal for gas sampling, or for pumping gas through apparatus
- two rubber valves with PVC connections

| safety bulb pipette filler (natural rubber) | cat. no. | capacity | length | OD | hole dia | pkg. | color | |
|--|------------------|----------|--------|-------|------------|------|-------|---|
| | 82024-574 | 15 mL | 106 mm | 41 mm | 4.5-6.5 mm | 3/ea | red | ■ |

- features a white ABS attachment containing a silicone conical adapter that seals around pipettes without twisting
- elastic band for hanging on reagent bottle or wrist
- filler fits standard size pipettes



Pipette Stands

Convenient and safe storage

3, 4, and 6-place pipette stands are designed to hold most major pipette brands.

- acrylic stands are transparent/clear; ABS stands are solid opaque colors

| cat. no. | material | holds | slot upper | slot lower | color | |
|------------------|----------|-------|------------|------------|-------|---|
| 82024-540 | acrylic | three | 25 mm | 20 mm | clear | □ |
| 82024-542 | acrylic | four | 25 mm | 20 mm | clear | □ |
| 82010-756 | acrylic | four | 25 mm | 20 mm | clear | □ |
| 82024-544 | acrylic | six | 25 mm | 20 mm | clear | □ |

Premium Microscope Slide Boxes

Slide boxes are perfect for storage or transportation of standard size slides (25 x 75 mm and 1 x 3 inch)

- molded, rigid boxes securely hold slides in grooved slots
- durable ABS plastic boxes are available with either cork or foam lining
- both the lining and the slots are numbered to correspond to the slide inventory sheet



25 place slide box

| easy access slides with friction-fit cover | cat. no. | capacity | lining | hinge | l x w x h | pkg. | color |
|--|-----------|-----------|--------|--------------|------------------|-------|---------|
| | 82003-418 | 25 slides | cork | friction fit | 141 x 88 x 35 mm | 72/cs | blue ■ |
| | 82003-420 | 25 slides | cork | friction fit | 141 x 88 x 35 mm | 72/cs | green ■ |
| | 82003-422 | 25 slides | cork | friction fit | 141 x 88 x 35 mm | 72/cs | red ■ |
| | 82003-426 | 25 slides | cork | friction fit | 141 x 88 x 35 mm | 72/cs | white □ |



50 place slide box

| thumb-latch lock mechanism | cat. no. | capacity | lining | hinge | l x w x h | pkg. | color |
|----------------------------|-----------|-----------|--------|-------------|------------------|-------|---------|
| | 82024-580 | 50 slides | cork | thumb-latch | 209 x 86 x 35 mm | 72/cs | blue ■ |
| | 82024-582 | 50 slides | cork | thumb-latch | 209 x 86 x 35 mm | 72/cs | red ■ |
| | 82024-584 | 50 slides | cork | thumb-latch | 209 x 86 x 35 mm | 72/cs | white □ |
| | 82024-588 | 50 slides | cork | thumb-latch | 209 x 86 x 35 mm | 72/cs | green ■ |

available in cork lining only



100 place slide box

| closes securely with a rust-resistant, nickel-plated clasp and pin hinges | cat. no. | capacity | lining | hinge | l x w x h | pkg. | color |
|---|-----------|------------|--------|--------------|-------------------|-------|----------|
| | 82003-406 | 100 slides | cork | nickle clasp | 208 x 175 x 34 mm | 72/cs | blue ■ |
| | 82003-408 | 100 slides | cork | nickle clasp | 208 x 175 x 34 mm | 72/cs | green ■ |
| | 82003-410 | 100 slides | cork | nickle clasp | 208 x 175 x 34 mm | 72/cs | red ■ |
| | 82003-412 | 100 slides | cork | nickle clasp | 208 x 175 x 34 mm | 72/cs | yellow ■ |
| | 82003-414 | 100 slides | cork | nickle clasp | 208 x 175 x 34 mm | 72/cs | white □ |
| | 82003-416 | 100 slides | cork | nickle clasp | 208 x 175 x 34 mm | 72/cs | gray ■ |



100 place tall slide box

| stackable with 100-place premium slide boxes | cat. no. | capacity | lining | hinge | l x w x h | pkg. | color |
|--|-----------|------------|--------|--------------|-------------------|-------|--------|
| | 89140-280 | 100 slides | cork | nickle clasp | 208 x 175 x 60 mm | 72/cs | blue ■ |

available in cork lining only

additional space for larger 3 x 2 in slides



Cardboard Slide Tray

Cardboard slide trays hold standard slides (25 x 75 mm and 1 x 3 inches)

- horizontal compartments are recessed to protect the specimens when the cover is closed
- features thumb cut for easy slide removal
- available in 2 sizes with 5 color options for the 20-place trays
- **disposable**

| cat. no. | capacity | l x w x h | qty. | trim/imprinting |
|------------------|-----------|-------------------|-------|-----------------|
| 82020-913 | 20 slides | 208 x 341 x 9 mm | 72/cs | black ■ |
| 82020-915 | 30 slides | 281 x 328 x 10 mm | 40/cs | black ■ |



20 and 30-place

Slide Mailers

Disposable, safe, secure transport of microscope slides

cardboard slide mailers

| heavy cardboard mailers for safe transport | cat. no. | capacity | l x w x h | qty. | color |
|--|------------------|----------|-----------------|-------|-----------|
| | 82020-908 | 1 place | 95 x 42 x 5 mm | 25/pk | natural ■ |
| | 82020-911 | 2 place | 103 x 79 x 5 mm | 36/pk | natural ■ |

- thumb grooves easy slide removal
- 1-and-2-place slide mailer versions
- each will hold 25 x 75 mm (1 x 3 in) and or 26 x 76 mm sized slides
- **disposable**



1-place

2-place

plastic 2-place slide mailers

| ideal for secure shipping and storing of slides | cat. no. | capacity | l x w x h | qty. | color |
|---|------------------|----------|----------------|-------|-----------|
| | 82007-222 | 2 place | 93 x 68 x 5 mm | 25/pk | ■ ■ ■ ■ * |
| | 82007-224 | 2 place | 93 x 68 x 5 mm | 25/pk | natural □ |

*assorted colors (red, yellow, blue, purple, orange)

- polypropylene protectors hold two standard (25 x 75 mm and 1 x 3 in) slides securely
- large, easy to-secure snaps, provide safe storage for valuable slides
- variety of colors for coding
- **disposable**



plastic 5-place slide mailers

| ideal for secure shipping, storing, or staining slides | cat. no. | capacity | l x w x h | qty. | color |
|--|------------------|----------|-----------------|-------|-----------|
| | 82024-598 | 5 place | 81 x 19 x 31 mm | 25/pk | natural □ |
| | 82003-431 | 5 place | 43 x 24 x 88 mm | 25/pk | ■ ■ ■ ■ * |
| | 82003-433 | 5 place | 43 x 24 x 88 mm | 25/pk | natural □ |

*assorted colors (red, yellow, blue, purple, orange)

- polypropylene containers hold five standard slides (25 x 75 mm and 1 x 3 in)
- inside wall features slotted channels to keep slides separated and to make removal easier
- lid closes securely to keep slides safe
- **disposable**



82024-598

side opening

end opening



20-Place Slide Tray Racks

Search, identify, and easily locate the specimen slide you need

store multiple slide trays for easy retrieval—ideal for microbiology, histology, and pathology fields

- ABS rack holds five 20-place slide trays (*not included see above*)
- eliminates the worry of breakage during retrieval—easy access to each of the trays
- matching clasps provide area for identification labels and doubles as a handle
- portable and refrigerator safe for cold sample storage
- stack multiple units to create your own library

| cat. no. | capacity | l x w x h | color |
|------------------|---------------|--------------------|--------|
| 89140-274 | 5 slide trays | 226 x 340 x 125 mm | blue ■ |



Microscope Slide Dispenser

Designed for easy load and slide dispensing

- works effortlessly with 25 x 75 mm (1 x 3 in) and 26 x 76 mm slides
- turn either side wheel to dispense slides one at a time


| cat. no. | capacity | l x w x h | color |
|------------------|-----------------|--------------------|--------|
| 82024-522 | up to 91 slides | 113 x 123 x 163 mm | blue ■ |

Slide Holders

Hold slides during evaluation, cleaning, and drying

- made from chemical-resistant polyethylene
- solid, strong construction with rubber feet for bench top stability

| cat. no. | capacity | l x w x h |
|------------------|-----------|-------------------|
| 82024-526 | 40 slides | 298 x 105 x 24 mm |
| 82024-528 | 78 slides | 299 x 178 x 21 mm |
| 82024-524 | 15 slides | 153 x 102 x 21 mm |



Standard Microscope Slide Boxes

Constructed of durable polypropylene

slide boxes hold either 25 or 100 standard slides (25 x 75 mm and 1 x 3 inches)

- hinged lids are held securely closed with a latch lock
- stacking feet maximize stability
- slide identification card

| unlined and stackable | cat. no. | capacity | hinge | l x w x h | color |
|-----------------------|------------------|------------|------------|-------------------|---------|
| | 82024-614 | 25 slides | latch lock | 141 x 92 x 30 mm | blue ■ |
| | 82024-616 | 25 slides | latch lock | 141 x 92 x 30 mm | red ■ |
| | 82024-612 | 25 slides | latch lock | 141 x 92 x 30 mm | white □ |
| | 82024-608 | 100 slides | latch lock | 210 x 169 x 37 mm | blue ■ |
| | 82024-610 | 100 slides | latch lock | 210 x 169 x 37 mm | red ■ |
| | 82024-606 | 100 slides | latch lock | 210 x 169 x 37 mm | white □ |





Storage Rack for Microscope Slide Boxes

Create a file system for 50 and 100 place slide boxes

storage rack makes organization of samples easy

- clear acrylic construction allows the labels to be seen on the slide boxes
- place boxes in rows or columns by color to create an easy coding system

| cat. no. | bins | l x w x h |
|------------------|---------|-----------------|
| 82024-620 | 15 bins | 711 x 171 x 228 |



Freezer Racks

Offering a variety of racks to safely store samples in your freezer

choose from a wide selection of corrosion resistant stainless steel racks that are compatible with major brands of upright and chest freezers

- all have handles for easy use
- racks store boxes up to a 135 x 135 mm footprint and accommodate 2" tall boxes (up to 5 cm) as well as 3" tall boxes (up to 7.6 cm) depending on the rack model

upright freezer racks 2 inch boxes

| cat. no. | config | box size | l x w x h |
|------------------|--------|---------------|---|
| 82024-520 | 4 x 4 | 2 in (5.3 cm) | 21.8 x 5.5 x 8.9 in (55.6 x 14 x 22.5 cm) |

upright freezer racks 3 inch boxes

| cat. no. | config | box size | l x w x h |
|------------------|--------|---------------|---|
| 82024-518 | 3 x 4 | 3 in (7.6 cm) | 22 x 5.5 x 8.9 in (55.6 x 14.0 x 22.5 cm) |

Weighing Boats

Flexible polystyrene boats feature smooth surfaces to provide accurate pour-outs with minimal sample loss

select white, black, or blue for best color contrast to ensure sample is completely removed

biologically inert and resistant to dilute and weak acids, aqueous solutions, alcohols, and bases

temperature range -10 to 70°C (14 to 158°F)

- can be used with microwave moisture analyzers
- **all are FDA/food-grade compliant**
- **disposable**

smooth surfaces
provide accurate
pour-outs

| antistatic square weighing boats | antistatic material cat. no. | size | l x w x h | pkg. | color |
|---|------------------------------|--------|-------------------|--------|---|
| | 10770-442 | small | 46 x 46 x 8 mm | 500/pk | white <input type="checkbox"/> |
| | 10770-448 | medium | 85 x 85 x 24 mm | 500/pk | white <input type="checkbox"/> |
| | 10770-450 | large | 140 x 140 x 22 mm | 500/pk | white <input type="checkbox"/> |
| | 10770-406 | medium | 85 x 85 x 24 mm | 500/pk | blue <input checked="" type="checkbox"/> |
| | 10803-152 | large | 140 x 140 x 22 mm | 500/pk | blue <input checked="" type="checkbox"/> |
| | 10803-162 | small | 46 x 46 x 8 mm | 500/pk | blue <input checked="" type="checkbox"/> |
| | 10803-164 | large | 140 x 140 x 22 mm | 500/pk | black <input checked="" type="checkbox"/> |
| | 10803-174 | medium | 85 x 85 x 24 mm | 500/pk | black <input checked="" type="checkbox"/> |
| | 10803-176 | small | 46 x 46 x 8 mm | 500/pk | black <input checked="" type="checkbox"/> |

| square weighing boats | material cat. no. | size | l x w x h | pkg. | color |
|-----------------------------|-------------------|--------|-------------------|--------|---|
| | 10803-116 | small | 46 x 46 x 8 mm | 500/pk | blue <input checked="" type="checkbox"/> |
| | 10803-136 | small | 46 x 46 x 8 mm | 500/pk | white <input type="checkbox"/> |
| | 10803-138 | small | 46 x 46 x 8 mm | 500/pk | black <input checked="" type="checkbox"/> |
| | 10803-146 | large | 140 x 140 x 22 mm | 500/pk | black <input checked="" type="checkbox"/> |
| | 10803-148 | large | 140 x 140 x 22 mm | 500/pk | white <input type="checkbox"/> |
| | 10803-150 | large | 140 x 140 x 22 mm | 500/pk | blue <input checked="" type="checkbox"/> |
| | 10803-146 | large | 140 x 140 x 22 mm | 500/pk | black <input checked="" type="checkbox"/> |
| | 10803-168 | medium | 64 x 51 x 19 mm | 500/pk | blue <input checked="" type="checkbox"/> |
| | 10803-170 | medium | 64 x 51 x 19 mm | 500/pk | white <input type="checkbox"/> |
| | 10803-172 | medium | 64 x 51 x 19 mm | 500/pk | black <input checked="" type="checkbox"/> |

- boats are shallow with a wide, flat bottom to resist tipping
- rounded corners simplify transfer

| hexagonal weighing boats | antistatic material cat. no. | size | l x w x h | pkg. | color |
|--------------------------------|------------------------------|---------|------------------|--------|--------------------------------|
| | 10770-452 | small | 38 x 25 x 9 mm | 500/pk | white <input type="checkbox"/> |
| | 10770-454 | medium | 64 x 51 x 19 mm | 500/pk | white <input type="checkbox"/> |
| | 10770-456 | large | 102 x 76 x 22 mm | 500/pk | white <input type="checkbox"/> |
| | 10770-458 | x-large | 127 x 89 x 31 mm | 500/pk | white <input type="checkbox"/> |

- bend flexible boats at points to form pour spouts
- flat bottoms and sloping sides provide accurate transfer

| pour-boat weighing dishes | antistatic material cat. no. | size | l x w x h | pkg. | color |
|---------------------------------|------------------------------|--------|-------------------|--------|--------------------------------|
| | 10770-408 | small | 43 x 58 x 13 mm | 500/pk | white <input type="checkbox"/> |
| | 10770-438 | medium | 83 x 132 x 26 mm | 500/pk | white <input type="checkbox"/> |
| | 10770-440 | large | 108 x 183 x 26 mm | 500/pk | white <input type="checkbox"/> |

- pour-spout design easily dispenses samples without flexing the dish.
- flat-bottom dishes are made from 12-mil thick polystyrene for better strength, stability, and handling





| Diamond weighing boats | regular material cat. no. | antistatic material cat. no. | size | l x w x h | pkg. | color | |
|------------------------|---------------------------|------------------------------|--------|------------------|--------|-------|-------------------------------------|
| | 10803-114 | 10803-134 | 5 mL | 55 x 35 x 6 mm | 500/pk | White | <input type="checkbox"/> |
| | 10803-124 | 10803-142 | 25 mL | 71 x 46 x 14 mm | 500/pk | White | <input type="checkbox"/> |
| | 10803-132 | 10803-144 | 30 mL | 78 x 56 x 14 mm | 500/pk | White | <input type="checkbox"/> |
| | 10803-156 | 10803-126 | 100 mL | 119 x 90 x 19 mm | 500/pk | White | <input type="checkbox"/> |
| | 10803-110 | 10803-112 | 5 mL | 55 x 35 x 6 mm | 500/pk | Black | <input checked="" type="checkbox"/> |
| | 10803-128 | 10803-120 | 25 mL | 71 x 46 x 14 mm | 500/pk | Black | <input checked="" type="checkbox"/> |
| | 10803-130 | 10803-122 | 30 mL | 78 x 56 x 14 mm | 500/pk | Black | <input checked="" type="checkbox"/> |
| | 10803-158 | 10803-118 | 100 mL | 119 x 90 x 19 mm | 500/pk | Black | <input checked="" type="checkbox"/> |

- Diamond shape that easily flexes for accurate pouring
- Wide, flat bottom to resist tipping
- Flexible boats bend to form pour spout



Grabbit™ Temp Mitts

Protect your hands from hot or cold labware

the silicone rubber mitt is designed to comfortably fit your hand and they have internal ribs in the pockets for added dexterity

- oversized surfaces with nubs for secure grip
- mitts withstand temperatures from -70 to 500°F (-57 to 260°C)
- **autoclavable**

| cat. no. | description | l x w x h |
|-----------|----------------|-------------------|
| 89026-680 | full-hand mitt | 177 x 100 x 32 mm |

Wire Petri Dish Holder

- holds up to 6 x 100 mm petri dishes
- epoxy-coated steel wire
- **autoclavable**

| cat. no. | l x w x h |
|-----------|-------------------|
| 82008-018 | 127 x 86 x 107 mm |





Drawer Organizer

Keeps supplies organized and easily accessible

polystyrene organizer is perfect for storing and organizing small laboratory supplies

- protect small and fragile items
- fit in most standard drawers or store on benches and shelves

| cat. no. | description | overall l x w x h |
|-----------|-------------------------------------|-------------------|
| 82024-548 | 3-Compartment Tray: | 460 x 307 x 63 mm |
| | 2 compartments measure 430 x 80 mm | |
| | 1 compartment measures 430 x 110 mm | |



Acrylic Dispenser for Parafilm® Sealing Film

Heavy-duty construction for increased durability

Store, dispense, and cut Parafilm®, tape, labels, and Tough-Tags®.

- Strong, crystal-clear acrylic dispenser makes your film supply visible and features a serrated cutting blade
- Dispenser holds one or two 50 mm (2-inch) rolls, or one 100 mm (4-inch) roll

| cat. no. | W x H x D |
|-----------|--------------------|
| 82024-547 | 119 x 168 x 175 mm |



Glove Box Holders

Keep gloves organized and ready for immediate use

- free up bench top space by mounting holders on wall for easy dispensing of gloves
- all holders fit most major brands of boxes
- mounting screws included

| acrylic glove box holders | cat. no. | description | material | w x h x d |
|---------------------------|-----------|--------------|----------|--------------------|
| | 89011-720 | 1-box holder | acrylic | 141 x 260 x 110 mm |
| | 89011-722 | 2-box holder | acrylic | 295 x 260 x 110 mm |
| | 89011-724 | 3-box holder | acrylic | 295 x 385 x 110 mm |

- includes a polyester foam insert to securely hold a range of glove box sizes—snug fit ensures easy removal of gloves

| wire glove box holders | cat. no. | description | material | w x h x d |
|------------------------|-----------|--------------|----------|--------------------|
| | 89011-726 | 1-box holder | wire | 140 x 209 x 108 mm |
| | 89011-728 | 3-box holder | wire | 250 x 454 x 108 mm |

- epoxy-coated steel wire racks are an economical solution



Mini Bin™ Foil Dispenser

Safely cut foil to any size

the mini bin foil dispenser is a great alternative to wasteful giant-sized foil rolls and expensive precut foil squares

- compact to fit in limited work space areas and easy to transport
- safe to use—load the foil roll, press the cover to activate the safety guard
- foam strip prevents foil retracting back into unit
- for use with 102 mm or 153 mm (4 or 6 inch) length foil rolls (sold separately)

| cat. no. | w x h x d |
|------------------|--------------------|
| 89234-754 | 190 x 100 x 112 mm |



Benchtop Dispensing Bins

Work great with gloved hands, bright colors

ergonomically designed with large, angled openings and oversized lip

- lid remains open when tilted back
- nonskid rubber feet and pre-drilled holes for wall mounting
- green bin includes a removable divider
- mounting screws included

| cat. no. | compartments | w x h x d | colors |
|------------------|--------------|--------------------|--------------|
| 82024-532 | dual | 180 x 155 x 165 mm | neon green ■ |
| 82024-530 | single | 155 x 155 x 170 mm | neon blue ■ |



Foil Rolls

Designed for easy use and minimal waste

- custom-sized pieces of foil can easily be created for covering small surfaces such as beakers, flasks, tubes and petri dishes.
- ideal for autoclaving, freezing, incubating, protecting labware, lining, sealing, and general storage
- two sizes available. designed to fit the Mini Bin™ foil dispenser
- **disposable**

| cat. no. | description | roll length |
|------------------|--------------------------------|-------------|
| 89234-756 | 4 in (102 mm) length foil roll | 152.4 m |
| 89234-758 | 6 in (153 mm) length foil roll | 152.4 m |



Laboratory Notebooks

For recording research data, experimental procedures, and design ideas, use these notebooks to keep information safe And well organized

- numbered, acid-free, pages include spaces for title, project number, book number, date, and witness signature
- lined notebooks feature college-ruled pages for easy report writing
- covers for all versions are waterproof and chemical-resistant
- all notebooks contain periodic and conversion tables in the back

| cat. no. | pages | page format | l x w | qty. | color |
|------------------|-------|-------------|--------------|-------|---------|
| 82007-219 | 100 | grid | 218 x 281 mm | 12/cs | black ■ |
| 82007-233 | 200 | grid | 218 x 281 mm | 12/cs | blue ■ |
| 82026-407 | 200 | lined | 218 x 281 mm | 12/cs | red ■ |

Notes

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.



A series of horizontal lines for writing, spanning the width of the page. There are 20 lines in total, evenly spaced, providing a template for text entry.



A series of horizontal lines for writing, spanning the width of the page. There are 20 lines in total, evenly spaced, providing a template for text entry.



1.800.932.5000 | vwr.com

Prices and product details are current when published; subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks, of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owner(s). | Visit vwr.com to view our privacy policy, trademark owners and additional disclaimers. ©2013 VWR International, LLC. All rights reserved.