

# Certificate of Quality

VWR Catalogue Number	See List
Description	VWR® Self-Seal Sterilization Pouches
Country of Origin	Manufactured in Canada
Date of Issue (yyyy-mm-dd)	2021-09-01

VWR Catalogue Number	Dimensions
89140-800	13.3 × 25.4 cm (5 1/4 × 10")
89140-802	19.1 × 33 cm (7 1/2 × 13")
89140-804	30.5 × 45.7 cm (12 × 18")

## Quality System Compliance

Products are manufactured under the [ISO 13485:2016](#) standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP

## QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements

## Product Specifications

Material	Pouch paper	Medical grade paper 60 GSM
	Laminate film	CPP/48ga, PET
	Adhesive tape	Made of polyester film coated with medical-grade acrylic adhesive.
Sterilization	Not applicable.	
ATP Assay	Not applicable.	
DNase & RNase Free	This product is free of any detectable DNase/RNase contamination.	

<b>BSE/TSE</b>	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
<b>Cytotoxicity</b>	This product has been tested for cytotoxicity and have been shown to be non-toxic.
<b>Non-Pyrogenic Statement</b>	Not applicable.
<b>Latex Statement</b>	This product is not made from natural rubber latex.
<b>BPA Statement</b>	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
<b>DEHP Statement</b>	Not applicable.
<b>RoHS</b>	Not applicable.
<b>REACH Statement</b>	Not applicable.
<b>Storage Conditions</b>	Store at room temperature of 15 °C (59 °F) and 25 °C (77 °F) with relative humidity of 40% - 70%. Storage location should be clean, dust free, and away from direct sunlight, fluorescent or ultraviolet light.
<b>Shelf Life</b>	Shelf life of 3 years from the date of manufacture.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.  
This document has been produced electronically and is valid without a signature.