



Certificate of Conformance

This certification is provided as full assurance that the following product code and lot number was manufactured in accordance with prescribed procedures and specifications.

Product Part Number:	CR100-34WHIR
Lot Number:	007242-000
Lot Quantity:	20 CASE(S)
Date of manufacture:	05/06/2018
Shelf Life:	24 MONTHS
Description:	WHITE VINYL TAPE, IRR'D (48/CS)
Width:	3/4"
Length:	36 YARDS
Irradiation Reference Number:	0001283-00

A handwritten signature in black ink, appearing to read "John A. Smith", written over a horizontal line.

05/08/2018

Authorized Signature
Quality Control Department

Date

CC: With Product
Customer File
C of C File

3431 WEST LOMITA BOULEVARD • TORRANCE, CA 90505-5010
Tel. (310) 784-6990 • FAX (310) 784-6980

© COPYRIGHT 2015
Form No. GEN-1034
REV. No. 01

Certificate Of Processing

Prepared for MICRONOVA MFG INC

STERIS

Gamma Process Run ID 108578A

Product Code

CR100-34WHIR

Product Lot Number

007242-000

Quantity

20

UOM

CS

Other Information:

0001283-00

Processing Run Start Date/Time: 06-May-2018 04:39:30 am

Approx. Downtime (hours): 0.95

Processing Run End Date/Time: 06-May-2018 09:36:38 am

Minimum Specified Dose (kGy): 40.0

Minimum Delivered Dose (kGy): 43.1

Maximum Specified Dose (kGy): 60.0

Maximum Delivered Dose (kGy): 54.9

Product meets Customer specifications; zero nonconformities occurred during this irradiation run.

Signature Manifest

Reviewed and E-Signed By

Adrian Martinez (QS & RC Technician)

Document Content Revision: 1

Signed On 5/7/2018 at 5:01 AM

UTC / GMT Offset (hh:mm): -7:00

Processing Location:

STERIS

1000 South Sarah Place

Ontario, CA 91761

Phone: 909-390-9942

Fax: 909-390-7854

Operating facilities are in compliance with applicable state and federal regulations (FDA, NRC, EPA, and OSHA) and provide services under a quality system which meets the requirements of FDA QSR, EN/ISO 13485, and in alignment with EN ANSI/AAMI/ISO 11137. STERIS certifies that these processed items received the indicated doses within the precision and accuracy of the dosimetry system used.