



Anachemia
A VWR Company

Product Specification

Material	74796-300
Material description	Potassium persulfate
Grade	ACS
CAS Number	7727-21-1
Molecular formula	K ₂ S ₂ O ₈
Molecular mass	270,32

Characteristics	Specifications
Assay	Min. 99,0 %
Chlorine compounds (as Cl)	Max. 0,001 %
Heavy metals (as Pb)	Max. 0,001 %
Insoluble matter	Max. 0,005 %
Fe (Iron)	Max. 5 ppm
Mn (Manganese)	Max. 2 ppm

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Michel Boisvert, Head of the Laboratory - Lachine
VWR Anachemia Canada, 255 rue Norman, Lachine, H8R 1A3,
QC, Canada