



Anachemia
A VWR Company

Product Specification

Material	52716-140
Material description	Lithium nitrate, crystallized
Grade	Reagent
CAS Number	7790-69-4
Molecular formula	LiNO ₃
Molecular mass	68.95

Characteristics	Specifications
Assay (calculated on dried substance)	Min. 99.0 %
Heavy metals (as Pb)	Max. 10 ppm
Insoluble matter	Max. 50 ppm
Cl (Chloride)	Max. 20 ppm
SO ₄ (Sulphate)	Max. 100 ppm
Fe (Iron)	Max. 10 ppm

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Michel Boisvert, Head of the Laboratory - Lachine
VWR Anachemia Canada, 255 rue Norman, Lachine, H8R 1A3,
QC, Canada

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.