

Biology

From the highest quality living and preserved specimens to the widest selection of media, microscope slides, and more, we have more options for your biology classroom than anywhere else.



Image Guided Injection Trainer for	Lumbar Spine	Thoracic Spine	Cervical Spine
Product number	470311-930	470318-134	470318-132
Skeletal Anatomy for imaging and palpation	Sacrum, S1-S3 (with sacral hiatus canal)	Vertebrae T3-T8	Occipital bone (cranial to external occipital protuberance)
	Ilium, bilateral without ischial bone and hip joint	Ribs 3-8	Vertebrae C1-T2
	Vertebrae T12-L5		Ribs 1 and 2
	Coccyx		
Interventional spine procedures	Transforaminal Epidural Steroid Injections (TFSI)	Interlaminar Epidural Steroid Injection	Greater Occipital Nerve (GON) Infiltration
	Intralaminar Epidural Injection	Thoracic Transforaminal Injection	Transforaminal Epidural Steroid Injections (TFSI)
	Facet Blocks	Thoracic Zygapophysial Joint	Intralaminar Epidural Injection
	Medial Branch Block (MBB/RF)	Nerve (Medial Branch) Injection	Medial Branch Block
	S1-block	Intraarticular Injection	Facet radiofrequency denervation
	Sacroiliac Joint Injection (SIJ-injection)	Thoracic Zygapophysial Joint	C1/C2 intraarticular
		Intercostal Nerve Block (ICNB)	Cervical facet intraarticular