

Product Name:	PF021WWW/0M VWR #76271-332
Description:	<p>White Diamond Series Countertop Freezer Freestanding</p> <p>1.7 Cu. Ft. Capacity</p> <p>One solid swing door, right hinged*</p> <p>Microprocessor temperature controller</p> <p>Digital temperature display</p> <p>Audible and visual high and low temperature alarms</p> <p>Remote alarm contacts</p> <p>2 probes; 1 in air and 1 in glycol bottle</p> <p>Adjustable operating temperature range: -15°C – -25°C</p> <p>2 adjustable shelves with guard rail to prevent product from being pushed against back wall</p> <p>Bottom footed shelf to allow air circulation</p> <p>Keyed door lock</p> <p>Magnetic door gasket for positive seal</p> <p>Leveling legs</p> <p>LED interior light</p> <p>Integrated door handle</p> <p>Probe access port (3/8")</p> <p>Powder coat exterior</p> <p>High density urethane foam cabinet and door insulation</p> <p>HFC-free refrigerant (R600a), foam insulation and packaging</p> <p>Three year parts and labor warranty, plus an additional two year compressor parts only warranty</p> <p>Pyxis®, Omnicell® and AcuDose RX® compatible</p> <p>Exterior dimensions: 23 ¾ "W x 24" D x 20 3/8"H</p> <p><i>Note: This unit must have 4" clearance on sides and back for adequate ventilation</i></p> <p>Manual defrost</p> <p>115V, 60 Hz, 1 Amp, 1/8 HP</p> <p>UL/C-UL listed</p> <p>Shipping weight: 106 lbs.</p> <p>*This product is available left hinged, product number PF021WWWLH/0M</p>